 <b>GAUTENG PROVINCE</b> PROVINCIAL TREASURY REPUBLIC OF SOUTH AFRICA	<b>PROVINCIAL SUPPLY CHAIN MANAGEMENT</b>	
	<b>PRICING SCHEDULE – NON - FIRM PRICES (PURCHASES)</b>	<b>Page: 1 of 4</b>

**NOTE: PRICE ADJUSTMENTS WILL BE ALLOWED AT THE PERIODS AND TIMES SPECIFIED IN THE BIDDING DOCUMENTS.**

**IN CASES WHERE DIFFERENT DELIVERY POINTS INFLUENCE THE PRICING, A SEPARATE PRICING SCHEDULE MUST BE SUBMITTED FOR EACH DELIVERY POINT**


<b>NAME OF THE BIDDER</b>		<b>BID NUMBER</b>	
---------------------------	--	-------------------	--

<b>CLOSING TIME</b>		<b>CLOSING DATE</b>	
---------------------	--	---------------------	--

<b>OFFER TO BE VALID FOR</b>		<b>DAYS FROM THE CLOSING DATE OF BID</b>
------------------------------	--	--

ITEM NO	QUANTITY	DESCRIPTION	BID PRICE IN RSA CURRENCY ** (ALL APPLICABLE TAXES INCLUDED)

<b>REQUIRED BY:</b>	
<b>AT:</b>	
<b>BRAND MODEL:</b>	
<b>COUNTRY OF ORIGIN:</b>	

 <b>GAUTENG PROVINCE</b> PROVINCIAL TREASURY REPUBLIC OF SOUTH AFRICA	<b>PROVINCIAL SUPPLY CHAIN MANAGEMENT</b>	
	<b>PRICING SCHEDULE – NON - FIRM PRICES (PURCHASES)</b>	<b>Page: 2 of 4</b>


Does the offer comply with the specification(s)?	YES		NO	
--	-----	--	----	--

If not to specification, indicate deviation(s)

<b>PERIOD REQUIRED FOR DELIVERY</b>	
-------------------------------------	--

<b>DELIVERY</b>	<b>FIRM</b>		<b>NOT FIRM</b>	
-----------------	-------------	--	-----------------	--

**“All applicable taxes” includes value- added tax, pay as you earn, income tax, unemployment insurance fund contributions and skills development levies.**

 <b>GAUTENG PROVINCE</b> PROVINCIAL TREASURY REPUBLIC OF SOUTH AFRICA	<b>PROVINCIAL SUPPLY CHAIN MANAGEMENT</b>	
	<b>PRICING SCHEDULE – NON - FIRM PRICES (PURCHASES)</b>	<b>Page: 3 of 4</b>

## PRICE ADJUSTMENTS

### A NON-FIRM PRICES SUBJECT TO ESCALATION

1. IN CASES OF PERIOD CONTRACTS, NON-FIRM PRICES WILL BE ADJUSTED (LOADED) WITH THE ASSESSED CONTRACT PRICE ADJUSTMENTS IMPLICIT IN NON-FIRM PRICES WHEN CALCULATING THE COMPARATIVE PRICES
2. IN THIS CATEGORY PRICE ESCALATIONS WILL ONLY BE CONSIDERED IN TERMS OF THE FOLLOWING FORMULA:


$$Pa = (1 - V)Pt \left( D1 \frac{R1t}{R1o} + D2 \frac{R2t}{R2o} + D3 \frac{R3t}{R3o} + D4 \frac{R4t}{R4o} \right) + VPt$$

Where:

Pa	=	The new escalated price to be calculated.
(1-V) Pt	=	85% of the original bid price. <b>Note that Pt must always be the original bid price and not an escalated price.</b>
D1, D2..	=	Each factor of the bid price eg. labour, transport, clothing, footwear, etc. The total of the various factors D1, D2...etc. must add up to 100%.
R1t, R2t.....	=	Index figure obtained from new index (depends on the number of factors used).
R1o, R2o	=	Index figure at time of bidding.
VPt	=	15% of the original bid price. This portion of the bid price remains firm i.e. it is not subject to any price escalations.

3. THE FOLLOWING INDEX/INDICES MUST BE USED TO CALCULATE YOUR BID PRICE:

INDEX	DATED	INDEX	DATED	INDEX	DATE

 <b>GAUTENG PROVINCE</b> PROVINCIAL TREASURY REPUBLIC OF SOUTH AFRICA	<b>PROVINCIAL SUPPLY CHAIN MANAGEMENT</b>	
	<b>PRICING SCHEDULE – NON - FIRM PRICES (PURCHASES)</b>	<b>Page: 4 of 4</b>

4. FURNISH A BREAKDOWN OF YOUR PRICE IN TERMS OF ABOVE-MENTIONED FORMULA. THE TOTAL OF THE VARIOUS FACTORS MUST ADD UP TO 100%.

FACTOR (DA, D2, etc. eg. Labour, Transport etc.)	PERCENTAGE OF BID PRICE

#### B. PRICES SUBJECT TO RATE OF EXCHANGE VARIATIONS

1. Please furnish full particulars of your financial institution, state the currencies used in the conversion of the prices of the items to South African currency, which portion of the price is subject to rate of exchange variations and the amounts remitted abroad.

PARTICULARS OF FINANCIAL INSTITUTION	ITEM NO	PRICE	CURRENCY	RATE	PORTION OF PRICE SUBJECT TO ROE	AMOUNT IN FOREIGN CURRENCY REMITTED ABROAD
				ZAR=		
				ZAR=		
				ZAR=		

2. Adjustments for rate of exchange variations during the contract period will be calculated by using the average monthly exchange rates as issued by your commercial bank for the periods indicated hereunder: (Proof from bank required)

AVERAGE MONTHLY EXCHANGE RATES FOR THE PERIOD	DATE DOCUMENTATION MUST BE SUBMITTED TO THIS OFFICE	DATE FROM WHICH NEW CALCULATED PRICES WILL BECOME EFFECTIVE	DATE UNTIL WHICH NEW CALCULATED PRICE WILL BE EFFECTIVE

GT/GDH/036/2025 THE SUPPLEMENTARY FOR GT/GDH/002/2024 – APPOINTMENT OF SERVICE PROVIDER/S TO SUPPLY AND DELIVER REHABILITATION CONSUMABLES TO THE GDOH INSTITUTIONS FOR A PERIOD ENDING 15 NOVEMBER 2027																			
											PROVIDE A COST BREAKDOWN AS A PERCENTAGE OF TOTAL UNIT PRICE								IS PRODUCT IMPORTED OR MANUFACTURED LOCALLY?
ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
A: - PHYSIOTHERAPY																			
ITEM 1	Therapeutic Hand Putty																		
ITEM-1.1	158919	<b>Therapeutic hand Putty</b> Soft resistance • Used for therapeutic exercises for the hand or finger/s • Offers soft resistance for strengthening of muscles of the finger/s or hand • Non-Toxic • Non-Oily Smooth and pliable <b>Dimensions:</b> • 200g quantities <b>Packaging requirements:</b> • Packaged in re-sealable airtight durable containers <b>Estimated usage:</b> 0	Yes, representative sample for items 1.1 and 1.2																
ITEM-1.2	177796	<b>Therapeutic hand Putty</b> Soft resistance • Used for therapeutic exercises for the hand or finger/s • Offers soft resistance for strengthening of muscles of the finger/s or hand • Non-Toxic • Non-Oily Smooth and pliable <b>Dimensions:</b> •500g quantities <b>Packaging requirements:</b> • Packaged in re-sealable airtight durable containers <b>Estimated usage:</b> 22	Yes, representative sample for items 1.1 and 1.2																
ITEM-1.3	159010	<b>Therapeutic hand putty</b> Medium resistance • Used for therapeutic exercises for the hand or finger/s • Offers medium resistance for strengthening of muscles of the finger/s or hand • Non-toxic • Non-Oily Smooth and pliable <b>Dimensions:</b> 500g quantities <b>Packaging requirements:</b> Packaged in re-sealable airtight durable containers <b>Estimated usage:</b> 7	Yes, representative sample for items 1.3 and 1.4		Each														
ITEM-1.4	177797	<b>Therapeutic hand putty</b> Medium resistance • Used for therapeutic exercises for the hand or finger/s • Offers medium resistance for strengthening of muscles of the finger/s or hand • Non-toxic • Non-Oily Smooth and pliable <b>Dimensions:</b> 200g quantities <b>Packaging requirements:</b> Packaged in re-sealable airtight durable containers <b>Estimated usage:</b> 1	Yes, representative sample for items 1.3 and 1.4		Each														
ITEM-1.5	177798	<b>Therapeutic hand putty</b> Firm resistance • Used for therapeutic exercises for the hand or finger/s • Offers firm resistance for strengthening of muscles of the finger/s or hand • Non-toxic • Non-Oily Smooth and pliable <b>Dimensions:</b> 500g quantities <b>Packaging requirements:</b> Packaged in re-sealable airtight durable containers <b>Estimated usage:</b> 1	Yes, representative sample for items 1.5 and 1.6		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-1.6	159011	<b>Therapeutic hand putty</b> Firm resistance • Used for therapeutic exercises for the hand or finger/s • Offers firm resistance for strengthening of muscles of the finger/s or hand • Non-toxic • Non-Oily Smooth and pliable <b>Dimensions:</b> 200g quantities <b>Packaging requirements:</b> Packaged in re-sealable airtight durable containers <b>Estimated usage:</b> 1	Yes, representative sample for items 1.5 and 1.6		Each														
ITEM 2	<b>Resistive Exercise Band</b>																		
ITEM-2.1	159012	<b>Resistive Exercise Band</b> – medium resistance • Resistive exercise band used for therapeutic strengthening of muscles • Offers medium resistance to stretch • Made of elastomer • No odour, Latex free <b>Dimensions:</b> Width = 10-15cm, Length of roll = 25m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 16	Yes, representative sample for items 2.1 and 2.2		Per Roll														
ITEM-2.2	177807	<b>Resistive Exercise Band</b> – medium resistance • Resistive exercise band used for therapeutic strengthening of muscles • Offers medium resistance to stretch • Made of elastomer • No odour, Latex free <b>Dimensions:</b> Width = 10-15cm, Length of roll = 45m - 50m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 2.1 and 2.2		Per Roll														
ITEM-2.3	159013	<b>Resistive Exercise Band</b> – heavy resistance • Resistive exercise band used for therapeutic strengthening of muscles • Offers heavy resistance to stretch • Made of elastomer • No odour, Latex free <b>Dimensions:</b> Width = 10-15cm, Length of roll = 25m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 2.3 and 2.4		Per Roll														
ITEM-2.4	178716	<b>Resistive Exercise Band</b> – heavy resistance • Resistive exercise band used for therapeutic strengthening of muscles • Offers heavy resistance to stretch • Made of elastomer • No odour, Latex free <b>Dimensions:</b> Width = 10-15cm, Length of roll = 45m - 50m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 2.3 and 2.4		Per Roll														
ITEM-2.5	159014	<b>Resistive Exercise Band</b> – extra heavy resistance • Resistive exercise band used for therapeutic strengthening of muscles • Offers extra heavy resistance to stretch • Made of elastomer • No odour, Latex free <b>Dimensions:</b> Width = 10-15cm, Length of roll = 25m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 18	Yes, representative sample for items 2.5 and 2.6		Per Roll														
ITEM-2.6	178717	<b>Resistive Exercise Band</b> – extra heavy resistance • Resistive exercise band used for therapeutic strengthening of muscles • Offers extra heavy resistance to stretch • Made of elastomer •No odour, Latex free <b>Dimensions:</b> Width = 10-15cm, Length of roll = 45m - 50m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 29	Yes, representative sample for items 2.5 and 2.6		Per Roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-2.7	159016	<b>Resistive Exercise Band</b> – double extra heavy resistance <ul style="list-style-type: none"><li>• Resistive exercise band used for therapeutic strengthening of muscles</li><li>• Offers double extra heavy resistance to stretch</li><li>• Made of elastomer</li><li>• No odour, Latex free</li></ul> <b>Dimensions:</b> Width = 10-15cm, Length of roll = 25m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 7	Yes		Per Roll														
ITEM-2.8	113421	<b>Resistive Exercise Band</b> – light resistance <ul style="list-style-type: none"><li>• Resistive exercise band used for therapeutic strengthening of muscles</li><li>• Offers light resistance to stretch</li><li>• Made of elastomer</li><li>• No odour, Latex free</li></ul> <b>Dimensions:</b> Width = 10-15cm, Length of roll = 25m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 2.8 and 2.9		Per Roll														
ITEM-2.9	156602	<b>Resistive Exercise Band</b> – light resistance <ul style="list-style-type: none"><li>• Resistive exercise band used for therapeutic strengthening of muscles</li><li>• Offers light resistance to stretch</li><li>• Made of elastomer</li><li>• No odour, Latex free</li></ul> <b>Dimensions:</b> Width = 10-15cm, Length of roll = 45m - 50m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 7	Yes, representative sample for items 2.8 and 2.9		Per Roll														
ITEM 3	<b>Exercise Band/Tubing Handles</b>																		
ITEM-3.1	159480	<b>Exercise Band/Tubing Handles</b> <ul style="list-style-type: none"><li>• Handle for exercise band/tubing</li><li>• Easy grip handle</li><li>• Exercise band/tubing can be easily attached to handle</li></ul> <b>Dimensions:</b> To fit size of adult hand <b>Packaging requirements:</b> One Pair <b>Estimated usage:</b> 4	Yes		Pair														
ITEM-3.2	174316	<b>Medium Resistant Exercise Tubing</b> <ul style="list-style-type: none"><li>• Resistive exercise tubing used for therapeutic strengthening of muscles</li><li>• Offers medium resistance to stretch</li></ul> <b>Dimensions:</b> Length of Roll: 25m <b>Packaging requirements:</b> Single Box <b>Estimated usage:</b> 1	Yes		Box														
ITEM-3.3	174312	<b>Heavy Resistant Exercise Tubing</b> <ul style="list-style-type: none"><li>• Resistive exercise tubing used for therapeutic strengthening of muscles</li><li>• Offers heavy resistance to stretch</li></ul> <b>Dimensions:</b> Length of Roll: 25m <b>Packaging requirements:</b> Single Box <b>Estimated usage:</b> 2	Yes		Box														
ITEM-3.4	174406	<b>Extra Heavy Resistant Exercise Tubing</b> <ul style="list-style-type: none"><li>• Resistive exercise tubing used for therapeutic strengthening of muscles</li><li>• Offers extra heavy resistance to stretch</li></ul> <b>Dimensions:</b> Length of Roll: 25m <b>Packaging requirements:</b> Single Box <b>Estimated usage:</b> 1	Yes		Box														
ITEM 4	<b>Strapping</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-4.1	159481	<b>Strapping – underwrap</b> <ul style="list-style-type: none"><li>• Used in conjunction with strapping</li><li>• Made of non-woven, breathable (porous) material which allows some stretch</li><li>• Adhesive-should be able to remain on the skin for up to one week</li><li>• Hypoallergenic</li><li>• Printed grid on backing sheet</li></ul> <b>Dimensions:</b> Width = 5cm, Length of roll ±10m <b>Packaging requirements:</b> Single roll in a Box <b>Estimated usage:</b> 142	Yes		Per Roll in a box														
ITEM-4.2	159482	<b>Strapping – rigid tape</b> <ul style="list-style-type: none"><li>• Used in conjunction with underwrap</li><li>• Used for therapeutic strapping of muscles and or joints</li><li>• Aggressively adhesive rigid strapping tape</li><li>• Non-stretch</li><li>• Made of strong rayon material, spread with zinc oxide material</li><li>• Porous</li><li>• Should be able to remain on the skin for up to one week</li><li>• Flesh coloured</li><li>• Hypoallergenic</li></ul> <b>Dimensions:</b> Width 3.8cm, Length of roll ±13.7m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 49	Yes		Per Roll in a box														
ITEM-4.3	159484	<b>Strapping – sports tape</b> <ul style="list-style-type: none"><li>• Used for therapeutic strapping of muscles and or joints especially ankles</li><li>• Aggressively adhesive rigid strapping tape</li><li>• Non-stretch</li><li>• Made of strong rayon material, spread with zinc oxide material</li><li>• Porous</li><li>• Remains in place during vigorous activities</li><li>• Hypoallergenic</li><li>• Hand-tearable</li></ul> <b>Dimensions:</b> Width 3.8cm, Length of roll ±13.7m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 40	Yes		Per Roll in a box														
ITEM 5	Kinesiotape																		
ITEM-5.1	118875	<b>Kinesiotape, 5cm wide</b> <ul style="list-style-type: none"><li>- Therapeutic elastic adhesive strapping tape</li><li>- latex free</li><li>- anti-allergy/hypoallergenic</li><li>- wave patterned adhesive design</li><li>- water resistant</li><li>- stretchable</li><li>- adhesive design that lifts skin to maintain flexibility, increase circulation, decrease pain, re-educate muscles</li><li>- Known as Kinesiotape or Kinesiology Tape</li></ul> <b>Dimensions:</b> Width 5cm, Length 5m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 5.1 and 5.2		Per Roll in a box														



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-5.2	124249	<b>Kinesiotape, 7.5cm wide</b> - Therapeutic elastic adhesive strapping tape - latex free - anti-allergy/hypoallergenic - wave patterned adhesive design - water resistant - stretchable - adhesive design that lifts skin to maintain flexibility, increase circulation, decrease pain, re-educate muscles - Known as Kinesiotape or Kinesiology Tape <b>Dimensions:</b> Width 7.5cm, Length 5m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 5.1 and 5.2		Per Roll in a box														
ITEM 6	<b>Dynamic tape</b>																		
ITEM-6.1	New	<b>Dynamic tape, 5cm</b> • Extremely soft, highly conformable, breathable and water resistant • Allows full range of movement without limitation • Four way stretch <b>Dimensions:</b> 5cm x 5m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 6.1 and 6.2		Each														
ITEM-6.2	New	<b>Dynamic tape, 7.5cm</b> • Extremely soft, highly conformable, breathable and water resistant • Allows full range of movement without limitation • Four way stretch <b>Dimensions:</b> 7.5cm x 5m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 6.1 and 6.2		Each														
ITEM 7	<b>Compression bandage short stretch</b>																		
ITEM-7.1	New	<b>Compression bandage short stretch, 6cm</b> • Short stretch textile elastic • Permeable, latex free, 100% cotton • Flexible, high wearing comfort and re-usable, can used for management of lymphedema <b>Dimensions:</b> +- 6cm x 5m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 7.1 to 7.4		Each														
ITEM-7.2	New	<b>Compression bandage short stretch, 8cm</b> • Short stretch textile elastic • Permeable, latex free, 100% cotton • Flexible, high wearing comfort and re-usable, can used for management of lymphedema <b>Dimensions:</b> 8cm x 5m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 7.1 to 7.4		Each														
ITEM-7.3	New	<b>Compression bandage short stretch, 10cm</b> • Short stretch textile elastic • Permeable, latex free, 100% cotton • Flexible, high wearing comfort and re-usable, can used for management of lymphedema <b>Dimensions:</b> 10cm x 5m <b>Packaging requirements:</b> Single roll in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 7.1 to 7.4		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-7.4	New	<b>Compression bandage short stretch, 12cm</b> <ul style="list-style-type: none"><li>• Short stretch textile elastic</li><li>• Permeable, latex free, 100% cotton</li><li>• Flexible, high wearing comfort and re-usable, can used for management of lymphedema</li></ul> <b>Dimensions:</b> 12cm x 5m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 7.1 to 7.4		Each														
ITEM-7.5	New	<b>Open cell foam bandage for padding</b> <ul style="list-style-type: none"><li>• Bandage for padding under compression bandaging</li><li>• 100% polyurethane, latex free</li><li>• Air and moisture permeable.</li></ul> Non-sliding and tear resistant. Re-usable with same patient <b>Dimensions:</b> 10cm x 2.5m x 0.4cm <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes		Each														
ITEM 8	<b>Elastic Tubular Support Bandage</b>																		
ITEM-8.1	New	<b>Elastic tubular support bandage, 4.2cm</b> <ul style="list-style-type: none"><li>• Elasticated tubular bandage that provides firm and comfortable support, 100% cotton</li><li>• Reduces persistent swelling</li><li>• Easy to apply, re-usable with same patient, easy to wash and latex free.</li></ul> <b>Dimensions:</b> 4.2cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 8.1 to 8.8		Each														
ITEM-8.2	New	<b>Elastic tubular support bandage, 6.25cm</b> <ul style="list-style-type: none"><li>• Elasticated tubular bandage that provides firm and comfortable support, 100%</li><li>• Reduces persistent swelling</li><li>• Easy to apply, re-usable with same patient, easy to wash and latex free.</li></ul> <b>Dimensions:</b> 6.25cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 8.1 to 8.8		Each														
ITEM-8.3	New	<b>Elastic tubular support bandage, 6.75cm</b> <ul style="list-style-type: none"><li>• Elasticated tubular bandage that provides firm and comfortable support, 100% cotton</li><li>• Reduces persistent swelling</li><li>• Easy to apply, re-usable with same patient, easy to wash and latex free.</li></ul> <b>Dimensions:</b> 6.75cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 8.1 to 8.8		Each														
ITEM-8.4	New	<b>Elastic tubular support bandage, 7.5cm</b> <ul style="list-style-type: none"><li>• Elasticated tubular bandage that provides firm and comfortable support, 100% cotton</li><li>• Reduces persistent swelling</li><li>• Easy to apply, re-usable with same patient, easy to wash and latex free.</li></ul> <b>Dimensions:</b> 7.5cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 8.1 to 8.8		Each														
ITEM-8.5	New	<b>Elastic tubular support bandage, 8.75cm</b> <ul style="list-style-type: none"><li>• Elasticated tubular bandage that provides firm and comfortable support, 100% cotton</li><li>• Reduces persistent swelling</li><li>• Easy to apply, re-usable with same patient, easy to wash and latex free.</li></ul> <b>Dimensions:</b> 8.75cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage:</b> 1	Yes, representative sample for items 8.1 to 8.8		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-8.6	New	<b>Elastic tubular support bandage, 10cm</b> • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <b>Dimensions:</b> 10cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage: 0</b>	Yes, representative sample for items 8.1 to 8.8		Each														
ITEM-8.7	New	<b>Elastic tubular support bandage, 12cm</b> • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <b>Dimensions:</b> 12cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage:0</b>	Yes, representative sample for items 8.1 to 8.8		Each														
ITEM-8.8	New	<b>Elastic tubular support bandage, 14cm</b> • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <b>Dimensions:</b> 14cm x 10m <b>Packaging requirements:</b> Single roll <b>Estimated usage:0</b>	Yes, representative sample for items 8.1 to 8.8		Each														
<b>ITEM 9</b>																			
ITEM-9.1	159541	<b>Paraffin wax refill, beads</b> For therapeutic wax bath • Wax beads • Scent free • Paraffin wax mixed with paraffin oil <b>Dimensions:</b> 500g <b>Packaging requirements:</b> Sealed packaging <b>Estimated usage: 7</b>	Yes, representative sample for items 9.1 and 9.2		Pkt														
ITEM-9.2	113522	<b>Paraffin wax refill, beads</b> For therapeutic wax bath • Wax beads • Scent free • Paraffin wax mixed with paraffin oil <b>Dimensions:</b> 1kg <b>Packaging requirements:</b> Sealed packaging <b>Estimated usage: 7</b>	Yes, representative sample for items 9.1 and 9.2		Pkt														
ITEM-9.3	159541	<b>Paraffin wax refill, block</b> For therapeutic wax bath • Wax block • Scent free • Paraffin wax mixed with paraffin oil <b>Dimensions:</b> 500g <b>Packaging requirements:</b> Sealed packaging <b>Estimated usage: 7</b>	Yes, representative sample for items 9.3 and 9.4		Pkt														
ITEM-9.4	136424	<b>Paraffin wax refill, block</b> For therapeutic wax bath • Wax Block • Scent free • Paraffin wax mixed with paraffin oil <b>Dimensions:</b> 1kg <b>Packaging requirements:</b> Sealed packaging <b>Estimated usage: 7</b>	Yes, representative sample for items 9.3 and 9.4		Pkt														
<b>ITEM 10</b>																			
Arnica Oil																			

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-10.1	159542	<b>Arnica Oil – 250 ml</b> <ul style="list-style-type: none"><li>• Oil containing Arnica concentrate</li><li>• Contain essential oils which promotes muscular relaxation</li><li>• Easily removable</li></ul> <b>Dimensions:</b> 250 ml <b>Packaging:</b> Per Bottle <b>Estimated usage:</b> 6	Yes, representative sample for items 10.1 to 10.3		Per bottle														
ITEM-10.2	159543	<b>Arnica Oil – 500 ml</b> <ul style="list-style-type: none"><li>• Massage Oil containing Arnica concentrate</li><li>• Contain essential oils which promotes muscular relaxation</li><li>• Easily removable</li></ul> <b>Dimensions:</b> 500 ml <b>Packaging:</b> Per Bottle <b>Estimated usage:</b> 11	Yes, representative sample for items 10.1 to 10.3		Per bottle														
ITEM-10.3	159544	<b>Arnica Oil – 1L</b> <ul style="list-style-type: none"><li>• Massage Oil containing Arnica concentrate</li><li>• Contain essential oils which promotes muscular relaxation</li><li>• Easily removable</li></ul> <b>Dimensions:</b> 1L <b>Packaging:</b> Per Bottle <b>Estimated usage:</b> 13	Yes, representative sample for items 10.1 to 10.3		Per bottle														
ITEM 11		<b>Acupuncture needles</b>																	
ITEM-11.1	159546	<b>Acupuncture needles – 13 mm</b> <ul style="list-style-type: none"><li>• Copper handled needles with guide tubes</li><li>• Disposable: for single use only</li></ul> <b>Dimensions:</b> Length = 13mm, Gauge = 0.22mm <b>Packaging requirements:</b> Sterile individually packed, 100 needles per box <b>Estimated usage:</b> 31	Yes		Per box														
ITEM-11.2	159611	<b>Acupuncture needles – 25mm</b> <ul style="list-style-type: none"><li>• Copper handled needles with guide tubes</li><li>• Disposable: for single use only</li></ul> <b>Dimensions:</b> Length = 25mm, Gauge = 0.25mm <b>Packaging requirements:</b> Sterile individually packed, 100 needles per box <b>Estimated usage:</b> 38	Yes		Per box														
ITEM-11.3	162854	<b>Acupuncture needles – 40 mm</b> <ul style="list-style-type: none"><li>• Copper handled needles with guide tubes</li><li>• Disposable: for single use only</li></ul> <b>Dimensions:</b> Length = 40mm, Gauge = 0.30mm <b>Packaging requirements:</b> Sterile individually packed, 100 needles per box <b>Estimated usage:</b> 19	Yes		Per box														
ITEM-11.4	159612	<b>Acupuncture needles – 50 mm</b> <ul style="list-style-type: none"><li>• Copper handled needles with guide tubes</li><li>• Disposable: for single use only</li></ul> <b>Dimensions:</b> Length = 50mm, Gauge = 0.30mm <b>Packaging requirements:</b> Sterile individually packed, 100 needles per box <b>Estimated usage:</b> 28	Yes		Per box														
ITEM-11.5	136425	<b>Acupuncture needles – 30 mm</b> <ul style="list-style-type: none"><li>• Copper handled needles with guide tubes</li><li>• Disposable: for single use only</li></ul> <b>Dimensions:</b> Length = 0,25mm, Gauge = 0.30-35mm <b>Packaging requirements:</b> Sterile individually packed, 100 needles per box <b>Estimated usage:</b> 39	Yes		Per box														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-11.6	159613	<b>Acupuncture needles – 75 mm</b> • Copper handled needles with guide tubes • Disposable: for single use only <b>Dimensions:</b> Length = 75mm, Gauge = 0.30-35mm <b>Packaging requirements:</b> Sterile individually packed, 100 needles per box <b>Estimated usage:</b> 39	Yes		Per box														
ITEM-11.7	159614	<b>Acupuncture needles –100 mm</b> • Copper handled needles with guide tubes • Disposable: for single use only <b>Dimensions:</b> Length = 100mm, Gauge = 0.30-35mm <b>Packaging requirements:</b> Sterile individually packed, 100 needles per box <b>Estimated usage:</b> 28	Yes		Per box														
ITEM 12		<b>Flutter / Broncho-vibe</b>																	
ITEM-12.1	160068	<b>Flutter / Broncho-vibe</b> • A device which generates controlled oscillating positive pressure during expiration • Aids in mobilizing secretions • Hard, plastic cylindrical mouth piece around a circular core of hardened material • Contains a high density stainless steel ball • Has a perforated protective cover • <b>Washable</b> <b>Dimensions:</b> Length ± 15cm, Diameter ± 3cm <b>Packaging requirements:</b> Each in a sealed pack <b>Estimated usage:</b> 92	Yes		Each														
ITEM 13		<b>Luft Respiratory Device</b>																	
ITEM-13.1	New	<b>Luft Respiratory Device</b> • Resistance training for the lungs • Assists with breathing • Hand held device • Re-usable with same patient <b>Dimensions:</b> None <b>Packaging requirements:</b> Each in box <b>Estimated usage:</b> 1	Yes		each														
ITEM 14		<b>Positive Expiratory Pressure</b>																	
ITEM-14.1	137067	<b>Positive Expiratory Pressure (PEP)</b> - Made of transparent hard plastic - Comes with removable mouth-piece - Small volume nebulizer - Collapsible reservoir tubing which must be able to attach to the nebulizer - It must be delivered with all accessories - It must be autoclavable <b>Packaging:</b> Each in a sealed plastic pack <b>Estimated usage :</b> 1	Yes		Each														
ITEM15		<b>Respi-flow Respiratory Exerciser</b>																	
ITEM-15.1	137068	<b>Respi-flow Respiratory Exerciser</b> - Must have 3 chambers - 3 balls - Disposed after a single patient use - Made of transparent plastic with removable re-usable mouth piece - Flow maximum inspiration - Help to develop & improve respiratory fitness <b>Packaging:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 16		<b>Transparent cling wrap</b>																	

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-16.1	160125	<b>Transparent cling wrap</b> • Transparent • Self attaching • Non breathable • Used for wrapping the body or parts thereof during dermatology treatments to facilitate uptake of preparations into the skin <b>Dimensions:</b> Length ± 40m, Width ± 30cm <b>Diameter:</b> ± 3cm <b>Packaging requirements:</b> Sealed packaging, Single item <b>Estimated usage:</b> 177	Yes		Each														
ITEM 17	Rectangular sponge sleeves																		
ITEM-17.1	160126	<b>Rectangular sponge sleeves, medium</b> For carbon rubber electrodes of interferential machine • Durable, high quality sponge 3 Layers thick <b>Dimensions:</b> To fit rubber electrode of 60 x 80 mm <b>Packaging requirements:</b> Sets of 4 <b>Estimated usage:</b> 30	Yes, representative sample for items 17.1 and 17.2		Per set														
ITEM-17.2	160129	<b>Rectangular sponge sleeves, large</b> For carbon rubber electrodes of interferential machine • Durable, high quality sponge 3 Layers thick <b>Dimensions:</b> To fit rubber electrode of 80 x 120 mm <b>Packaging requirements:</b> Sets of 4 <b>Estimated usage:</b> 1	Yes, representative sample for items 17.1 and 17.2		Per set														
ITEM 18	Electrodes																		
ITEM-18.1	160131	<b>Reusable round/oval gel electrodes – 50mm</b> • Solid gel with cover • Suitable for multiple applications for 10-15 days • Self adhering • Low profile with pigtail pin connection <b>Dimensions:</b> = ± 50mm <b>Packaging requirements:</b> Set of 4 <b>Estimated usage:</b> 28	Yes		Per set														
ITEM-18.2	160428	<b>Reusable rectangular gel electrodes – 50 x 90mm</b> • Solid gel with cover • Suitable for multiple applications for 10-15 days • Self adhering • Low profile with pigtail pin connection <b>Dimensions:</b> = 50mm x ± 90mm <b>Packaging requirements:</b> Set of 4 <b>Estimated usage:</b> 56	Yes		Per set														
ITEM-18.3	160431	<b>Carbon electrode for multiple use – 50 x 50mm</b> • Carbon electrodes for use with TENS machines • Suitable for multiple use <b>Dimensions:</b> 50 x 50mm <b>Packaging requirements:</b> Set of 4 <b>Estimated usage:</b> 27	Yes		Per set														
ITEM 19	Protective eyewear																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-19.1	160434	<b>Protective Eyewear for Therapeutic Laser machine of ± 830 nm wavelength</b> <ul style="list-style-type: none"><li>• Sturdy ergonomic design</li><li>• Suitable for emissions of wavelength of ± 830 nm</li></ul> <b>Dimensions:</b> Suitable for most face types <b>Packaging requirements:</b> Single pair <b>Estimated usage:</b> 1	Yes		Each														
ITEM-19.2	178075	<b>Protective Eyewear for UV Phototherapy</b> <ul style="list-style-type: none"><li>• For use by patients when undergoing phototherapy</li><li>• Blocks 100% UVA / UVB</li><li>• Totally opaque to UV Rays</li></ul> <b>Dimensions:</b> Rubber eye pieces with adjustable nose bridge for comfortable fit <b>Packaging requirements:</b> Single pair <b>Estimated usage:</b> 1	Yes		Each														
ITEM-19.3	178076	<b>Protective Eyewear for UV Phototherapy</b> <ul style="list-style-type: none"><li>• For use by patients after UVA / UVB light treatment</li><li>• Blocks 100% UVA / UVB</li><li>• Totally opaque to UV rays</li></ul> <b>Dimensions:</b> UV blocking wraparound glasses <b>Packaging requirements:</b> Single pair <b>Estimated usage:</b> 1	Yes		Each														
ITEM 20	<b>Positioner: Preformed Resting Foot Drop Splint</b>																		
ITEM-20.1	New	<b>Positioner: Preformed Resting Foot Drop Splint - Small Left</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes 3-6 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 20.1 to 20.6		Each														
ITEM-20.2	New	<b>Positioner: Preformed Resting Foot Drop Splint - Small Right</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes 3-6 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 20.1 to 20.6		Each														
ITEM-20.3	New	<b>Positioner: Preformed Resting Foot Drop Splint - Medium Left</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 7-9 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 20.1 to 20.6		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-20.4	New	<b>Positioner: Preformed Resting Foot Drop Splint - Medium Right</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 7-9 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 20.1 to 20.6		Each														
ITEM-20.5	New	<b>Positioner: Preformed Resting Foot Drop Splint - Large Left</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 10-12 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 20.1 to 20.6		Each														
ITEM-20.6	New	<b>Positioner: Preformed Resting Foot Drop Splint - Large Right</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 10-12 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 20.1 to 20.6		Each														
ITEM 21		<b>Paediatric Positioner: Preformed Resting Foot Drop Splint</b>																	
ITEM-21.1	New	<b>Paediatric Positioner: Preformed Resting Foot Drop Splint - Small Left</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 4-6 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 21.1 to 21.6		Each														
ITEM-21.2	New	<b>Paediatric Positioner: Preformed Resting Foot Drop Splint - Small Right</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 4-6 <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 21.1 to 21.6		Each														



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-21.3	New	<b>Paediatric Positioner: Preformed Resting Foot Drop Splint - Medium Left</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes 7-9 <b>Packaging requirements:</b> Single item <b>Estimated usage:1</b>	Yes, representative sample for items 21.1 to 21.6		Each														
ITEM-21.4	New	<b>Paediatric Positioner: Preformed Resting Foot Drop Splint - Medium Right</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes 7-9 <b>Packaging requirements:</b> Single item <b>Estimated usage: 1</b>	Yes, representative sample for items 21.1 to 21.6		Each														
ITEM-21.5	New	<b>Paediatric Positioner: Preformed Resting Foot Drop Splint - Large Left</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 10-12 <b>Packaging requirements:</b> Single item <b>Estimated usage: 1</b>	Yes, representative sample for items 21.1 to 21.6		Each														
ITEM-21.6	New	<b>Paediatric Positioner: Preformed Resting Foot Drop Splint - Large Right</b> <ul style="list-style-type: none"><li>• Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles</li><li>• Helps protect heel area by preventing excessive pressure and lessening friction</li><li>• Velcro closure system</li><li>• Washable/ autoclavable</li></ul> <b>Dimensions:</b> Shoe sizes ± 10-12 <b>Packaging requirements:</b> Single item <b>Estimated usage: 1</b>	Yes, representative sample for items 21.1 to 21.6		Each														
ITEM 22	<b>Positioner: Foam Resting Foot Drop Splint</b>																		
ITEM 22.1	New	<b>Positioner: Foam Resting Foot drop splint - right foot</b> <ul style="list-style-type: none"><li>• Foam or foam coated with washable vinyl orthosis</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles and to reduce pressure on the heel</li><li>• Anatomically contoured right</li><li>• Side opening for cross ventilation and visual contact</li><li>• Secure fastening system</li><li>• Washable or autoclavable</li></ul> <b>Dimensions:</b> ± 40cm long from heel to calf, for right foot <b>Packaging requirements:</b> Single item <b>Estimated usage: 1</b>	Yes, representative sample for items 22.1 and 22.2		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-22.2	New	<b>Positioner: Foam Resting Foot drop splint - left foot</b> <ul style="list-style-type: none"><li>• Foam or foam coated with washable vinyl orthosis</li><li>• Not for ambulation; used on bed-ridden clients</li><li>• Used to prevent muscle and tendon shortening of the calf muscles and to reduce pressure on the heel</li><li>• Anatomically contoured left</li><li>• Side opening for cross ventilation and visual contact</li><li>• Secure fastening system</li><li>• Washable or autoclavable</li></ul> <b>Dimensions:</b> ± 40cm long from heel to calf, for left foot <b>Packaging requirements:</b> Single item <b>Estimated usage:</b> 1	Yes, representative sample for items 22.1 and 22.2		Each														
ITEM 23	<b>Hot pack moist heat hydro</b>																		
ITEM-23.1	172223	<b>Hot pack moist heat hydro, half size 13x30cm</b> Hotpack moist head hydrocollator; half size, 13cm width; 30cm length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-23.2	172221	<b>Hot pack moist heat hydro, oversize 38x61cm</b> Hotpack moist head hydrocollator; oversize 38cm width; 61cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-23.3	172236	<b>Hot pack moist heat hydro, spinal 25 x 46cm</b> Hotpack moist head hydrocollator; spinal 25cm width; 46cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-23.4	172238	<b>Hot pack moist heat hydro, spinal 25 x 61cm</b> Hotpack moist head hydrocollator; spinal 25cm width; 61cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-23.5	172226	<b>Hot pack moist heat, knee-shoulder 25x50cm</b> Hotpack moist head hydrocollator; spinal 25cm width; 50cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-23.6	136433	<b>Hot pack moist heat hydro, cervical 60cm x 15cm</b> Hotpack moist head hydrocollator; cervical 15cm width; 60cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. To fit the contours of the neck. Apply moist heat for a minimum of 30mins. Pack are re-usable. <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 0	Yes		Each														
ITEM-23.7	172220	<b>Hot pack moist heat hydro, standard 25cm x 30cm</b> Hotpack moist head hydrocollator; standard 25cm width; 30cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable. <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 24	<b>Washable hot pack covers</b>																		
ITEM-24.1	117731	<b>Washable hot pack covers, standard 48x69cm</b> Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. <b>Size:</b> standard 48x69cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-24.2	117733	<b>Washable hot pack covers, oversize 61x91cm</b> Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. <b>Size:</b> oversize 61x91cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-24.3	117734	<b>Washable hot pack covers, cervical 64x46cm</b> Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot packs covers. <b>Size:</b> cervical 64x46cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-24.4	117805	<b>Washable hot pack covers, to fit knee/shoulder hotpack 25x51cm</b> Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. <b>Size:</b> knee/shoulder 25x51cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-24.5	117806	<b>Washable hot pack covers, to fit half size 13x30cm</b> Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. <b>Size:</b> half size 13x30cm. <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 25	<b>Hemi-Arm Sling</b>																		
ITEM-25.1	136434	<b>Hemi-Arm Sling for Paeds - Small</b> Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. <b>Size:</b> small <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 25.1 to 25.3		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-25.2	136435	<b>Hemi-Arm Sling for Paeds - Medium</b> Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. <b>Size:</b> medium <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 25.1 to 25.3		Each														
ITEM-25.3	136436	<b>Hemi-Arm Sling for Paeds - Large</b> Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. <b>Size:</b> large <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 25.1 to 25.3		Each														
ITEM-25.4	117813	<b>Hemi-Arm Sling for Adult - Small</b> Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. <b>Size:</b> small <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 25.4 to 25.6		Each														
ITEM-25.5	117814	<b>Hemi-Arm Sling for Adult - Medium</b> Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. <b>Size:</b> medium <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 25.4 to 25.6		Each														
ITEM-25.6	117815	<b>Hemi-Arm Sling for Adult - Large</b> Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. <b>Size:</b> large <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 25.4 to 25.6		Each														
ITEM 26	Probes																		
ITEM-26.1	136437	<b>Intra vaginal probe</b> Can be used with most muscle stimulators and biofeedback equipment. Single patient use. Medical grade surgical stainless steel. Unique rectangular section resist lateral movement. Countered electrode provides comfortable stimulation and accurate biofeedback. Smooth clean shape. Long integrated cables. Push fit indicator wand with extension pieces. <b>Dimensions:</b> length 77mm, conductive part width 14mm with a 20mm width of the probe <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-26.2	136438	<b>Intra-anal and small intra vaginal probe</b> When inserted, the probe allows the anal canal sphincter to close around the narrow neck of the probe. Revolutionary shape for repeatable positioning and ambulatory use. Can be used with most muscle stimulators and biofeedback equipment. Single patient use. Medical grade surgical stainless steel. Unique rectangular section resist lateral movement. Countered electrode provides comfortable stimulation. Smooth clean shape. Long integrated cables. Always sits in the same position comparable EMR readings can be taken over different sessions. Allows for comfortable muscle stimulation and accurate biofeedback. <b>Dimensions:</b> length 83mm, conductive part diameter 10mm with a 14mm diameter of the probe <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 27	<b>Pelvic floor exercise system</b>																		
ITEM-27.1	136439	<b>Pelvic floor exercise system</b> Exercise system consists of two cones and a set of 4 weight or more sizes (5g,10,2x20g). Single patient use. <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
<b>B:- OCCUPATIONAL THERAPY</b>																			
ITEM 28	<b>Thermoplastic splinting material with cohesive &amp; adhesive properties</b>																		
ITEM-28.1	156745	<b>PERFORATED thermoplastic splinting material with cohesive &amp; adhesive properties -3.2mm</b> • Low temperature thermoplastic material used to make splints • Micro-perforated 13% • Has the following properties: - Resistance to stretch: Moderate - Drape/Conformity: Maximum - Memory/Recovery: 100% - Rigidity: minimum - Bonding: cohesive & adhesive properties - Surface finish: Fingerprint resist • Colours: Standard (fleshy pink) <b>Dimensions:</b> 46cm x 61cm x 3.2mm sheet <b>Packaging requirements:</b> Pack of 4 sheets <b>Estimated usage:</b> 1	Yes		Pack of 4 sheet														
ITEM-28.2	178107	<b>NON PERFORATED thermoplastic splinting material with cohesive &amp; adhesive properties - 1,6mm</b> • Low temperature thermoplastic material used to make splints • Has the following properties: - Resistance to stretch: Moderate - Drape/Conformity: Maximum - Memory/Recovery: 100% - Rigidity: minimum - Bonding: cohesive & adhesive properties - Surface finish: Fingerprint resist • Colours: Standard (fleshy pink) <b>Dimensions:</b> 46cm x 61cm x 1.6mm sheet <b>Packaging requirements:</b> Pack of 4 sheets <b>Estimated usage:</b> 1	Yes		Pack of 4 sheet														
ITEM 29	<b>Thermoplastic splinting material with temporary bonding properties</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-29.1	178109	<b>NON-PERFORATED thermoplastic splinting material with temporary bonding properties - 2,4 mm</b> <ul style="list-style-type: none"><li>• Low temperature thermoplastic material used to make splints</li><li>• Non-perforated</li><li>• Has the following properties:<ul style="list-style-type: none"><li>- Resistance to stretch: Maximum</li><li>- Drape/Conformity: Moderate</li><li>- Memory/Recovery: 100%</li><li>- Rigidity: Maximum</li><li>- Bonding: Temporary bonding qualities only</li><li>- Surface finish: Fingerprint resist; transparent when heated; non-stick coating</li></ul></li><li>• Colour: White</li></ul> <b>Dimensions:</b> 46cm x 61 cm x 2.4mm sheet <b>Packaging requirements:</b> Pack of 4 sheets <b>Estimated usage:</b> 24	Yes		Sheet														
ITEM-29.2	156750	<b>PERFORATED thermoplastic splinting material with temporary bonding properties - 2,4mm</b> <ul style="list-style-type: none"><li>• Low temperature thermoplastic material used to make splints</li><li>• 1% perforated</li><li>• Has the following properties:<ul style="list-style-type: none"><li>- Resistance to stretch: Maximum</li><li>- Drape/Conformity: Moderate</li><li>- Memory/Recovery: 100%</li><li>- Rigidity: Maximum</li><li>- Bonding: Temporary bonding qualities only</li><li>- Surface finish: Fingerprint resist, transparent when heated, non-stick coating</li></ul></li><li>• Colour: White</li></ul> <b>Dimensions:</b> 46cm x 61 cm x 2.4 mm sheet <b>Packaging requirements:</b> Pack of 4 sheets <b>Estimated usage:</b> 1	Yes		Sheet														
ITEM 30	<b>Non-Perforated thermo-plastic splinting material -Equivalent to OMEGA MAX</b>																		
ITEM-30.1	New	<b>Non-Perforated thermo-plastic splinting material - 3.2mm</b> <ul style="list-style-type: none"><li>• Low temperature thermoplastic material used to make splints</li><li>• Has the following properties:<ul style="list-style-type: none"><li>- Resistance to stretch: Moderate</li><li>- Drape/Conformity: Moderate to maximum</li><li>- Memory/Recovery: 100%</li><li>- Rigidity: Moderate to maximum</li><li>- Bonding: Has a powder coating to give a slightly tacky surface that holds material in place. - When cut gives smooth self sealed edges.</li><li>- Working time of 3 minutes to form and harden after softening.</li><li>- Similar or equivalent to OMEGA MAX material.</li><li>- Surface finish: Fingerprint resist</li></ul></li><li>• Colour: Beige</li></ul> <b>Dimensions:</b> 45cm x 60cm x 3.2mm sheet <b>Packaging requirements:</b> Pack of 4 sheets <b>Estimated usage:</b> 1	Yes		each														



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-31.4	New	<b>Paediatric Wrist Brace with closure system - Small, Left</b> <ul style="list-style-type: none"><li>Breathable; durable elastic material</li><li>Latex sealed edging</li><li>Malleable aluminium strip covered with leather or vinyl</li><li>Removable volar metal stay</li><li>Palmar bar should not extend past palmar crease</li><li>Length of gutter should be a minimum of half of forearm length as appropriate for size of splint . Narrow strap through 1st webspace</li></ul> <b>Dimensions:</b> wrist circumference of 13-14cm, left hand <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 10	Yes, representative sample for items 31.1 to 31.8		Each														
ITEM-31.5	New	<b>Paediatric Wrist Brace with closure system - Medium, Left</b> <ul style="list-style-type: none"><li>Breathable; durable elastic material</li><li>Latex sealed edging</li><li>Malleable aluminium strip covered with leather or vinyl</li><li>Removable volar metal stay</li><li>Palmar bar should not extend past palmar crease</li><li>Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace.</li></ul> <b>Dimensions:</b> wrist circumference of 14-15cm, left hand <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 10	Yes, representative sample for items 31.1 to 31.8		Each														
ITEM-31.6	New	<b>Paediatric Wrist Brace with closure system - Medium, Right</b> <ul style="list-style-type: none"><li>Breathable; durable elastic material</li><li>Latex sealed edging</li><li>Malleable aluminium strip covered with leather or vinyl</li><li>Removable volar metal stay</li><li>Palmar bar should not extend past palmar crease</li><li>Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace.</li></ul> <b>Dimensions:</b> wrist circumference of 14-15cm, right hand <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 10	Yes, representative sample for items 31.1 to 31.8		Each														
ITEM-31.7	New	<b>Paediatric Wrist Brace with closure system - Large, Right</b> <ul style="list-style-type: none"><li>Breathable; durable elastic material</li><li>Latex sealed edging</li><li>Malleable aluminium strip covered with leather or vinyl</li><li>Removable volar metal stay</li><li>Palmar bar should not extend past palmar crease</li><li>Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace.</li></ul> <b>Dimensions:</b> wrist circumference of 15-17cm, right hand <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 10	Yes, representative sample for items 31.1 to 31.8		Each														
ITEM-31.8	New	<b>Paediatric Wrist Brace with closure system - Large, Left</b> <ul style="list-style-type: none"><li>Breathable; durable elastic material</li><li>Latex sealed edging</li><li>Malleable aluminium strip covered with leather or vinyl</li><li>Removable volar metal stay</li><li>Palmar bar should not extend past palmar crease</li><li>Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace.</li></ul> <b>Dimensions:</b> wrist circumference of 15-75cm, left hand <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 10	Yes, representative sample for items 31.1 to 31.8		Each														
ITEM 32	Wrist and Thumb Splint																		



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-32.1	New	<b>Wrist and Thumb Splint - Extra Small Right</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit right hand circumference of 13cm to 14cm <b>Packaging requirements:</b> Per item <b>Estimated usage:1</b>	Yes		Each														
ITEM-32.2	New	<b>Wrist and Thumb Splint - Extra Small Left</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit left hand circumference of 13cm to 14cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes		Each														
ITEM-32.3	New	<b>Wrist and Thumb Splint - Small Right</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit right hand circumference of 14cm to 16cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes		Each														
ITEM-32.4	New	<b>Wrist and Thumb Splint - Small Left</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit left hand circumference of 14cm to 16cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-32.5	New	<b>Wrist and Thumb Splint - Medium Right</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint</li></ul> . Must include thumb support <b>Dimensions:</b> To fit right hand circumference of 16cm to 19cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes, representative sample for items 32.1 to 32.10		Each														
ITEM-32.6	New	<b>Wrist and Thumb Splint - Medium Left</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit left hand circumference of 16cm to 19cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes, representative sample for items 32.1 to 32.10		Each														
ITEM-32.7	New	<b>Wrist and Thumb Splint - Large Right</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit right hand circumference of 19cm to 20cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes, representative sample for items 32.1 to 32.10		Each														
ITEM-32.8	New	<b>Wrist and Thumb Splint - Large Left</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit left hand circumference of 19cm to 20cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes, representative sample for items 32.1 to 32.10		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-32.9	New	<b>Wrist and Thumb Splint - Extra Large Right</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support</li></ul> <b>Dimensions:</b> To fit right hand circumference of 20cm to 22cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes, representative sample for items 32.1 to 32.10		Each														
ITEM-32.10	New	<b>Wrist and Thumb Splint - Extra Large Left</b> <ul style="list-style-type: none"><li>• Used to prevent thumb MP joint motion,</li><li>• Made of breathable material with D-Ring or hook and loop straps</li><li>• Malleable aluminium strip covered with leather or vinyl</li><li>• Removable volar metal stay</li><li>• Palmar bar should not extend past palmar crease</li><li>• Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support.</li></ul> <b>Dimensions:</b> To fit left hand circumference of 20cm to 22cm <b>Packaging requirements:</b> Per item <b>Estimated usage: 1</b>	Yes, representative sample for items 32.1 to 32.10		Each														
ITEM 33	<b>Pressure garment material</b>																		
ITEM-33.1	156777	<b>Pressure garment material - 50m, beige</b> <ul style="list-style-type: none"><li>• Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation</li><li>• Highly durable</li><li>• Is resistant to stretch - has 100% memory</li><li>• Has extended shelf life</li><li>• Colours: Beige</li></ul> <b>Dimensions:</b> Width of material: 180cm <b>Packaging requirements:</b> 50 metre roll <b>Estimated usage: 13</b>	Yes, representative sample for items 33.1 to 33.6		Per roll														
ITEM-33.2	New	<b>Pressure garment material - 25m, beige</b> <ul style="list-style-type: none"><li>• Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation</li><li>• Highly durable</li><li>• Is resistant to stretch - has 100% memory</li><li>• Has extended shelf life</li><li>• Colours: Beige</li></ul> <b>Dimensions:</b> Width of material: 180cm <b>Packaging requirements:</b> 25 metre roll <b>Estimated usage: 1</b>	Yes, representative sample for items 33.1 to 33.6		Per roll														
ITEM-33.3	New	<b>Pressure garment material - 10m, beige</b> <ul style="list-style-type: none"><li>• Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation</li><li>• Highly durable</li><li>• Is resistant to stretch - has 100% memory</li><li>• Has extended shelf life</li><li>• Colours: Beige</li></ul> <b>Dimensions:</b> Width of material: 180cm <b>Packaging requirements:</b> 10 metre roll <b>Estimated usage: 1</b>	Yes, representative sample for items 33.1 to 33.6		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-33.4	156776	<b>Pressure garment material - 50m, brown</b> <ul style="list-style-type: none"><li>• Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation</li><li>• Highly durable</li><li>• Is resistant to stretch - has 100% memory</li><li>• Has extended shelf life</li><li>• Colours: Brown</li></ul> <b>Dimensions:</b> Width of material: 180cm <b>Packaging requirements:</b> 50 metre roll <b>Estimated usage:</b> 13	Yes, representative sample for items 33.1 to 33.6		Per roll														
ITEM-33.5	New	<b>Pressure garment material - 25m, brown</b> <ul style="list-style-type: none"><li>• Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation</li><li>• Highly durable</li><li>• Is resistant to stretch - has 100% memory</li><li>• Has extended shelf life</li><li>• Colours: Brown</li></ul> <b>Dimensions:</b> Width of material: 180cm <b>Packaging requirements:</b> 25 metre roll <b>Estimated usage:</b> 1	Yes, representative sample for items 33.1 to 33.6		Per roll														
ITEM-33.6	New	<b>Pressure garment material - 10m, brown</b> <ul style="list-style-type: none"><li>• Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation</li><li>• Highly durable</li><li>• Is resistant to stretch - has 100% memory</li><li>• Has extended shelf life</li><li>• Colours: Brown</li></ul> <b>Dimensions:</b> Width of material: 180cm <b>Packaging requirements:</b> 10 metre roll <b>Estimated usage:</b> 1	Yes, representative sample for items 33.1 to 33.6		Per roll														
ITEM 34	Zipper																		
ITEM-34.1	New	<b>Zipper Beige 30cm long</b> <ul style="list-style-type: none"><li>• Can be used to fabricate pressure garments</li><li>• Highly durable</li><li>• Plastic zipper</li><li>• Must have top and bottom stop</li><li>• Colours: Beige</li></ul> <b>Dimensions:</b> • 30cm long <b>Packaging requirements:</b> • 10 in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 34.1 to 34.6		Each														
ITEM-34.2	New	<b>Zipper Beige 20cm Long</b> <ul style="list-style-type: none"><li>• Can be used to fabricate pressure garments</li><li>• Highly durable</li><li>• Plastic zipper</li><li>• Must have top and bottom stop</li><li>• Colours: Beige</li></ul> <b>Dimensions:</b> 20cm long <b>Packaging requirements:</b> 10 in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 34.1 to 34.6		Each														
ITEM-34.3	New	<b>Zipper Beige 15cm Long</b> <ul style="list-style-type: none"><li>• Can be used to fabricate pressure garments</li><li>• Highly durable</li><li>• Plastic zipper</li><li>• Must have top and bottom stop</li><li>• Colours: Beige</li></ul> <b>Dimensions:</b> 15cm long <b>Packaging requirements:</b> 10 in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 34.1 to 34.6		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-34.4	New	<b>Zipper Brown 30cm Long</b> <ul style="list-style-type: none"><li>• Can be used to fabricate pressure garments</li><li>• Highly durable</li><li>• Plastic zipper</li><li>• Must have top and bottom stop</li><li>• Colours: Brown</li></ul> <b>Dimensions:</b> 30cm long <b>Packaging requirements:</b> 10 in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 34.1 to 34.6		Each														
ITEM-34.5	New	<b>Zipper Brown 20cm Long</b> <ul style="list-style-type: none"><li>• Can be used to fabricate pressure garments</li><li>• Highly durable</li><li>• Plastic zipper</li><li>• Must have top and bottom stop</li><li>• Colours: Brown</li></ul> <b>Dimensions:</b> 20cm long <b>Packaging requirements:</b> 10 in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 34.1 to 34.6		Each														
ITEM-34.6	New	<b>Zipper Brown 15cm Long</b> <ul style="list-style-type: none"><li>• Can be used to fabricate pressure garments</li><li>• Highly durable</li><li>• Plastic zipper</li><li>• Must have top and bottom stop</li><li>• Colours: Brown</li></ul> <b>Dimensions:</b> 15cm long <b>Packaging requirements:</b> 10 in a box <b>Estimated usage:</b> 1	Yes, representative sample for items 34.1 to 34.6		Each														
ITEM 35	Neoprene																		
ITEM-35.1	156779	<b>Neoprene 3.2mm thick Non-Perforated</b> <ul style="list-style-type: none"><li>• Nylon lining on one side</li><li>• Pile material on the other over a rubber core,</li><li>• Has a degree of elasticity resulting in comfortable and secure strapping.</li><li>• Adhesive to male velcro</li><li>• Latex free</li></ul> <b>Dimensions:</b> 46cm x 61cm x 3,2mm sheet <b>Packaging requirements:</b> Single sheet <b>Estimated usage:</b> 29	Yes, representative sample for items 35.1 and 35.2		Per sheet														
ITEM-35.2	156781	<b>Neoprene 1.5mm thick Non-Perforated</b> <ul style="list-style-type: none"><li>• Nylon lining on one side</li><li>• Pile material on the other over a rubber core,</li><li>• Has a degree of elasticity resulting in comfortable and secure strapping.</li><li>• Adhesive to male velcro</li><li>• Latex free</li></ul> <b>Dimensions:</b> 46cm x 61 cm x 1.5mm sheet <b>Packaging requirements:</b> Single sheet <b>Estimated usage:</b> 61	Yes, representative sample for items 35.1 and 35.2		Per sheet														
ITEM-35.3	New	<b>Neoprene 3.2mm thick Perforated</b> <ul style="list-style-type: none"><li>• Nylon lining on one side</li><li>• Pile material on the other over a rubber core,</li><li>• Has a degree of elasticity resulting in comfortable and secure strapping.</li><li>• Adhesive to male velcro</li><li>• Latex free</li></ul> <b>Dimensions:</b> 46cm x 61cm x 3,2mm sheet <b>Packaging requirements:</b> Single sheet <b>Estimated usage:</b> 29	Yes, representative sample for items 35.3 and 35.4		Per sheet														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-35.4	New	<b>Neoprene 1.5mm thick Perforated</b> <ul style="list-style-type: none"><li>• Nylon lining on one side</li><li>• Pile material on the other over a rubber core,</li><li>• Has a degree of elasticity resulting in comfortable and secure strapping.</li><li>. Adhesive to male velcro</li><li>• Latex free</li></ul> <b>Dimensions:</b> 46cm x 61 cm x 1.5mm sheet <b>Packaging requirements:</b> Single sheet <b>Estimated usage:</b> 61	Yes, representative sample for items 35.3 and 35.4		Per sheet														
ITEM 36	<b>Neoprene like wrap on a ROLL</b>																		
ITEM-36.1	156782	<b>Neoprene like wrap on a ROLL; 7.6cm wide</b> <ul style="list-style-type: none"><li>• Piled fabric on one side</li><li>• Woven blend of nylon &amp; lycra outer layer</li><li>• Elastomeric foam inner layer.</li><li>• Medium conformability</li><li>• Min to Moderate elasticity</li><li>. Adhesive to male velcro</li><li>• Mostly linear stretch</li><li>• Non fraying</li></ul> <b>Dimensions:</b> 7.6cm x 4.5m roll x ± 2mm thick <b>Packaging requirements:</b> Per roll <b>Estimated usage:</b> 1	Yes, representative sample for items 36.1 to 36.4		Per roll														
ITEM-36.2	156787	<b>Neoprene like wrap on a ROLL; 15cm wide</b> <ul style="list-style-type: none"><li>• Piled fabric on one side</li><li>• Woven blend of nylon &amp; lycra outer layer</li><li>• Elastomeric foam inner layer.</li><li>• Medium conformability</li><li>• Min to Moderate elasticity</li><li>. Adhesive to male velcro</li><li>• Mostly linear stretch</li><li>• Non fraying</li></ul> <b>Dimensions:</b> 15cm x 4.5m roll x ± 2mm thick <b>Packaging requirements:</b> Per roll <b>Estimated usage:</b> 1	Yes, representative sample for items 36.1 to 36.4		Per roll														
ITEM-36.3	156788	<b>Neoprene like wrap on a ROLL; 6.4cm wide</b> <ul style="list-style-type: none"><li>• Piled fabric on one side</li><li>• Woven blend of nylon &amp; lycra outer layer</li><li>• Elastomeric foam inner layer.</li><li>• Moderate elasticity</li><li>• Multi-directional stretch</li><li>. Adhesive to male velcro</li><li>• Non fraying</li><li>• Medium conformability</li></ul> <b>Dimensions:</b> 6.4cm x 4.5m roll x 1.5mm <b>Packaging requirements:</b> Per roll <b>Estimated usage:</b> 2	Yes, representative sample for items 36.1 to 36.4		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-36.4	156789	<b>Neoprene like wrap on a ROLL; 15cm wide</b> <ul style="list-style-type: none"><li>• Piled fabric on one side</li><li>• Woven blend of nylon &amp; lycra outer layer</li><li>• Elastomeric foam inner layer.</li><li>• Moderate elasticity</li><li>• Multi-directional stretch</li><li>• Adhesive to male velcro</li><li>• Non fraying</li><li>• Medium conformability</li></ul> <b>Dimensions:</b> 15cm x 4.5m roll x 1.5mm <b>Packaging requirements:</b> Per roll <b>Estimated usage:</b> 37	Yes, representative sample for items 36.1 to 36.4		Per roll														
ITEM-36.5	156415	<b>Neoprene-like-wrap in a SHEET form - Black</b> <ul style="list-style-type: none"><li>• Piled fabric on one side</li><li>• Woven blend of nylon &amp; lycra outer layer</li><li>• Elastomeric foam inner layer.</li><li>• High conformability</li><li>• Moderate elasticity</li><li>• Multi-directional stretch</li><li>• Adhesive to male velcro</li><li>• Non fraying</li><li>• Colour: Black</li></ul> <b>Dimensions:</b> ±43.5cm x ±56.4cm x ±1.5mm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 1	Yes, representative sample for items 36.5 and 36.6		Sheet														
ITEM-36.6	156429	<b>Neoprene-like-wrap in a SHEET form - Dark blue</b> <ul style="list-style-type: none"><li>• Piled fabric on one side</li><li>• Woven blend of nylon &amp; lycra outer layer</li><li>• Elastomeric foam inner layer.</li><li>• High conformability</li><li>• Moderate elasticity</li><li>• Multi-directional stretch</li><li>• Adhesive to male velcro</li><li>• Non fraying</li><li>• Colour: Dark blue</li></ul> <b>Dimensions:</b> ±43.5cm x ±56.4cm ±1.5mm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 2	Yes, representative sample for items 36.5 and 36.6		Sheet														
ITEM 37	Padding																		
ITEM-37.1	New	<b>Pressure relief padding</b> <ul style="list-style-type: none"><li>• Insulates against shock and vibration</li><li>• Use to cushion and protects hand splints</li><li>• Conforms to all body contours</li><li>• Effectively eliminating abrasion and shearing</li><li>• Latex free, cleanable, 2mm thick gel pad with self adhesive backing</li></ul> <b>Dimensions:</b> 20cm x 91cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 4	Yes		per sheet														
ITEM-37.2	156443	<b>Padding material for thermoplastic splinting – adhesive, towelling</b> <ul style="list-style-type: none"><li>• Self adhesive</li><li>• Towelling material</li><li>• Non-stretch</li><li>• Washable</li></ul> <b>Dimensions:</b> 50cm x 100cm x 3mm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 25	Yes		Sheet														
ITEM 38	Leather suede (for making wheelchair gloves)																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-38.1	174261	<b>Leather suede (for making wheelchair gloves)</b> • suede leather . Colour: Black <b>Dimensions:</b> ± 2 mm thick <b>Packaging requirements:</b> Per square metre <b>Estimated usage:</b> 1	Yes, representative sample for items 38.1 and 38.2		M²														
ITEM-38.2	178113	<b>Leather suede (for making wheelchair gloves)</b> • suede leather . Colour: Tan <b>Dimensions:</b> ± 2 mm thick <b>Packaging requirements:</b> Per square metre <b>Estimated usage:</b> 1	Yes, representative sample for items 38.1 and 38.2		M²														
ITEM 39	<b>Towelling</b>																		
ITEM-39.1	173609	<b>Towelling</b> • Piled on both sides of material • White <b>Dimensions:</b> Minimum 150cm wide <b>Packaging requirements:</b> 10 metre <b>Estimated usage:</b> 1	Yes, representative sample for items 39.1 and 39.2		Per meter														
ITEM-39.2	156497	<b>Towelling</b> • Piled on both sides of material • Dark colours <b>Dimensions:</b> Minimum 150cm wide <b>Packaging requirements:</b> 10 metre <b>Estimated usage:</b> 1	Yes, representative sample for items 39.1 and 39.2		Per meter														
ITEM 40	<b>Calico</b>																		
ITEM-40.1	156498	<b>Calico White</b> • Natural cotton • Pre-shrunk . Colour: White <b>Dimensions:</b> Minimum 115 cm wide <b>Packaging requirements:</b> 10 metre <b>Estimated usage:</b> 1	Yes, representative sample for items 40.1 and 40.2		Per meter														
ITEM-40.2	173610	<b>Calico Cream</b> • Natural cotton • Pre-shrunk . Colour: Cream <b>Dimensions:</b> Minimum 115 cm wide <b>Packaging requirements:</b> 10 metre <b>Estimated usage:</b> 1	Yes, representative sample for items 40.1 and 40.2		Per meter														
ITEM 41	<b>Foam</b>																		
ITEM-41.1	156499	<b>Foam, medium density sheet, 5cm thick</b> • Open cell foam • Medium density Non flammable <b>Dimensions:</b> Minimum 180cm x 91cm x 5cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 1	Yes, representative sample for items 41.1 and 41.2		Sheet														
ITEM-41.2	161592	<b>Foam, medium density sheet, 10cm thick</b> • Open cell foam • Medium density Non flammable <b>Dimensions:</b> Minimum 180cm x 91cm x 10cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 1	Yes, representative sample for items 41.1 and 41.2		Sheet														



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-41.3	161593	<b>Foam, high density sheets, 5cm thick</b> • Open cell foam • High density Non flammable <b>Dimensions:</b> Minimum 180cm x 91cm x 5cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage: 1</b>	Yes, representative sample for items 41.3 and 41.4		Sheet														
ITEM-41.4	156500	<b>Foam, high density sheets, 10cm thick</b> • Open cell foam • High density Non flammable <b>Dimensions:</b> Minimum 180cm x 91cm x 10cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage: 1</b>	Yes, representative sample for items 41.3 and 41.4		Sheet														
ITEM 42	<b>Memory Foam</b>																		
ITEM-42.1	New	<b>Memory Foam, medium density sheets, 2cm thick</b> • 16-80kg/m3 , Visco elastic polymer • Medium density Non flammable <b>Dimensions:</b> Minimum 100cm x 180cm x 2cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage: 1</b>	Yes, representative sample for items 42.1 and 42.2		Per sheet														
ITEM-42.2	New	<b>Memory Foam, medium density sheets, 5cm thick</b> • 16-80kg/m3, Visco elastic polymer • Medium density Non flammable <b>Dimensions:</b> Minimum 100cm x 180cm x 5cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage: 1</b>	Yes, representative sample for items 42.1 and 42.2		Per sheet														
ITEM 43	<b>Wood, SA pine</b>																		
ITEM-43.1	156506	<b>Wood, SA Pine, 22mm thickness</b> • South African grown pine • Solid • Planed • To be cut in meter square blocks <b>Dimensions:</b> •22 mm thick <b>Packaging requirements:</b> • Per plank; cut as specified <b>Estimated usage: 1</b>	Yes		M²														
ITEM 44	<b>Wood, Superwood</b>																		
ITEM-44.1	156513	<b>Wood, Superwood, 3mm thick</b> • Smooth on both sides • Durable • To be cut in meter square blocks <b>Dimensions:</b> .Thickness = 3mm <b>Packaging requirements:</b> Per plank; cut as specified <b>Estimated usage: 1</b>	Yes, representative sample for items 44.1 to 44.7		M²														
ITEM-44.2	156509	<b>Wood, Superwood, 6mm thick</b> • Smooth on both sides • Durable • To be cut in meter square blocks <b>Packaging requirements:</b> Per plank; cut as specified <b>Estimated usage: 1</b>	Yes, representative sample for items 44.1 to 44.7		M²														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-44.3	156510	<b>Wood, Superwood, 9mm thick</b> <ul style="list-style-type: none"><li>• Smooth on both sides</li><li>• Durable</li><li>• To be cut in meter square blocks</li></ul> <b>Dimensions:</b> .Thickness = 9mm <b>Packaging requirements:</b> Per plank; cut as specified <b>Estimated usage: 1</b>	Yes, representative sample for items 44.1 to 44.7		M²														
ITEM-44.4	156511	<b>Wood, Superwood, 12mm thick</b> <ul style="list-style-type: none"><li>• Smooth on both sides</li><li>• Durable</li><li>• To be cut in meter square blocks</li></ul> <b>Dimensions:</b> .Thickness = 12mm <b>Packaging requirements:</b> Per plank; cut as specified <b>Estimated usage:1</b>	Yes, representative sample for items 44.1 to 44.7		M²														
ITEM-44.5	156515	<b>Wood, Superwood, 18mm thick</b> <ul style="list-style-type: none"><li>• Smooth on both sides</li><li>• Durable</li><li>• To be cut in meter square blocks</li></ul> <b>Dimensions:</b> .Thickness =18mm <b>Packaging requirements:</b> Per plank; cut as specified <b>Estimated usage: 1</b>	Yes, representative sample for items 44.1 to 44.7		M²														
ITEM-44.6	156517	<b>Wood, Superwood, 22mm thick</b> <ul style="list-style-type: none"><li>• Smooth on both sides</li><li>• Durable</li><li>• To be cut in meter square blocks</li></ul> <b>Dimensions:</b> .Thickness = 22mm <b>Packaging requirements:</b> Per plank; cut as specified <b>Estimated usage: 1</b>	Yes, representative sample for items 44.1 to 44.7		M²														
ITEM-44.7	156522	<b>Wood, Superwood, 32mm thick</b> <ul style="list-style-type: none"><li>• Smooth on both sides</li><li>• Durable</li><li>• To be cut in meter square blocks</li></ul> <b>Dimensions:</b> .Thickness = 32mm <b>Packaging requirements:</b> Per plank; cut as specified <b>Estimated usage: 1</b>	Yes, representative sample for items 44.1 to 44.7		M²														
ITEM 45	<b>Stretch loop, non-adhesive (equivalent to velcro)</b>																		
ITEM-45.1	156601	<b>Stretch loop, 2.5cm, non-adhesive (equivalent to velcro)</b> <ul style="list-style-type: none"><li>• Stretch loop</li><li>• Elastic loop provides even tension</li><li>• Made of rubber with soft edges</li><li>• Moderate resistance to stretch for securing splints</li><li>• Snug fit and secure closure</li><li>• Mates with all Velcro type hook strapping</li><li>• Colours: white</li></ul> <b>Dimensions:</b> 25m x 2.5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 3</b>	Yes, representative sample for items 45.1 to 45.1		Per roll														
ITEM 46	<b>Soft loop, non-adhesive (equivalent to velcro)</b>																		
ITEM-46.1	174266	<b>Soft loop, 2.5cm. non-adhesive (equivalent to velcro)</b> <ul style="list-style-type: none"><li>• Loop – 2.5cm wide</li><li>• Soft</li><li>• Non-abrasive loop Velcro</li><li>• Highly durable</li><li>• Maximum bonding capability</li><li>• Colours: Dark Brown</li></ul> <b>Dimensions:</b> 25m x 2,5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 1</b>	Yes, representative sample for items 46.1 to 46.2		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-46.2	174269	<b>Soft loop, 5cm, non-adhesive (equivalent to velcro)</b> Loop – 5cm wide • Soft • Non-abrasive loop Velcro • Highly durable • Maximum bonding capability • Colours: Dark Brown <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> 1	Yes, representative sample for items 46.1 to 46.2		Per roll														
ITEM 47	<b>Loop, non-adhesive (equivalent to velcro)</b>																		
ITEM-47.1	164587	<b>Loop, non-adhesive, 2.5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> White <b>Dimensions:</b> 50m x 2,5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.2	164588	<b>Loop, non-adhesive, 5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> White <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.3	118877	<b>Loop, non-adhesive, 2cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> White <b>Dimensions:</b> 50m x 2cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.4	118265	<b>Loop, non-adhesive, 2.5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> Beige <b>Dimensions:</b> 50m x 2,5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.5	118264	<b>Loop, non-adhesive, 5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> Beige <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.6	173614	<b>Loop, non-adhesive, 2cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> Brown <b>Dimensions:</b> 50m x 2cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.7	171336	<b>Loop, non-adhesive, 2.5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> Brown <b>Dimensions:</b> 50m x 2,5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-47.8	171326	<b>Loop, non-adhesive, 5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> Brown <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.9	164567	<b>Loop, non-adhesive, 2.5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> Black <b>Dimensions:</b> 50m x 2,5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM-47.10	164569	<b>Loop, non-adhesive, 5cm (equivalent to velcro)</b> Loop, non-adhesive, nylon <b>Colours:</b> Black <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 90 rolls over range of sizes	Yes, representative sample for items 47.1 to 47.10		Per roll														
ITEM 48	<b>Hook, non-adhesive (equivalent to velcro)</b>																		
ITEM-48.1	156528	<b>Hook, non-adhesive, 2.5cm (equivalent to velcro)</b> Hook, non-adhesive, nylon Colours: Beige <b>Dimensions:</b> 25m x 2.5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 48.1 to 48.7		Per roll														
ITEM-48.2	156530	<b>Hook, non-adhesive, 5cm (equivalent to velcro)</b> Hook, non-adhesive, nylon Colours: Beige <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 48.1 to 48.7		Per roll														
ITEM-48.3	171299	<b>Hook, non-adhesive, 2.5cm (equivalent to velcro)</b> Hook, non-adhesive, nylon Colours: Brown <b>Dimensions:</b> 50m x 2.5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 48.1 to 48.7		Per roll														
ITEM-48.4	171311	<b>Hook, non-adhesive, 5cm (equivalent to velcro)</b> Hook, non-adhesive, nylon Colours: Brown <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 48.1 to 48.7		Per roll														
ITEM-48.5	173612	<b>Hook, non-adhesive, 2.5cm (equivalent to velcro)</b> Hook, non-adhesive, nylon Colours: Dark Brown <b>Dimensions:</b> 25m x 2.5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 48.1 to 48.7		Per roll														
ITEM-48.6	173613	<b>Hook, non-adhesive, 5cm (equivalent to velcro)</b> Hook, non-adhesive, nylon Colours: Dark Brown <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 48.1 to 48.7		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-48.7	164560	<b>Hook, non-adhesive, 5cm (equivalent to velcro)</b> Hook non adhesive nylon Colours: Black <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 48.1 to 48.7		Per roll														
ITEM 49	<b>Self adhesive loop (equivalent to velcro)</b>																		
ITEM-49.1	173615	<b>Self adhesive loop 2.5cm (equivalent to velcro)</b> Loop – 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: white • Highly durable <b>Dimensions:</b> 25m x 2.5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> 8	Yes, representative sample for items 49.1 to 49.8		Per roll														
ITEM-49.2	171273	<b>Self adhesive loop 5cm (equivalent to velcro)</b> Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: white • Highly durable <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 45 rolls over range of sizes	Yes, representative sample for items 49.1 to 49.8		Per roll														
ITEM-49.3	118878	<b>Self adhesive loop 2.5cm (equivalent to velcro)</b> Loop - 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: beige • Highly durable <b>Dimensions:</b> 50m x 2.5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 45 rolls over range of sizes	Yes, representative sample for items 49.1 to 49.8		Per roll														
ITEM-49.4	118879	<b>Self adhesive loop 5cm (equivalent to velcro)</b> Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: beige • Highly durable <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 45 rolls over range of sizes	Yes, representative sample for items 49.1 to 49.8		Per roll														
ITEM-49.5	171336	<b>Self adhesive loop 2.5cm (equivalent to velcro)</b> Loop - 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: brown • Highly durable <b>Dimensions:</b> 50m x 2.5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 45 rolls over range of sizes	Yes, representative sample for items 49.1 to 49.8		Per roll														
ITEM-49.6	171279	<b>Self adhesive loop 5 cm (equivalent to velcro)</b> Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: brown • Highly durable <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 45 rolls over range of sizes	Yes, representative sample for items 49.1 to 49.8		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORT ED	DIRECT MATER IAL: % LOCAL CONTE NT	%DIRE CT LABOU R	% TRAN SPOR T	%OTHE R OVER HEAD S AND %PRO FITS	%TOT AL	RATE OF EXCH ANGE	ADJ UST ME NT FRE QU EN CY	IMPO RTE D	LOCALL Y MANUF ACTUR ED
ITEM-49.7	171284	<b>Self adhesive loop 2.5 cm (equivalent to velcro)</b> Loop - 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: black • Highly durable <b>Dimensions:</b> 50m x 2.5cm <b>Packaging requirements:</b> Per roll <b>Estimated requirements:</b> 45 rolls over range of sizes	Yes, representative sample for items 49.1 to 49.8		Per roll														
ITEM-49.8	171289	<b>Self adhesive loop 5 cm (equivalent to velcro)</b> Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: black • Highly durable <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 45 rolls over range of sizes	Yes, representative sample for items 49.1 to 49.8		Per roll														
ITEM 50	<b>Hook self adhesive (equivalent to velcro)</b>																		
ITEM-50.1	156531	<b>Hook self adhesive 2cm (equivalent to velcro)</b> Hook – 2cm self adhesive • Maximum bonding capability • Colours: white • Highly durable <b>Dimensions:</b> 25m x 2cm <b>Packaging requirements:</b> • 25m reel <b>Estimated usage: 0</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.2	156533	<b>Hook self adhesive 2,5cm (equivalent to velcro)</b> Hook – 2,5cm self adhesive • Maximum bonding capability • Colours: white • Highly durable <b>Dimensions:</b> 25m x 2,5cm <b>Packaging requirements:</b> • 25m reel <b>Estimated usage: 42</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.3	171270	<b>Hook self adhesive 2,5cm (equivalent to velcro)</b> Hook - 2,5cm self adhesive Maximum bonding capability • Colours: white • Highly durable <b>Dimensions:</b> 50m x 2,5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 45 rolls over range of sizes	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.4	156535	<b>Hook self adhesive 5cm (equivalent to velcro)</b> Hook – 5cm self adhesive • Maximum bonding capability • Colours: white • Highly durable <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 11</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.5	164583	<b>Hook self adhesive 5cm (equivalent to velcro)</b> Hook - 5cm self adhesive. • Maximum bonding capability • Colours: white • Highly durable <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 10 rolls over range of sizes	Yes, representative sample for items 50.1 to 50.15		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-50.6	156532	<b>Hook self adhesive 2cm (equivalent to velcro)</b> Hook – 2cm self adhesive • Maximum bonding capability • Colours: beige • Highly durable <b>Dimensions:</b> 25m x 2cm <b>Packaging requirements:</b> • 25m reel <b>Estimated usage: 1</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.7	156534	<b>Hook self adhesive 2,5cm (equivalent to velcro)</b> Hook – 2,5cm self adhesive • Maximum bonding capability • Colours: beige • Highly durable <b>Dimensions:</b> 25m x 2,5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 1</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.8	156536	<b>Hook self adhesive 5cm equivalent to velcro</b> Hook – 5cm self adhesive • Maximum bonding capability • Colours: beige • Highly durable <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 9</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.9	173614	<b>Hook self adhesive 2cm (equivalent to velcro)</b> Hook – 2cm self adhesive • Maximum bonding capability • Colours: Dark Brown • Highly durable <b>Dimensions:</b> 25m x 2cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 1</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.10	174272	<b>Hook self-adhesive 2,5cm (equivalent to velcro)</b> Hook – 2,5cm self adhesive • Maximum bonding capability • Colours: Dark Brown • Highly durable <b>Dimensions:</b> 25m x 2,5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 1</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.11	173616	<b>Hook self adhesive 5cm (equivalent to velcro)</b> Hook – 5cm self adhesive • Maximum bonding capability • Colours: Dark Brown • Highly durable <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 1</b>	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.12	164584	<b>Hook self adhesive 2,5cm (equivalent to velcro)</b> Hook - 2,5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable <b>Dimensions:</b> 50m x 2,5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 10 rolls over range of sizes	Yes, representative sample for items 50.1 to 50.15		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-50.13	164586	<b>Hook self adhesive 5cm (equivalent to velcro)</b> Hook - 5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated requirements:</b> 10 rolls over range of sizes	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.14	164557	<b>Hook self adhesive 2,5cm (equivalent to velcro)</b> Hook - 2,5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable <b>Dimensions:</b> 50m x 2,5cm <b>Packaging requirements:</b> Per roll <b>Estimated requirements:</b> 10 rolls over range of sizes	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM-50.15	171268	<b>Hook self adhesive 5cm (equivalent to velcro)</b> Hook - 5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable <b>Dimensions:</b> 50m x 5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 10 rolls over range of sizes	Yes, representative sample for items 50.1 to 50.15		Per roll														
ITEM 51	<b>Hook maximum bonding (equivalent to velcro)</b>																		
ITEM-51.1	156528	<b>Hook, maximum bonding, 2.5cm (equivalent to velcro)</b> Hook – 2.5cm wide • Highly durable • Maximum bonding capability • Colours: beige <b>Dimensions:</b> 25m x 2.5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 1</b>	Yes, representative sample for items 51.1 to 51.5		Per roll														
ITEM-51.2	156530	<b>Hook, maximum bonding, 5cm (equivalent to velcro)</b> Hook – 5cm wide • Highly durable • Maximum bonding capability • Colours: beige <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage: 1</b>	Yes, representative sample for items 51.1 to 51.5		Per roll														
ITEM-51.3	171299	<b>Hook, maximum bonding, 2.5cm (equivalent to velcro)</b> Hook - 2.5cm wide • Highly durable • Maximum bonding capability • Colours: brown <b>Dimensions:</b> 50m x 2.5cm <b>Packaging requirements:</b> Per roll <b>Estimated quantity:</b> 25 rolls over range of sizes	Yes, representative sample for items 51.1 to 51.5		Per roll														



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-51.4	173612	<b>Hook, maximum bonding, 2.5cm (equivalent to velcro)</b> Hook – 2.5cm wide • Highly durable • Maximum bonding capability • Colours: Dark Brown <b>Dimensions:</b> 25m x 2.5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> 1	Yes, representative sample for items 51.1 to 51.5		Per roll														
ITEM-51.5	173613	<b>Hook, maximum bonding, 5cm (equivalent to velcro)</b> Hook – 5cm wide • Highly durable • Maximum bonding capability • Colours: Dark Brown <b>Dimensions:</b> 25m x 5cm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> 1	Yes, representative sample for items 51.1 to 51.5		Per roll														
ITEM 52	<b>Self gripping (equivalent to velcro)</b>																		
ITEM-52.1	New	<b>Self gripping 25mm (equivalent to velcro)</b> Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. <b>Colour:</b> Tan/beige <b>Dimensions:</b> 25m x 25mm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> New	Yes, representative sample for items 52.1 to 52.6		Per roll														
ITEM-52.2	New	<b>Self gripping 50mm (equivalent to velcro)</b> Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. <b>Colour:</b> Tan/beige <b>Dimensions:</b> 25m x 50mm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> New	Yes, representative sample for items 52.1 to 52.6		Per roll														
ITEM-52.3	New	<b>Self gripping 25mm (equivalent to velcro)</b> Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. <b>Colour:</b> Black <b>Dimensions:</b> 25m x 25mm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> New	Yes, representative sample for items 52.1 to 52.6		Per roll														
ITEM-52.4	New	<b>Self gripping 50mm (equivalent to velcro)</b> Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. <b>Colour:</b> Black <b>Dimensions:</b> 25m x 50mm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> New	Yes, representative sample for items 52.1 to 52.6		Per roll														
ITEM-52.5	New	<b>Self gripping 25mm (equivalent to velcro)</b> Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. <b>Colour:</b> White <b>Dimensions:</b> 25m x 25mm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> New	Yes, representative sample for items 52.1 to 52.6		Per roll														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-52.6	New	<b>Self gripping 50mm (equivalent to velcro)</b> Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. <b>Colour:</b> White <b>Dimensions:</b> 25m x 50mm <b>Packaging requirements:</b> 25m roll <b>Estimated usage:</b> New	Yes, representative sample for items 52.1 to 52.6		Per roll														
ITEM 53	<b>Silicon Gel Sheets</b>																		
ITEM-53.1	118267	<b>Silicon Gel Sheets</b> Reduce the appearance of scars; Re-usable; Self- adhesive; Dimensions: 12 x 15cm <b>Packaging Requirements:</b> Pack of 10 sheets <b>Estimated usage:</b> 1	Yes		Per pack														
ITEM 54	<b>Silicon</b>																		
ITEM-54.1	New	<b>Silicon Gel 15ml</b> Reduce the appearance of scars; improves pliability and mobility Self drying silicon topical gel that is absorbable by the skin and can be applied to microporous tape Ingredients: Silicon gel, centella, asiatica <b>Dimensions:</b> 15ml tube <b>Packaging Requirements:</b> Pack of 5 <b>Estimated usage:</b> 1	Yes, representative sample for items 54.1 and 54.2		Per pack														
ITEM-54.2	New	<b>Silicon Gel 50ml</b> Reduce the appearance of scars; improves pliability and mobility Self drying silicon topical gel that is absorbable by the skin and can be applied to microporous tape Ingredients: Silicon gel, centella, asiatica <b>Dimensions:</b> 50ml tube <b>Packaging Requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 54.1 and 54.2		Each														
ITEM-54.3	New	<b>Silicon Spray</b> Treatment of scars to improve pliability, motility and appearance. Self drying silicon spray Ingredients: Silicon dioxide, polysiloxanes. <b>Dimensions:</b> 100ml Aerosol Can <b>Packaging Requirements:</b> Box of 6 <b>Estimated usage:</b> 1	Yes		Per box														
ITEM 55	<b>Squeeze Ball</b>																		
ITEM-55.1	118263	<b>Squeeze Ball Hand Exerciser 6 cm</b> Palm size ball; Returns to original diameter; <b>Dimensions:</b> 6cm diameter <b>Packaging Requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 56	<b>Pen and Pencil holder</b>																		
ITEM-56.1	118194	<b>Pen and Pencil holder</b> Soft PVC moulding that fits onto a normal sized pen/pencil to make them easier to grip and control writing. Dimensions: Internal Diameter 8mm. <b>Packaging Requirements:</b> Pack of 3 <b>Estimated usage:</b> 1	Yes		Per pack														
ITEM 57	<b>Long-Handled sponge</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-57.1	118193	<b>Long-Handled sponge</b> Mounted on plastic mouldable handle; Dimensions: Length 57.2cm. <b>Packaging Requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 57.1 and 57.2		Each														
ITEM-57.2	New	<b>Long-Handled sponge</b> Mounted on plastic mouldable handle; Dimensions: Length of more than 60cm <b>Packaging Requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 57.1 and 57.2		Each														
ITEM 58	<b>Sock Aid</b>																		
ITEM-58.1	118192	<b>Sock Aid</b> Cord attached assists in pulling up sock over foot; Cord style; Dimensions: Width 10cm; Cord style 25.5cm. <b>Packaging Requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 59	<b>Button hook and Zip puller</b>																		
ITEM-59.1	118191	<b>Button hook and Zip puller</b> Assist in tying buttons and pulling up zips with one hand. <b>Packaging Requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 60	<b>Plate guard</b>																		
ITEM-60.1	118741	<b>Plate guard</b> Fits around ordinary plate to assist with one handed eating; Microwave and Dishwasher safe. Dimensions: 19 to 25cm <b>Packaging Requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 61	<b>Finger coil spring</b>																		
ITEM-61.1	118190	<b>Finger coil spring</b> stainless steel spring; Gauge 0,7mm; Dimensions: Coil Diameter of 5mm. <b>Packaging Requirements:</b> Pack of 10 <b>Estimated usage:</b> 1	Yes		Per pack														
ITEM 62	<b>Knuckle bender coil spring</b>																		
ITEM-62.1	118188	<b>Knuckle bender coil spring</b> Stainless steel spring; Designed to assist intrinsic muscle weakness; Gauge 1mm. Torque force at full finger extension: ca.1000g.cm Dimensions: Coil Diameter 1cm. <b>Packaging Requirements:</b> Pack of 10 <b>Estimated usage:</b> 1	Yes		Per pack														
ITEM 63	<b>Finger extension assist</b>																		
ITEM-63.1	118187	<b>Finger extension assist</b> Stainless steel spring; Assist in extending fingers. Gauge 1mm. At total finger flexion the spring is tightened up to a torque force of : ca.2000g.cm Dimensions: Coil Diameter 4mm. <b>Packaging Requirements:</b> Pack of 10 <b>Estimated usage:</b> 1	Yes		Per pack														
ITEM 64	<b>Wrist extension coil spring</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-64.1	118180	<b>Wrist extension coil spring</b> Stainless steel spring; Two toes - left and right. Gauge 1,5 mm. Torque force when tightened over an angle of 130 degrees: 5000g.cm Diameter: 15mm. <b>Packaging Requirements:</b> Pack of 10 (5 left and 5 right) <b>Estimated usage: 1</b>	Yes		Per pack														
ITEM 65	<b>Plastic shoehorn</b>																		
ITEM 65.1	124872	<b>Plastic shoehorn</b> To assist in putting on shoes. Dimensions: Length 43cm; Weight 97g <b>Packaging Requirements:</b> Each <b>Estimated usage: 1</b>	Yes		Each														
ITM 66	<b>Webbing</b>																		
ITEM-66.1	156599	<b>Webbing 2.5cm</b> • Cotton • Non-elastic <b>Dimensions:</b> 2.5cm wide on a roll <b>Packaging requirements:</b> Single roll <b>Estimated usage: 6</b>	Yes, representative sample for items 66.1 and 66.2		Per roll														
ITEM-66.2	156600	<b>Webbing 5cm</b> • Cotton • Non-elastic <b>Dimensions:</b> 5cm wide on a roll <b>Packaging requirements:</b> Single roll <b>Estimated usage: 6</b>	Yes, representative sample for items 66.1 and 66.2		Per roll														
item 67	<b>Soft Cast</b>																		
ITEM-67.1	178631	<b>Soft Cast</b> • Semi-rigid casting tape <b>Dimensions:</b> 5cm width <b>Packaging requirements:</b> Box of 10 <b>Estimated usage: 1</b>	Yes, representative sample for items 67.1 to 67.3		Per box														
ITEM-67.2	113434	<b>Soft Cast</b> • Semi-rigid casting tape <b>Dimensions:</b> 7,5cm width <b>Packaging requirements:</b> Box of 10 <b>Estimated usage: 1</b>	Yes, representative sample for items 67.1 to 67.3		Per box														
ITEM-67.3	117664	<b>Soft Cast</b> • Semi-rigid casting tape <b>Dimensions:</b> 10cm width <b>Packaging requirements:</b> Box of 10 <b>Estimated usage: 1</b>	Yes, representative sample for items 67.1 to 67.3		Per box														
ITEM 68	<b>Outrigger Wire Copper Coated steel</b>																		
ITEM-68.1	113437	<b>Outrigger Wire Copper Coated steel, 1.6mm</b> Welding Rod 1.6mm Bends easily with pliers or with bender Length: 12 x 91cm <b>Packaging requirements:</b> Pack of 10 <b>Estimated usage: 1</b>	Yes, representative sample for items 68.1 to 68.3		Per pack														
ITEM-68.2	113723	<b>Outrigger Wire Copper Coated steel, 2.4mm</b> Welding Rod 2.4mm Bends easily with pliers or with bender Length: 12 x 91cm <b>Packaging requirements:</b> Pack of 10 <b>Estimated usage: 1</b>	Yes, representative sample for items 68.1 to 68.3		Per pack														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-68.3	113724	<b>Outrigger Wire Copper Coated steel, 3.2mm</b> Welding Rod 3.2mm Bends easily with pliers or with bender Length: 12 x 91cm <b>Packaging requirements:</b> Pack of 10 <b>Estimated usage:</b> 1	Yes, representative sample for items 68.1 to 68.3		Per pack														
ITEM 69	Tap turners																		
ITEM-69.1	156618	<b>Tap turners for cross type taps</b> • Strong moulded tap turner used to open a cross shaped tap head with leverage • Hooks over the top of a ‘+’ type tap, to make it easy to turn. • Supplied as a pair for hot and cold taps. <b>Dimensions:</b> Length ±140mm, Weight = 65g <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 54	Yes		Per pair														
ITEM-69.2	156619	<b>Tap turners for crystal type taps</b> • Highly durable tap turner which clamps around the top of cylindrical tap tops and provides the leverage to turn them easily • Barrel of the handle is spring loaded to allow the jaws to clamp down on the tap, which is gripped by soft non-slip rubber linings in the jaws. Jaws can also be locked in position. • Comfortable and easy to grip handle. • For use with taps up to 60 mm in diameters. <b>Dimensions:</b> Closed length = 190mm, Weight = 112g <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 70	Bath bench																		
ITEM-70.1	156628	<b>Bath bench</b> • To fit over top of bath to facilitate transfer into bath / enable client to sit above bath to wash with hand shower. • Must be adjustable to fit any bath and fit / secure safely to bath. <b>Dimensions:</b> To fit any bath <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 1	Yes		Each														
ITEM 71	Bath seat																		
ITEM-71.1	156629	<b>Bath seat</b> • Seat / bench which fits firmly into bath (seat ±½ way between bath rim and bottom of bath. • Adjustable to fit firmly into contour of various standard baths to enable client to sit safely above the bottom of the bath or transfer safely into the bottom of the bath. • Non-slip handles. <b>Dimensions:</b> To fit standard bath <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 4	Yes		Each														
ITEM 72	Rubber band post																		
ITEM-72.1	156605	<b>Rubber band post</b> • Used to anchor rubber bands <b>Dimensions:</b> 6,4mm x 13mm <b>Packaging requirements:</b> Pack of 10 <b>Estimated usage:</b> 10	Yes		Per pack														
ITEM 73	Closed cell foam tubing for assistive devices																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-73.1	156606	<b>Closed cell foam tubing for assistive devices, Size 1</b> <ul style="list-style-type: none"><li>• Closed cell foam tubing to build up tool &amp; utensil handles to give greater control to people who lack coordination, strength or gripping ability</li><li>• Slip resistant outer layer</li><li>• Dishwasher safe</li><li>• Can be cut to the required length</li><li>• Available in bright colours</li></ul> <b>Dimensions:</b> • Length = 300mm, • Inner diameter = 6mm, • Outer diameter 22mm <b>Packaging requirements:</b> Pack of 6 <b>Estimated usage:</b> 159	Yes, representative sample for items 73.1 to 73.3		Pack of 6														
ITEM-73.2	156607	<b>Closed cell foam tubing for assistive devices, Size 2</b> <ul style="list-style-type: none"><li>• Closed cell foam tubing to build up tool &amp; utensil handles to give greater control to people who lack coordination, strength or gripping ability</li><li>• Slip resistant outer layer</li><li>• Dishwasher safe</li><li>• Can be cut to the required length</li><li>• Available in bright colours</li></ul> <b>Dimensions:</b> • Length = 300mm, • Inner diameter = 10mm, • Outer diameter 28mm <b>Packaging requirements:</b> Pack of 6 <b>Estimated usage:</b> 158	Yes, representative sample for items 73.1 to 73.3		Pack of 6														
ITEM-73.3	156608	<b>Closed cell foam tubing for assistive devices, Size 3</b> <ul style="list-style-type: none"><li>• Closed cell foam tubing to build up tool &amp; utensil handles to give greater control to people who lack coordination, strength or gripping ability</li><li>• Slip resistant outer layer</li><li>• Dishwasher safe</li><li>• Can be cut to the required length</li><li>• Available in bright colours</li></ul> <b>Dimensions:</b> • Length = 300mm, • Inner diameter = 16mm, • Outer diameter 28mm <b>Packaging requirements:</b> Pack of 6 <b>Estimated usage:</b> 128	Yes, representative sample for items 73.1 to 73.3		Pack of 6														
ITEM 74	<b>Non slip rubber sheeting</b>																		
ITEM-74.1	156610	<b>Non slip rubber sheeting</b> <ul style="list-style-type: none"><li>• Non woven fabric, coated with non-slip properties.</li><li>• Suitable for moulding or wrapping around an object.</li></ul> <b>Dimensions:</b> 40cm x 9m x 1mm <b>Packaging requirements:</b> 9 meter roll <b>Estimated usage:</b> 3	Yes		Per roll														
ITEM 75	<b>Adapted knife with dagger handle</b>																		
ITEM-75.1	156615	<b>Adapted knife with dagger handle with serrated blade</b> <ul style="list-style-type: none"><li>• High quality stainless steel blade at ± 90° angle to handle to keep wrist in stress-free position.</li><li>• Serrated edge- ideal for slicing bread.</li><li>• Dishwasher safe.</li></ul> <b>Dimensions:</b> • Blade length = 254mm, • Weight 125g <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-75.2	156616	<b>Adapted knife with dagger handle with scalloped blade</b> <ul style="list-style-type: none"><li>• Stainless steel blade at ± 90° angle to handle to keep wrist in stress-free position.</li><li>• Permanently sharp blade</li><li>• Blade is ideal for cutting meat and vegetables</li><li>• Dishwasher safe.</li></ul> <b>Dimensions:</b> • Blade length = 203mm, • Weight 155g <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 76	<b>Long handled reacher</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-76.1	156617	<b>Long handled reacher</b> • Ultra light weight aluminium shaft • Assist those with restricted reach or mobility problems • Single moving jaw operated with a trigger action from the handle <b>Dimensions:</b> • Weight ±180 grams, • Length ± 600 mm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 9	Yes		Each														
ITEM 77	<b>Raised toilet seat</b>																		
ITEM-77.1	156625	<b>Raised toilet seat</b> • Light weight, plastic raised toilet seat with non-slip clips / brackets to attach firmly to / detach from most toilets easily. • Has a front cut-out to facilitate personal hygiene and a slight slope from back to front making sitting & rising easier. • Height-adjustable: can be raised 50mm, 100mm and 150mm above toilet rim • Moulded white polyethylene/ or similar compound that is chemically resilient, and easy to clean. <b>Dimensions:</b> • To fit standard toilet, • Height from toilet rim 50mm, 100mm and 150mm. <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 71	Yes		Each														
ITEM 78	<b>Nylon D rings</b>																		
ITEM-78.1	162856	<b>Nylon D rings, 25mm, White</b> • Nylon D-rings- rectangular • White <b>Dimensions:</b> 25 mm <b>Packaging requirements:</b> 50 per pack <b>Estimated usage:</b> 1	Yes, representative sample for items 78.1 and 78.2		Per packet														
ITEM-78.2	162857	<b>Nylon D rings, 50mm, White</b> • Nylon D-rings- rectangular • White <b>Dimensions:</b> 50mm <b>Packaging requirements:</b> 25 per pack <b>Estimated usage:</b> 1	Yes, representative sample for items 78.1 and 78.2		Per packet														
ITEM 79	<b>PATIENT TRANSFER ROLLBOARD</b>																		
ITEM-79.1	161149	<b>PATIENT TRANSFER ROLLBOARD, 170cm</b> A patient transfer roller-slide to provide a painless and comfortable transfer for patients with special needs. <input type="checkbox"/> Shape retaining sliding board covered with polyvinylchloride and PU-covered nylon <input type="checkbox"/> Able to bridge a gap of approx. 20cm, without bending, bearing a patient weighing up to 100kg <input type="checkbox"/> Allows for a patient to be transferred between to points even if there is a height difference <input type="checkbox"/> X-ray translucent <input type="checkbox"/> Foldable and lightweight <input type="checkbox"/> Removable machine washable cover, made of Ultra Low Friction material <input type="checkbox"/> Able to withstand surface disinfectants with 70% Ethanol, 1% Virkon or 45% Isopropyl alcohol <b>Dimensions:</b> Approx. 170cm x 50cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 10	Yes		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-79.2	161150	<b>PATIENT TRANSFER ROLLBOARD, 85cm</b> A patient transfer roller-slide to provide a painless and comfortable transfer for patients with special needs. <input type="checkbox"/> Shape retaining sliding board covered with polyvinylchloride and PU-covered nylon <input type="checkbox"/> Able to bridge a gap of approx. 20cm, without bending, bearing a patient weighing up to 100kg <input type="checkbox"/> Allows for a patient to be transferred between to points even if there is a height difference <input type="checkbox"/> X-ray translucent <input type="checkbox"/> Foldable and lightweight <input type="checkbox"/> Removable machine washable cover, made of Ultra Low Friction material <input type="checkbox"/> Able to withstand surface disinfectants with 70% Ethanol, 1% Virkon or 45% Isopropyl alcohol <b>Dimensions:</b> Approx. 85cm x 50cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 1	Yes		Each														
ITEM 80	<b>PATIENT TRANSFER SHEET / ROLL-OVER SLIDE TYPE</b>																		
ITEM-80.1	161174	<b>PATIENT TRANSFER SHEET / ROLL-OVER SLIDE TYPE, 50cm x 60cm</b> A patient transfer sheet to provide a painless and comfortable lateral transfer for patients, also to turn the patient. <input type="checkbox"/> A soft one piece padded roll bag which has flame retardant properties <input type="checkbox"/> Cotton outer material and inner Ultra low friction material <input type="checkbox"/> Has a polyurethane coated Nylon protective cover <input type="checkbox"/> Machine washable and or autoclavable <input type="checkbox"/> Able to withstand surface disinfectants with 45% Isopropyl alcohol <b>Dimensions:</b> Approx. 50cm x 60cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 0	Yes, representative sample for items 80.1 and 80.2		Each														
ITEM-80.2	161173	<b>PATIENT TRANSFER SHEET / ROLL-OVER SLIDE TYPE, 190cm x 60cm</b> A patient transfer sheet to provide a painless and comfortable lateral transfer for patients, also to turn the patient. <input type="checkbox"/> A soft one piece padded roll bag which has flame retardant properties <input type="checkbox"/> Cotton outer material and inner Ultra low friction material <input type="checkbox"/> Has a polyurethane coated Nylon protective cover <input type="checkbox"/> Machine washable and or autoclavable <input type="checkbox"/> Able to withstand surface disinfectants with 45% Isopropyl alcohol <b>Dimensions:</b> Approx. 190 x 60cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 1	Yes, representative sample for items 80.1 and 80.2		Each														
ITEM 81	<b>DISPOSABLE COVER FOR TRANSFER SHEET/ ROLL-OVER SLIDE TYPE</b>																		
ITEM-81.1	161175	<b>DISPOSABLE COVER FOR TRANSFER SHEET/ ROLL-OVER SLIDE TYPE, 190cm x 60cm</b> <input type="checkbox"/> Disposable covers shall be available for use with the transfer sheet <b>Dimensions:</b> Approx. 190 x 60cm <b>Packaging requirements:</b> Roll of 100 <b>Estimated usage:</b> 1	Yes, representative sample for items 81.1 and 81.2		Per roll														
ITEM-81.2	161176	<b>DISPOSABLE COVER FOR TRANSFER SHEET/ ROLL-OVER SLIDE TYPE, 50cm x 60cm</b> <input type="checkbox"/> Disposable covers shall be available for use with transfer sheet <b>Dimension:</b> 50cm X 60cm <b>Packaging requirements:</b> • Roll of 100 <b>Estimated usage:</b> 1	Yes, representative sample for items 81.1 and 81.2		Per roll														
ITEM 82	<b>PATIENT TRANNSFER SUPPORT AND MOVING AID</b>																		



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-82.1	161159	<b>PATIENT TRANSFER SUPPORT AND MOVING AID</b> A patient support and moving aid, assisting and supporting when raising and sitting the patient down, transfer to and from bed. <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface for patient comfort and safety <input type="checkbox"/> Two pairs of handles on either end, giving a choice of grip <b>Dimensions:</b> Approx. 50x 20cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 1	Yes		Each														
ITEM 83	<b>PATIENT SUPPORT AND MOVING BELT</b>																		
ITEM-83.1	161164	<b>PATIENT SUPPORT AND MOVING BELT, Small</b> Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles <input type="checkbox"/> Seven handles <b>Dimensions:</b> To accommodate waist sizes of 60 – 100cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 1	Yes, representative sample for items 83.1 to 83.4		Each														
ITEM-83.2	161165	<b>PATIENT SUPPORT AND MOVING BELT, Medium</b> Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles Nine handles <b>Dimensions:</b> To accommodate waist sizes of 72-120 cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 55	Yes, representative sample for items 83.1 to 83.4		Each														
ITEM-83.3	161166	<b>PATIENT SUPPORT AND MOVING BELT, Large</b> Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles Eleven handles <b>Dimensions:</b> To accommodate waist sizes of 100 – 160cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 5	Yes, representative sample for items 83.1 to 83.4		Each														
ITEM-83.4	161167	<b>PATIENT SUPPORT AND MOVING BELT, X-Large</b> Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles Thirteen handles <b>Dimensions:</b> To accommodate waist sizes of 115 – 175cm <b>Packaging requirements:</b> Single unit <b>Estimated usage:</b> 1	Yes, representative sample for items 83.1 to 83.4		Each	Each													
ITEM 84	<b>Patient glide board</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-84.1	161172	<b>Patient glide board</b> <input type="checkbox"/> For seating transfers <input type="checkbox"/> Manufactured from a Polyethylene and reinforced with non-rigid PVC <input type="checkbox"/> Flexible, allowing contouring to take place from different heights, during transfers from one surface to another <input type="checkbox"/> Elongated break allowing a side section to be bent upwards approx. 60 degree, to protect the patient during transfers <input type="checkbox"/> Anti-slip strips underneath <input type="checkbox"/> Handles at either end, for ease of use <input type="checkbox"/> Suitable for use during transfers, between bed and chair/wheelchair, between car and wheelchair and between wheelchair and toilet <input type="checkbox"/> Able to bridge gaps of approx. 10cm <input type="checkbox"/> Has no weight limit providing there is a stable surface, during transfer <input type="checkbox"/> Can be washed down, cleaned with alcohol or a surface disinfectant <b>Dimensions:</b> 75cm x 33cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 85	<b>Non-slip mat</b>																		
ITEM-85.1	New	<b>Non-slip mat</b> <input type="checkbox"/> Used on floor to prevent patients from slipping <input type="checkbox"/> Made of polyethylene Dimensions: Approx. 45 x 38cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 0	Yes		Each														
ITEM 86	<b>Transfer board</b>																		
ITEM-86.1	New	<b>Transfer board, wooden</b> <input type="checkbox"/> Tapered ends, light weight, smooth finish for easy sliding over board <input type="checkbox"/> Varnished Dimensions: 68cm x 30cm x 0.5cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 0	Yes		Each														
ITEM 87	<b>Conventional vacuum erection device</b>																		
ITEM-87.1	New	<b>Conventional vacuum erection device</b> <input type="checkbox"/> Kit must include pump, quick seal tension rings of different sizes, service oil, personal lubrication Dimensions: Standard <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 0	Yes		Each														
ITEM 88	<b>Wheelchair gloves (fingerless glove)</b>																		
ITEM-88.1	New	<b>Wheelchair gloves, heavy duty, extra small</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Extra Small <b>Packaging requirements:• per pair</b> <b>Estimated usage:</b> 1	Yes, representative sample for items 88.1 to 88.5		Per pair														
ITEM-88.2	New	<b>Wheelchair gloves, heavy duty, small</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Small <b>Packaging requirements:• per pair</b> <b>Estimated usage:</b> 1	Yes, representative sample for items 88.1 to 88.5		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-88.3	New	<b>Wheelchair gloves, heavy duty, medium</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Medium <b>Packaging requirements:• per pair</b> <b>Estimated usage: 1</b>	Yes, representative sample for items 88.1 to 88.5		Per pair														
ITEM-88.4	New	<b>Wheelchair gloves, heavy duty, large</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Large <b>Packaging requirements:• per pair</b> <b>Estimated usage: 1</b>	Yes, representative sample for items 88.1 to 88.5		Per pair														
ITEM-88.5	New	<b>Wheelchair gloves, heavy duty, X-Large</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: X-Large <b>Packaging requirements:• Per pair</b> <b>Estimated usage: 1</b>	Yes, representative sample for items 88.1 to 88.5		Per pair														
ITEM 89	<b>Wheelchair gloves (fully covered)</b>																		
ITEM-89.1	New	<b>Wheelchair gloves, heavy duty, extra small</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Extra Small <b>Packaging requirements:</b> Per pair <b>Estimated usage: 1</b>	Yes, representative sample for items 89.1 to 89.5		Per pair														
ITEM-89.2	New	<b>Wheelchair gloves, heavy duty, small</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Small <b>Packaging requirements:</b> Per pair <b>Estimated usage: 1</b>	Yes, representative sample for items 89.1 to 89.5		Per pair														
ITEM-89.3	New	<b>Wheelchair gloves, heavy duty, medium</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage: 1</b>	Yes, representative sample for items 89.1 to 89.5		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-89.4	New	<b>Wheelchair gloves, heavy duty, large</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 89.1 to 89.5		Per pair														
ITEM-89.5	New	<b>Wheelchair gloves, heavy duty, X-Large</b> <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: X-Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 89.1 to 89.5		Per pair														
ITEM 90	<b>Universal cuff</b>																		
ITEM-90.1	New	<b>Universal Cuff utensil pocket</b> <input type="checkbox"/> Made out of webbing material <input type="checkbox"/> Used to help patients with limited hand function to grasp eating utensils Has hook and loop closure on a D-Ring Includes utensil pocket Dimensions: can be adjusted to fit palms measuring 18cm to 23cm in circumference <b>Packaging requirements:</b> • Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-90.2	New	<b>Universal Cuff cylindrical foam</b> <input type="checkbox"/> Made out of webbing material <input type="checkbox"/> Used to help patients with limited hand function to grasp eating utensils Has hook and loop closure on a D-Ring Includes cylindrical foam to hold utensils Dimensions: can be adjusted to fit palms measuring 18cm to 23cm in circumference <b>Packaging requirements:</b> • Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 91	<b>Contact adhesive glue</b>																		
ITEM-91.1	New	<b>Contact adhesive glue</b> <input type="checkbox"/> Non - staining, quick drying, versatile on different material <input type="checkbox"/> Similar or equal to Genkem <b>Dimensions:</b> 500ml tin <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 0	Yes		Each														
ITEM 92	<b>Patient turn table</b>																		
ITEM-92.1	161674	<b>PATIENT TURN TABLE</b> Diameter. 38cm & 16mm thick <input type="checkbox"/> Used on floor to assist with turning patients that have difficulty in turning, or moving their feet. <input type="checkbox"/> Two discs manufactured from Polyethylene, fixed together in the centre on a axle. Bottom disc has non slip surface which prevents the patient from moving or slipping on it, Top disc rotates on axle, while bottom disc stays static. <input type="checkbox"/> Has carrying handle, <input type="checkbox"/> Easy to clean <b>Dimensions:</b> Diameter of 38cm, Thickness 12 mm with handle <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 93	<b>Cup specialized feeding</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-93.1	107265	<b>Cup specialized feeding transparent/semi</b> Cup, specialized feeding, should be transparent or semi-transparent with no sharp edge and the surface should be smooth, made of high grade polypropylene and be Bisphenol A (toxin) free,, must be able to withstand 125°C temp, must be calibrated in 5ml and must hold 50-60ml, must be autoclavable, the inner and outer surface of the cup must be reachable by hand and/or finger in order to clean, "Sinapi" equal or better. <b>Packaging requirements:</b> Per cup <b>Estimated usage:</b> 4,500 per year	Yes		Per cup														
ITEM 94	<b>Cup feeding system</b>																		
ITEM-94.1	107270	<b>Cup feeding system</b> Cup feeding system, consists of main holder, lid for the main holder and the feeding cup, main holder and feeding cup should be from virgin high-grade polypropylene material, be able to withstand a temperature of 120°C, the main holder must hold 200ml and have volume calibrations marked every 25ml, the feeding cup must hold 50 - 60ml, the feeding cup must be able to be stored inside the main holder with the lid in place, have no sharp edges and the surface texture should be smooth, the lid must fit watertight on the main holder, the inner and outer surface of the different components must be reachable by hand and/or finger in order to clean, the unit must be individually packaged and sealed in a plastic bag with an illustrated leaflet printed, "Sinapi" equal or better. <b>Packaging requirements:</b> Per system <b>Estimated usage:</b> 6,000 per year	Yes		Per system														
ITEM 95	<b>Breastmilk storage container, disposable</b>																		
ITEM-95.1	New	<b>Breastmilk storage container, disposable, 120-125ml</b> Breastmilk storage container, disposable, to store expressed breast milk. Manufactured from high density food grade polyethylene. BPA (Bisphenol-A) and DEHP (Di-Ethylhexyl Phthalate) free. Include foil cap. Include snap on cap made of polypropylene or polyethylene (water tight seal). Clean, filtered air used during blow moulding. Withstand pasteurisation up to 90°C. <b>Size:</b> 120-125ml capacity <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 5,000 per year	Yes, representative sample for items 95.1 and 95.2		Each														
ITEM-95.2	New	<b>Breastmilk storage container, disposable, 250ml</b> Breastmilk storage container, disposable, to store expressed breast milk. Manufactured from high density food grade polyethylene. BPA (Bisphenol-A) and DEHP (Di-Ethylhexyl Phthalate) free. Include foil cap. Include snap on cap made of polypropylene or polyethylene (water tight seal). Clean, filtered air used during blow moulding. Withstand pasteurisation up to 90°C. <b>Size:</b> 250ml capacity <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 10,000 per year	Yes, representative sample for items 95.1 and 95.2		Each														
<b>C:- PODIATRY</b>																			
ITEM 96	<b>Foam-o-felt</b>																		
ITEM-96.1	159280	<b>Foam-o-felt ,45x22.5cmx5mm</b> Foam-o-felt, mixture of 70% wool and 30% viscose (good for shock absorpion when natural fibrofatty pad has worn out), 45cm minimum length, 22.5cm width, 5mm thickness. <b>Dimensions:</b> 45cm x 22.5cm x 5mm <b>Packaging requirements:</b> Pack of 4 <b>Estimated usage:</b> 1	Yes		Per pack														
ITEM 98	<b>Poron</b>																		
ITEM-98.1	159282	<b>Poron, 31.5x22.5cmx3.2mm</b> Poron, soft supporting material, micro porous polyurethane with shock-absorbing properties, medical grade, 3.2mm thickness, 31.5cm length, 22.5cm width, per sheet, pack of 4 <b>Dimensions:</b> 31.5cm x 22.5cm x 3.2mm <b>Packaging requirements:</b> Pack of 4 <b>Estimated usage:</b> 72	Yes		Per pack														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM 99	PPT ultralux open cell																		
ITEM-99-1	159287	<b>PPT ultralux open cell, 152x30cmx1.6mm</b> PPT ultralux, durable open cell material that doesn't bottom out, 152cm length, 30cm width, 1.6mm thickness, good shock absorber, ideal covering material <b>Dimensions:</b> 152cm x 30cm x 1.6mm <b>Packaging requirements:</b> Pack of 2 <b>Estimated usage:</b> 1	Yes		Per pack														
ITEM 100	Cortex cork material																		
ITEM-100.1	159150	<b>Cortex cork material, 1mx600x1.33mm</b> Cortex, cork material, 1m length per sheet, 600mm width, 1.33mm thickness, used as base for corrective foot devices <b>Dimensions:</b> 1m x 600mm x 1.33mm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 1	Yes, representative sample for items 100.1 and 100.2		Per sheet														
ITEM-100.2	159151	<b>Cortex cork material, 1mx600x1.5mm</b> Cortex, cork material, 1m length per sheet, 600mm width, 1.5mm thickness, used as base for corrective foot devices <b>Dimensions:</b> 1m x 600mm x 1.5mm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 50	Yes, representative sample for items 100.1 and 100.2		Per sheet														
ITEM 101	Neoprene material																		
ITEM-101.1	178709	<b>Neoprene material, 1mx3mm</b> Neoprene material 1m length, 3mm thick <b>Dimensions:</b> 1m x 0,5cm <b>Packaging requirements:</b> Per sheet <b>Estimated usage:</b> 1	Yes		Per sheet														
ITEM 102	Adhesive contact liquid																		
ITEM-102.1	136798	<b>Adhesive contact liquid 1L</b> Adhesive contact, used for laminating all components of corrective foot devices, liquid, 1L <b>Size:</b> 1L <b>Packaging requirements:</b> Per 1L <b>Estimated usage:</b> 1	Yes, representative sample for items 102.1 and 102.2		Per L														
ITEM-102.2	159162	<b>Adhesive contact liquid 5L</b> Adhesive contact, used for laminating all components of corrective foot devices, liquid, 5L <b>Size:</b> 5L <b>Packaging requirements:</b> Per 5L <b>Estimated usage:</b> 10	Yes, representative sample for items 102.1 and 102.2		Per L														
ITEM 103	Remover solvent																		
ITEM-103.1	159163	<b>Remover solvent to dissolve adhesive, 5L</b> Remover solvent, chemical to dissolve contact adhesive, 5L <b>Size:</b> 5L <b>Packaging requirements:</b> Per 5L <b>Estimated usage:</b> 10	Yes		Per L														
ITEM 104	Fleecy web, cotton																		
ITEM-104.1	159283	<b>Fleecy web, cotton, 40 x 22.5cm</b> Fleecy web, 100% cotton padding with a raised fleecy surface, flexible application, 40cm length, 22.5cm width, pack of 4 <b>Dimensions:</b> 40 x 22.5cm <b>Packaging requirements:</b> Pack of 4 <b>Estimated usage:</b> 1	Yes		per pack														
ITEM 105	Felt, semi-compressed																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-105.1	114120	<b>Felt, semi-compressed, 3mm</b> Felt, semi-compressed, 100% pure wool, used as a highly versatile padding which can be easily cut and shaped, providing excellent pressure absorption. <b>Size:</b> 3mm <b>Packaging requirements:</b> Pack of 4 <b>Estimated usage:</b> 1	Yes, representative sample for items 105.1 and 105.2		per pack														
ITEM-105.2	114121	<b>Felt, semi-compressed, 5mm</b> Felt, semi-compressed, 100% pure wool, used as a highly versatile padding which can be easily cut and shaped, providing excellent pressure absorption. <b>Size:</b> 5mm <b>Packaging requirements:</b> Pack of 4 <b>Estimated usage:</b> 1	Yes, representative sample for items 105.1 and 105.2		Per pack														
ITEM 106	Tubifoam sponge tubular																		
ITEM-106.1	159284	<b>Tubifoam, sponge tubular, 18mmx25cm</b> Tubifoam, sponge, tubular, 18mm width, 25cm length, requires no strapping or tape, protects toes from friction, per tube (pack of 10 tubes) <b>Dimensions:</b> 18mm x 25cm <b>Packaging requirements:</b> Pack of 10 tubes <b>Estimated usage:</b> 1	Yes, representative sample for items 106.1 to 106.3		Per pack														
ITEM-106.2	159285	<b>Tubifoam, sponge tubular, 21mmx25cm</b> Tubifoam, sponge, tubular 21mm width, 25cm length, requires no strapping or tape, protects toes from friction, per tube (pack of 10 tubes) <b>Dimensions:</b> 21mm x 25cm <b>Packaging requirements:</b> Pack of 10 tubes <b>Estimated usage:</b> 1	Yes, representative sample for items 106.1 to 106.3		Per pack														
ITEM-106.3	159286	<b>Tubifoam, sponge tubular, 25mmx25cm</b> Tubifoam, sponge, tubular, 25mm width, 25cm length, requires no strapping or tape, protects toes from friction, per tube (pack of 10 tubes) <b>Dimensions:</b> 25mm x 25cm <b>Packaging requirements:</b> Pack of 10 tubes <b>Estimated usage:</b> 1	Yes, representative sample for items 106.1 to 106.3		Per pack														
ITEM 107	Heel lifts, PPT pad																		
ITEM-107.1	136777	<b>Heel lifts, PPT pad, medium</b> Heel lifts, PPT pad indicated for heel spurs, plantar fasciitis, very soft density to provide relief to sensitive areas. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 107.1 to 107.2		Per pair														
ITEM-107.2	136778	<b>Heel lifts, PPT pad, large</b> Heel lifts, PPT pad indicated for heel spurs, plantar fasciitis, very soft density to provide relief to sensitive areas. <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 107.1 to 107.2		Per pair														
ITEM 108	Heel raise																		
ITEM-108.1	159288	<b>Heel raise, small</b> Heel raise, indicated for leg length difference and heel pain, polyurethane, density 0.25-0.28 grams/cm <sup>3</sup> , hardness of 30-50° shore A, covering material is micro fibre material. <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 108.1 to 108.3		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-108.2	159289	<b>Heel raise, medium</b> Heel raise, indicated for leg length difference and heel pain, polyurethane, density 0.25-0.28 grams/cm <sup>3</sup> , hardness of 40-70° shore A, covering material is micro fibre material. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 108.1 to 108.3		Per pair														
ITEM-108.3	159300	<b>Heel raise, large</b> Heel raise, indicated for leg length difference and heel pain, polyurethane, density 0.25-0.28 grams/cm <sup>3</sup> , hardness of 40-70° shore A, covering material is micro fibre material. <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 108.1 to 108.3		Per pair														
ITEM 109	<b>Heel pad, horse shoe</b>																		
ITEM-109.1	114147	<b>Heel pad, horse shoe, small</b> Heel pad, horse shoe, indicated for off loading, latex sponge rubber, 8-10mm, with density of 0.25 - 0.28 grams/cm <sup>3</sup> , hardness of 40-70° shore A. <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 109.1 to 109.3		Per pair														
ITEM-109.2	114148	<b>Heel pad, horse shoe, medium</b> Heel pad, horse shoe, indicated for off loading, latex sponge rubber, 8-10mm, with density of 0.25 - 0.28 grams/cm <sup>3</sup> , hardness of 40-70° shore A. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 109.1 to 109.3		Per pair														
ITEM-109.3	114149	<b>Heel pad, horse shoe, large</b> Heel pad, horse shoe, indicated for off loading, latex sponge rubber, 8-10mm, with density of 0.25 - 0.28 grams/cm <sup>3</sup> , hardness of 40-70° shore A. <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 109.1 to 109.3		Per pair														
ITEM 110	<b>Heel pad horse shoe with extractable central core</b>																		
ITEM-110.1	114150	<b>Heel pad horse shoe with extractable central core</b> Heel pad, horse shoe, with extractable central core to relived calcaneal spurs, indicated for heel spurs and plantar fasciitis and off loading, latex sponge rubber, 5-10mm, with density of 0.25 - 0.28 grams/cm <sup>3</sup> , hardness of 15-20° shore A. <b>Size:</b> One size fits all <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes		Per pair														
ITEM 111	<b>Cobra pad with valgus pad extension</b>																		
ITEM-111.1	114153	<b>Cobra pad with valgus pad extension, small</b> Cobra pad with valgus pad extension, sponge rubber, density of 0.25 - 0.28 grams/cm <sup>3</sup> , hardness of 15-20° shore A. <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 111.1 to 111.3		Per pair														
ITEM-111.2	114152	<b>Cobra pad with valgus pad extension, medium</b> Cobra pad with valgus pad extension, sponge rubber, density of 0.25 - 0.28 grams/cm <sup>3</sup> , hardness of 15-20° shore A. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 111.1 to 111.3		Per pair														



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-111.3	114151	<b>Cobra pad with valgus pad extension, large</b> Cobra pad with valgus pad extension, sponge rubber, density of 0.25 - 0.28 grams/cm <sup>3</sup> , hardness of 15-20° shore A. <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 111.1 to 111.3		Per pair														
ITEM 112	<b>Metatarsal bar, PPT pad</b>																		
ITEM-112.1	159301	<b>Metatarsal bar, PPT pad, small</b> Metatarsal bar, PPT pad, small, allowing for extra protection under metatarsal heads (prevents growth of hard callouses and reduces pain and irritation). <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 112.1 to 112.3		Per pair														
ITEM-112.2	159302	<b>Metatarsal bar, PPT pad, medium</b> Metatarsal bar, PPT pad, medium, allowing for extra protection under metatarsal heads (prevent growth of hard areas and callouses and reduces pain and irritation). <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 112.1 to 112.3		Per pair														
ITEM-112.3	159303	<b>Metatarsal bar, PPT pad, large</b> Metatarsal bar, PPT pad, large, allowing for extra protection under metatarsal heads (prevents growth of hard areas and callouses and reduces pain and irritation) <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 50	Yes, representative sample for items 112.1 to 112.3		Per pair														
ITEM 113	<b>Metatarsal bar, synthetic latex pad with cut-out</b>																		
ITEM-113.1	114122	<b>Metatarsal bar, synthetic latex pad with cut-out, small</b> Metatarsal bar, synthetic latex pad, small, with a cut-out on the first and fifth metatarsal, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm <sup>3</sup> <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 113.1 to 113.3		Per pair														
ITEM-113.2	114123	<b>Metatarsal bar, synthetic latex pad with cut-out, medium</b> Metatarsal bar, synthetic latex pad, medium, with a cut-out on the first and fifth metatarsal, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm <sup>3</sup> <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 113.1 to 113.3		Per pair														
ITEM-113.3	114124	<b>Metatarsal bar, synthetic latex pad with cut-out, large</b> Metatarsal bar, synthetic latex pad, large, with a cut-out on the first and fifth metatarsal, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm <sup>3</sup> <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 113.1 to 113.3		Per pair														
ITEM 114	<b>Metatarsal bar, synthetic latex pad, straight</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-114.1	114125	<b>Metatarsal bar, synthetic latex pad, straight, small</b> Metatarsal bar, synthetic latex pad, small, straight across metatarsal heads, with no cut-out on the first and fifth metatarsal, med pad made of synthetic latex foam, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm <sup>3</sup> <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 114.1 and 114.2		Per pair														
ITEM-114.2	114126	<b>Metatarsal bar, synthetic latex pad, straight, large</b> Metatarsal bar, synthetic latex pad, large, straight across metatarsal heads, with no cut-out on the first and fifth metatarsal, med pad made of synthetic latex foam, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm <sup>3</sup> <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 114.1 and 114.2		Per pair														
ITEM 115	<b>Metatarsal pad with arch support, synthetic latex</b>																		
ITEM-115.1	114142	<b>Metatarsal pad with arch support, synthetic latex, small</b> Metatarsal pad with arch support, made of synthetic latex, with hardness of 15-20° shore A, density of 0.25gram/cm <sup>3</sup> <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 115.1 to 115.3		Per pair														
ITEM-115.2	114143	<b>Metatarsal pad with arch support, synthetic latex, medium</b> Metatarsal pad with arch support, made of synthetic latex, with hardness of 15-20° shore A, density of 0.25gram/cm <sup>3</sup> <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 115.1 to 115.3		Per pair														
ITEM-115.3	114144	<b>Metatarsal pad with arch support, synthetic latex, large</b> Metatarsal pad with arch support, made of synthetic latex, with hardness of 15-20° shore A, density of 0.25gram/cm <sup>3</sup> <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 115.1 to 115.3		Per pair														
ITEM 116	<b>Pad, metatarsal, foam rubber</b>																		
ITEM-116.1	159159	<b>Pad, metatarsal, foam rubber, small</b> Pad, metatarsal, soft foam rubber, made of foam rubber for its supportive and cushioning properties, small <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 116.1 to 116.3		Per pair														
ITEM-116.2	159160	<b>Pad, metatarsal foam, rubber, medium</b> Pad, metatarsal, soft foam rubber, made of foam rubber for its supportive and cushioning properties, small <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 6	Yes, representative sample for items 116.1 to 116.3		Per pair														
ITEM-116.3	159161	<b>Pad, metatarsal foam, rubber, large</b> Pad, metatarsal, soft foam rubber, made of foam rubber for its supportive and cushioning properties, small <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 6	Yes, representative sample for items 116.1 to 116.3		Per pair														
ITEM 117	<b>Pads, valgus arch support, EVA</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-117.1	159156	<b>Pads, valgus arch support, EVA, small</b> Pads valgus, arch support, made of EVA of soft density sponge rubber, small <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 117.1 to 117.3		Per pair														
ITEM-117.2	159157	<b>Pads, valgus arch support, EVA, medium</b> Pads valgus, arch support, made of EVA of soft density sponge rubber, small <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 6	Yes, representative sample for items 117.1 to 117.3		Per pair														
ITEM-117.3	159158	<b>Pads, valgus arch support, EVA, large</b> Pads valgus, arch support, made of EVA of soft density sponge rubber, small <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 6	Yes, representative sample for items 117.1 to 117.3		Per pair														
ITEM 118	<b>Pads, valgus arch support, synthetic latex, hard</b>																		
ITEM-118.1	114139	<b>Pads, valgus arch support, synthetic latex, hard, small</b> Pads, valgus arch support, made of synthetic latex foam, hard, density of 30-35° shore, hardness 1 grams/cm <sup>3</sup> <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 118.1 to 118.3		Per pair														
ITEM-118.2	114140	<b>Pads, valgus arch support, synthetic latex, hard, medium</b> Pads, valgus arch support, made of synthetic latex foam, hard, density of 30-35° shore, hardness 1 grams/cm <sup>3</sup> <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 118.1 to 118.3		Per pair														
ITEM-118.3	114138	<b>Pads, valgus arch support, synthetic latex, hard, large</b> Pads, valgus arch support, made of synthetic latex foam, hard, density of 30-35° shore, hardness 1 grams/cm <sup>3</sup> . <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 118.1 to 118.3		Per pair														
ITEM 119	<b>Pads, valgus arch support, synthetic latex, soft</b>																		
ITEM-119.1	114136	<b>Pads, valgus arch support, synthetic latex, soft, small</b> Pads, valgus arch support, made of synthetic latex foam, soft, density between 15-20° shore <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 119.1 to 119.3		Per pair														
ITEM-119.2	114137	<b>Pads, valgus arch support, synthetic latex, soft, medium</b> Pads, valgus arch support, made of synthetic latex foam, soft, density between 15-20° shore <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 119.1 to 119.3		Per pair														
ITEM-119.3	114141	<b>Pads, valgus arch support, synthetic latex, soft, large</b> Pads, valgus arch support, made of synthetic latex foam, soft, density between 15-20° shore <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 119.1 to 119.3		Per pair														
ITEM 120	<b>Digi-tube, ribbed elastic, strip</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-120.1	159305	<b>Digi-tube, ribbed, elastic strip, small</b> Digi-tube, ribbed elastic material, tubular strip, fully lined with silicone, small, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) <b>Size:</b> Small <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 30	Yes, representative sample for items 120.1 to 120.3		Each														
ITEM-120.2	159307	<b>Digi-tube, ribbed, elastic strip, medium</b> Digi-tube, ribbed elastic material, tubular strip, fully lined with silicone, medium, durable, washable which can be cut into different lengths. <b>Size:</b> Medium <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 20	Yes, representative sample for items 120.1 to 120.3		Each														
ITEM-120.3	159308	<b>Digi-tube, ribbed elastic, strip, large</b> Digi-tube, ribbed elastic material, tubular strip, fully lined with silicone, large, durable, washable which can be cut into different lengths. <b>Size:</b> Large <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 8	Yes, representative sample for items 120.1 to 120.3		Each														
ITEM 121	<b>Digi-tube, pure gel, strip</b>																		
ITEM-121.1	136779	<b>Digi-tube, pure gel, strip, small</b> Digi-tube, pure gel, tubular strip, small, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) <b>Size:</b> Small <b>Packaging requirements:</b> Each <b>Estimated usage:</b> New item	Yes, representative sample for items 121.1 to 121.3		Each														
ITEM-121.2	136780	<b>Digi-tube, pure gel, strip, medium</b> Digi-tube, pure gel, tubular strip, medium, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) <b>Size:</b> Medium <b>Packaging requirements:</b> Each <b>Estimated usage:</b> New item	Yes, representative sample for items 121.1 to 121.3		Each														
ITEM-121.3	136781	<b>Digi-tube, pure gel, strip, large</b> Digi-tube, pure gel, tubular strip, large, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) <b>Size:</b> Large <b>Packaging requirements:</b> Each <b>Estimated usage:</b> New item	Yes, representative sample for items 121.1 to 121.3		Each														
ITEM 122	<b>Toe loop</b>																		
ITEM-122.1	114157	<b>Toe loop, narrow</b> Toe loop, narrow loop, foam-lined material with velcro strap, ideal for repositioning the toes, strap with loops secures two or more toes together to stabilise deformed toes <b>Size:</b> Narrow loop, 1.5 x 7.6cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 122.1 and 122.2		Each														
ITEM-122.2	114158	<b>Toe loop, wider</b> Toe loop, wider loop, foam-lined material with velcro strap, ideal for repositioning the toes, strap with loops secures two or more toes together to stabilise deformed toes <b>Size:</b> Wider loop, 4 x 10cm <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes, representative sample for items 122.1 and 122.2		Each														
ITEM 123	<b>Toe bar</b>																		
ITEM-123.1	114154	<b>Toe bar, small</b> Toe bar, with adjustable elastic band that fits around toes, made of chamois leather, recommended for hammer toes. <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 123.1 to 123.3		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-123.2	114155	<b>Toe bar, medium</b> Toe bar, with adjustable elastic band that fits around toes, made of chamois leather, recommended for hammer toes. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 123.1 to 123.3		Per pair														
ITEM-123.3	114156	<b>Toe bar, large</b> Toe bar, with adjustable elastic band that fits around toes, made of chamois leather, recommended for hammer toes. <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 123.1 to 123.3		Per pair														
ITEM 124	<b>CAP digital, pure gel</b>																		
ITEM-124.1	159720	<b>CAP digital, pure gel, small</b> Cap digital, pure gel, small, eases pressure and friction caused by blisters, corns and it is closed in front, transparent, durable and washable, single <b>Size:</b> Small <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 1	Yes, representative sample for items 124.1 to 124.3		Each														
ITEM-124.2	159719	<b>CAP digital, pure gel, medium</b> Cap digital, pure gel, medium, eases pressure and friction caused by blisters, corns and it is closed in front, transparent, durable and washable, single <b>Size:</b> Medium <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 1	Yes, representative sample for items 124.1 to 124.3		Each														
ITEM-124.3	159722	<b>CAP digital, pure gel, large</b> Cap digital, pure gel, large, eases pressure and friction caused by blisters, corns and it is closed in front, transparent, durable and washable, single <b>Size:</b> Large <b>Packaging requirements:</b> Per item <b>Estimated usage:</b> 1	Yes, representative sample for items 124.1 to 124.3		Each														
ITEM 125	<b>Spreader toe, pure gel interdigital wedge</b>																		
ITEM-125.1	159797	<b>Spreader toe, pure gel interdigital wedge, small</b> Spreader toe, pure gel interdigital wedge, small, for deformities and overlapping toes, protection in helomas and inter-digital pressure, helps with re-alignment of toes <b>Size:</b> Small <b>Packaging item:</b> Per pack of 12 <b>Estimated usage:</b> 1	Yes, representative sample for items 125.1 to 125.3		Per packet														
ITEM-125.2	159798	<b>Spreader toe, pure gel interdigital wedge, medium</b> Spreader toe, pure gel interdigital wedge, medium, for deformities and overlapping toes, protection in helomas and inter-digital pressure, helps with re-alignment of toes. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pack of 12 <b>Estimated usage:</b> 1	Yes, representative sample for items 125.1 to 125.3		Per packet														
ITEM-125.3	159800	<b>Spreader toe, pure gel interdigital wedge, large</b> Spreader toe, pure gel interdigital wedge, large, for deformities and overlapping toes, protection in heloma and inter-digital, pressure, helps with re-alignment of toes. <b>Size:</b> Large <b>Packaging requirements:</b> Per pack of 12 <b>Estimated usage:</b> 1	Yes, representative sample for items 125.1 to 125.3		Per packet														
ITEM 126	<b>Spreader, hammer, made of soft polymer gel, big toe spreader</b>																		
ITEM-126.1	137019	<b>Spreader, hammer, made of soft polymer gel, big toe spreader, small</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Small <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 6	Yes, representative sample for items 126.1 to 126.3		Per packet														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-126.2	137020	<b>Spreader, hammer, made of soft polymer gel, big toe spreader, medium</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Medium <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 6	Yes, representative sample for items 126.1 to 126.3		Per packet														
ITEM-126.3	170844	<b>Spreader, hammer, made of soft polymer gel, big toe spreader, large</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Large <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 126.1 to 126.3		Per packet														
ITEM 127	<b>Spreader, hammer, made of soft polymer gel, advanced toe spreader</b>																		
ITEM-127.1	170848	<b>Spreader, hammer, made of soft polymer gel, advanced toe spreader, small</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Small <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 127.1 to 127.3		Per packet														
ITEM-127.2	170849	<b>Spreader, hammer, made of soft polymer gel, advanced toe spreader, medium</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Medium <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 127.1 to 127.3		Per packet														
ITEM-127.3	170860	<b>Spreader, hammer, made of soft polymer gel, advanced toe spreader, large</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Large <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 127.1 to 127.3		Per packet														
ITEM 128	<b>Spreader, hammer, made of soft polymer gel, universal toe spreader</b>																		
ITEM-128.1	170845	<b>Spreader, hammer, made of soft polymer gel, universal toe spreader, small</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Small <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 128.1 to 128.4		Per packet														
ITEM-128.2	170846	<b>Spreader, hammer, made of soft polymer gel, universal toe spreader, medium</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Medium <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 128.1 to 128.4		Per packet														
ITEM-128.3	170847	<b>Spreader, hammer, made of soft polymer gel, universal toe spreader, large</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> Large <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 128.1 to 128.4		Per packet														
ITEM-128.4	170839	<b>Spreader, hammer, made of soft polymer gel, universal toe spreader, X-Large</b> Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units <b>Size:</b> X-Large <b>Packaging requirements:</b> Pack of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 128.1 to 128.4		Per packet														
ITEM 129	<b>Pad, hammer, polymer gel</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-129.1	137021	<b>Pad, hammer, polymer gel</b> Hammer toes crest pad, made of polymer gel, soft, good stretch, relieves forefoot and toe pain caused by hammer, claw and mallet toes, toe crest provides maximum comfort and reduces pressure form the toes and metatarsals <b>Packaging requirements:</b> Pack of 10 units <b>Estimated usage:</b> 1	Yes		Per packet														
ITEM 130	<b>Separator toe, gel mould fitting</b>																		
ITEM-130.1	159803	<b>Separator toe, gel mould fitting, small</b> Separator toe, gel mould fitting, small, to reduce interdigital pressure between toes, protection from inter-digital pressure, helps with re-alignment of toes. <b>Size:</b> Small <b>Packaging requirements:</b> Per pack of 12 <b>Estimated usage:</b> 1	Yes, representative sample for items 130.1 to 130.3		Per packet														
ITEM-130.2	178702	<b>Separator toe, gel mould fitting, medium</b> Separator toe, gel mould fitting, medium, to reduce interdigital pressure between toes, protection from inter-digital pressure, helps with re-alignment of toes. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pack of 12 <b>Estimated usage:</b> 1	Yes, representative sample for items 130.1 to 130.3		Per packet														
ITEM-130.3	159805	<b>Separator toe, gel mould fitting, large</b> Separator toe, gel mould fitting, large, to reduce interdigital pressure between toes, protection from inter-digital pressure, helps with re-alignment of toes. <b>Size:</b> Large <b>Packaging requirements:</b> Per pack of 12 <b>Estimated usage:</b> 1	Yes, representative sample for items 130.1 to 130.3		Per packet														
ITEM 131	<b>Shield bunion, foam shell, gel disc</b>																		
ITEM-131.1	159712	<b>Shield bunion, foam shell, gel disc, small</b> Shield bunion, toes, soft foam shell for pressure absorption and gel disc for pressure positioning, small, protects the bunion joint. Reduces pain in the bone protuberance, thin silicone for better patient compliance. <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 131.1 to 131.3		Per pair														
ITEM-131.2	159713	<b>Shield bunion, foam shell, gel disc, medium</b> Shield bunion, toes, soft foam shell for pressure absorption and gel disc for pressure positioning, medium, protects the bunion joint. Reduces pain in bone protuberance, thin silicone for patient compliance. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 131.1 to 131.3		Per pair														
ITEM-131.3	159714	<b>Shield bunion, foam shell, gel disc, large</b> Shield bunion, toes, soft foam shell for pressure absorption and gel disc for pressure positioning, large, protects the bunion joint. Reduces pain in bone protuberance, thin silicone for better patient compliance. <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 131.1 to 131.3		Per pair														
ITEM 132	<b>Shield, bunion, polymer gel, soft</b>																		
ITEM-132.1	170827	<b>Shield, bunion, polymer gel, soft, small</b> Shield, bunion, toes, made of <u>soft</u> polymer gel, good stretch, in 2 levels of density, protector is slim, yet very effective and can be used with any shoe, protectors help guard a painful bunion from rubbing and pressure and for protection with loop around big toe. <b>Dimensions:</b> Length ± 2.5cm <b>Size:</b> Small <b>Packaging requirements:</b> Pack of 6 units <b>Estimated usage:</b> 1	Yes, representative sample for items 132.1 to 132.3		Per packet														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-132.2	137066	<b>Shield, bunion, polymer gel, soft, medium</b> Shield, bunion, toes, made of <u>soft</u> polymer gel, good stretch, in 2 levels of density, protector is slim, yet very effective and can be used with any shoe, protectors help guard a painful bunion from rubbing and pressure and for protection with loop around big toe. <b>Dimensions:</b> Length ± 2.5cm <b>Size:</b> Medium <b>Packaging requirements:</b> Pack of 6 units <b>Estimated usage:</b> 1	Yes, representative sample for items 132.1 to 132.3		Per packet														
ITEM-132.3	170829	<b>Shield, bunion, polymer gel, soft, large</b> Shield, bunion, toes, made of <u>soft</u> polymer gel, good stretch, in 2 levels of density, protector is slim, yet very effective and can be used with any shoe, protectors help guard a painful bunion from rubbing and pressure and for protection with loop around big toe. <b>Dimensions:</b> Length ± 2.5cm <b>Size:</b> Large <b>Packaging requirements:</b> Pack of 6 units <b>Estimated usage:</b> 1	Yes, representative sample for items 132.1 to 132.3		Per packet														
ITEM 133	<b>Shield bunion, silicone</b>																		
ITEM-133.1	178468	<b>Shield bunion, silicone, small</b> Shield bunion toes layer of silicone gel placed to protect prominent big toe joint, covered with elastic material to keep silicone gel in place, reduce pain in bone protuberance to relief pressure <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 133.1 to 133.3		Per pair														
ITEM-133.2	178470	<b>Shield bunion, silicone, medium</b> Shield bunion toes layer of silicone gel placed to protect prominent big toe joint, covered with elastic material to keep silicone gel in place, reduce pain in bone protuberance to relief pressure <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 133.1 to 133.3		Per pair														
ITEM-133.3	178472	<b>Shield bunion, silicone, large</b> Shield bunion toes layer of silicone gel placed to protect prominent big toe joint, covered with elastic material to keep silicone gel in place, reduce pain in bone protuberance to relief pressure <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 133.1 to 133.3		Per pair														
ITEM 134	<b>Splint, Bunion / Hallux Valgus</b>																		
ITEM-134.1	114145	<b>Splint, Bunion / Hallux Valgus</b> Splint, for correction of bunion / hallux valgus, keeps first radius of foot in correct position, thin enough to wear in shoes, non-slip, breathable foam, holds the correction even during weight-bearing. <b>Size:</b> One size fits all <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 135	<b>Arch lift plantar fasciitis</b>																		
ITEM-135.1	114146	<b>Arch lift plantar fasciitis</b> Arch lift, for plantar fasciitis, foam lined material with velcro straps, alleviates symptoms of plantar fasciitis and heel spurs, provides upward compression to the arch, supporting and stabilising the planter fascia ligaments and includes a built-in pad for additional support to medial arch, helping to prevent excessive pronation. <b>Size:</b> One size fits all <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 136	<b>CAP digital, polymer gel, soft</b>																		



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-136.1	137063	<b>CAP digital, polymer gel, soft, small</b> CAP digital textile fabric inside coated with polymer gel, <u>soft</u> , stretchable, protects toes against rubbing and pressure, transparent, durable, washable, due to its softness, elasticity and moistening qualities (medical oil + vaseline). It works to cushion apical areas, as digital amputation cover, in-growing nail care, as keloid scar prevention. <b>Size:</b> Small <b>Packaging requirements:</b> Packet of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 136.1 to 136.3		Per packet														
ITEM-136.2	137064	<b>CAP digital, polymer gel, soft, medium</b> CAP digital textile fabric inside coated with polymer gel, <u>soft</u> , stretchable, protects toes against rubbing and pressure, transparent, durable, washable, due to its softness, elasticity and moistening qualities (medical oil + vaseline). It works to cushion apical areas, as digital amputation cover, in-growing nail care, as keloid scar prevention. <b>Size:</b> Medium <b>Packaging requirements:</b> Packet of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 136.1 to 136.3		Per packet														
ITEM-136.3	137065	<b>CAP digital, polymer gel, soft, large</b> CAP digital textile fabric inside coated with polymer gel, <u>soft</u> , stretchable, protects toes against rubbing and pressure, transparent, durable, washable, due to its softness, elasticity and moistening qualities (medical oil + vaseline). It works to cushion apical areas, as digital amputation cover, in-growing nail care, as keloid scar prevention. <b>Size:</b> Large <b>Packaging requirements:</b> Packet of 12 units <b>Estimated usage:</b> 1	Yes, representative sample for items 136.1 to 136.3		Per packet														
ITEM 137	<b>Heel Spur, polymer gel, soft</b>																		
ITEM-137.1	170811	<b>Heel Spur, polymer gel, soft, small</b> Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. <b>Size:</b> Small <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 137.1 to 137.4		Per pair														
ITEM-137.2	137058	<b>Heel Spur, polymer gel, soft, medium</b> Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. <b>Size:</b> Medium <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 137.1 to 137.4		Per pair														
ITEM-137.3	137059	<b>Heel Spur, polymer gel, soft, large</b> Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. <b>Size:</b> Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 137.1 to 137.4		Per pair														
ITEM-137.4	174308	<b>Heel Spur, Polymer gel, soft, X-Large</b> Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. <b>Size:</b> X-Large <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 137.1 to 137.4		Per pair														
ITEM 138	<b>Rings, digital, textile, polymer gel</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-138.1	170822	<b>Rings, digital, textile, polymer gel, small</b> Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of <u>soft</u> polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. <b>Length:</b> ± 2.5cm <b>Size:</b> Small <b>Packaging requirements:</b> Pack of 10 rings <b>Estimated usage:</b> 1	Yes, representative sample for items 138.1 to 138.4		Per packet														
ITEM-138.2	170823	<b>Rings, digital, textile, polymer gel, medium</b> Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of <u>soft</u> polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. <b>Length:</b> ± 2.5cm <b>Size:</b> Medium <b>Packaging requirements:</b> Pack of 10 rings <b>Estimated usage:</b> 1	Yes, representative sample for items 138.1 to 138.4		Per packet														
ITEM-138.3	170825	<b>Rings, digital, textile, polymer gel, large</b> Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of <u>soft</u> polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. <b>Length:</b> ± 2.5cm <b>Size:</b> Large <b>Packaging requirements:</b> Pack of 10 rings <b>Estimated usage:</b> 1	Yes, representative sample for items 138.1 to 138.4		Per packet														
ITEM-138.4	170826	<b>Rings, digital, textile, polymer gel, X-Large</b> Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of <u>soft</u> polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. <b>Length:</b> ± 2.5cm <b>Size:</b> X-Large <b>Packaging requirements:</b> Pack of 10 rings <b>Estimated usage:</b> 1	Yes, representative sample for items 138.1 to 138.4		Per packet														
ITEM 139	<b>Silicone putty kit, soft mouldable for relief</b>																		
ITEM-139.1	159806	<b>Silicone putty kit, soft, mouldable for relief</b> Silicone putty kit, mouldable for pressure-relief and correcting digital deformities, normal density, shore density of between 12-14 , 250g base with 250g catalyst, mixed 1:1; it is soft, quick setting, self-mixing, pliable and washable. <b>Size:</b> Base with catalyst = 500g product <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> 1	Yes		Per kit														
ITEM 140	<b>Silicone putty kit, hard, mouldable for relief</b>																		
ITEM-140.1	114127	<b>Silicone putty kit, hard, mouldable for relief</b> Silicone putty kit, mouldable for pressure-relief and correcting digital deformities, normal density, shore density of 20-25, 250g base with 250g catalyst, mixed 1:1; it is hard, quick setting, self-mixing, pliable and washable. <b>Size:</b> Base with catalyst = 500g product <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> 1	Yes		Per kit														
ITEM 141	<b>Insole, shoe, thin moulded EVA</b>																		
ITEM-141.1	159695	<b>Insole, shoe size 4, thin moulded EVA</b> Insole, shoe size 4, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 4 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-141.2	159697	<b>Insole, shoe size 5, thin moulded EVA</b> Insole, shoe size 5, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 5 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM-141.3	136782	<b>Insole, shoe size 6, thin moulded EVA</b> Insole, shoe size 6, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 6 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM-141.4	136783	<b>Insole, shoe size 7, thin moulded EVA</b> Insole, shoe size 7, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 7 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM-141.5	136784	<b>Insole, shoe size 8, thin moulded EVA</b> Insole, shoe size 8, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 8 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM-141.6	159700	<b>Insole, shoe size 9, thin moulded EVA</b> Insole, shoe size 9, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 9 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM-141.7	136786	<b>Insole, shoe size 10, thin moulded EVA</b> Insole, shoe size 10, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 10 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM-141.8	159702	<b>Insole, shoe size 11, thin moulded EVA</b> Insole, shoe size 11, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 11 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM-141.9	159703	<b>Insole, shoe size 12, thin moulded EVA</b> Insole, shoe size 12, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. <b>Size:</b> Shoe size 12 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 141.1 to 141.9		Per pair														
ITEM 142	<b>Insole, shoe, thin cork</b>																		
ITEM-142.1	136787	<b>Insole, shoe size 4, thin cork</b> Insole, shoe size 4, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 4 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> New item	Yes, representative sample for items 142.1 to 142.9		Per pair														
ITEM-142.2	136788	<b>Insole, shoe size 5, thin cork</b> Insole, shoe size 5, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 5 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> New item	Yes, representative sample for items 142.1 to 142.9		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-142.3	136789	<b>Insole, shoe size 6, thin cork</b> Insole, shoe size 6, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 6 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> New item	Yes, representative sample for items 142.1 to 142.9		Per pair														
ITEM-142.4	136790	<b>Insole, shoe size 7, thin cork</b> Insole, shoe size 7, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 7 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> New item	Yes, representative sample for items 142.1 to 142.9		Per pair														
ITEM-142.5	136791	<b>Insole, shoe size 8, thin cork</b> Insole, shoe size 8, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 8 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> New item	Yes, representative sample for items 142.1 to 142.9		Per pair														
ITEM-142.6	136792	<b>Insole, shoe size 9, thin cork</b> Insole, shoe size 9, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 9 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 142.1 to 142.9		Per pair														
ITEM-142.7	136793	<b>Insole, shoe size 10, thin cork</b> Insole, shoe size 10, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 10 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> New item	Yes, representative sample for items 142.1 to 142.9		Per pair														
ITEM-142.8	136794	<b>Insole, shoe size 11, thin cork</b> Insole, shoe size 11, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 11 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 142.1 to 142.9		Per pair														
ITEM-142.9	136795	<b>Insole, shoe size 12, thin cork</b> Insole, shoe size 12, thin cork, pre-cut, 1-2mm thick. <b>Size:</b> Shoe size 12 <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 142.1 to 142.9		Per pair														
item 143	<b>Innersole, polyurethane</b>																		
ITEM-143.1	114128	<b>Innersole, polyurethane, shoe size 4, full length</b> Innersole, polyurethane foam, shoe size 4, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 4 full length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.2	114130	<b>Innersole, polyurethane, shoe size 5, full length</b> Innersole, polyurethane foam, shoe size 5, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 5 full length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-143.3	159696	<b>Innersole, polyurethane, shoe size 6, full length</b> Innersole, polyurethane foam, shoe size 6, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 6 full length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.4	159698	<b>Innersole, polyurethane, shoe size 7, full length</b> Innersole, polyurethane foam, shoe size 7, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 7 full length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.5	159699	<b>Innersole, polyurethane, shoe size 8, full length</b> Innersole, polyurethane foam, shoe size 8, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 8 full length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.6	114135	<b>Innersole, polyurethane, shoe size 9, full length</b> Innersole, polyurethane foam, shoe size 9, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 9 full length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.7	159701	<b>Innersole, polyurethane, shoe size 10, full length</b> Innersole, polyurethane foam, shoe size 10, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is black micro fibre material. <b>Size:</b> Shoe size 10 full length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.8	114129	<b>Innersole, polyurethane, shoe size 4, ¾ length</b> Innersole, polyurethane foam, shoe size 4, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 4 ¾ length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.9	114131	<b>Innersole, polyurethane, shoe size 5, ¾ length</b> Innersole, polyurethane foam, shoe size 5, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 5 ¾ length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-143.10	114132	<b>Innersole, polyurethane, shoe size 6, ¾ length</b> Innersole, polyurethane foam, shoe size 6, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 6 ¾ length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.11	114133	<b>Innersole, polyurethane, shoe size 7, ¾ length</b> Innersole, polyurethane foam, shoe size 7, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 7 ¾ length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM-143.12	114134	<b>Innersole, polyurethane, shoe size 8, ¾ length</b> Innersole, polyurethane foam, shoe size 8, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm <sup>3</sup> , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. <b>Size:</b> Shoe size 8 ¾ length <b>Packaging requirements:</b> Per pair <b>Estimated usage:</b> 1	Yes, representative sample for items 143.1 to 143.12		Per pair														
ITEM 144	<b>Orthotic kit, medium, adhesive</b>																		
ITEM-144.1	New	<b>Orthotic kit, medium, adhesive, shoe size 5, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 5 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 144.1 to 144.7		Per kit														
ITEM-144.2	New	<b>Orthotic kit, medium, adhesive, shoe size 6, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 6 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 144.1 to 144.7		Per kit														
ITEM-144.3	New	<b>Orthotic kit, medium, adhesive, shoe size 7, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 7 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 144.1 to 144.7		Per kit														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-144.4	New	<b>Orthotic kit, medium, adhesive, shoe size 8, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 8 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 144.1 to 144.7		Per kit														
ITEM-144.5	New	<b>Orthotic kit, medium, adhesive, shoe size 9, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 9 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 144.1 to 144.7		Per kit														
ITEM-144.6	New	<b>Orthotic kit, medium, adhesive, shoe size 10, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 10 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 144.1 to 144.7		Per kit														
ITEM-144.7	New	<b>Orthotic kit, medium, adhesive, shoe size 11, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 11 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 144.1 to 144.7		Per kit														
ITEM 145	<b>Orthotic kit, hard, adhesive</b>																		
ITEM-145.1	New	<b>Orthotic kit, hard, adhesive, shoe size 5, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 5 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 145.1 to 145.7		Per kit														
ITEM-145.2	New	<b>Orthotic kit, hard, adhesive, shoe size 6, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 6 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 145.1 to 145.7		Per kit														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-145.3	New	<b>Orthotic kit, hard, adhesive, shoe size 7, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 7 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 145.1 to 145.7		Per kit														
ITEM-145.4	New	<b>Orthotic kit, hard, adhesive, shoe size 8, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 8 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 145.1 to 145.7		Per kit														
ITEM-145.5	New	<b>Orthotic kit, hard, adhesive, shoe size 9, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 9 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 145.1 to 145.7		Per kit														
ITEM-145.6	New	<b>Orthotic kit, hard, adhesive, shoe size 10, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 10 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 145.1 to 145.7		Per kit														
ITEM-145.7	New	<b>Orthotic kit, hard, adhesive, shoe size 11, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 11 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 145.1 to 145.7		Per kit														
ITEM 146	<b>Orthotic kit, medium, non-adhesive</b>																		
ITEM-146.1	New	<b>Orthotic kit, medium, non-adhesive, shoe size 5, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 5 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 146.1 to 146.7		Per kit														



ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-146.2	New	<b>Orthotic kit, medium, non-adhesive, shoe size 6, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 6 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 146.1 to 146.7		Per kit														
ITEM-146.3	New	<b>Orthotic kit, medium, non-adhesive, shoe size 7, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 7 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 146.1 to 146.7		Per kit														
ITEM-146.4	New	<b>Orthotic kit, medium, non-adhesive, shoe size 8, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 8 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 146.1 to 146.7		Per kit														
ITEM-146.5	New	<b>Orthotic kit, medium, non-adhesive, shoe size 9, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 9 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 146.1 to 146.7		Per kit														
ITEM-146.6	New	<b>Orthotic kit, medium, non-adhesive, shoe size, 10 full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 10 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 146.1 to 146.7		Per kit														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-146.7	New	<b>Orthotic kit, medium, non-adhesive, shoe size 11, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>medium</b> . <b>Size:</b> Shoe size 11 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 146.1 to 146.7		Per kit														
ITEM-147	<b>Orthotic kit, hard. non-adhesive</b>																		
ITEM-147.1	New	<b>Orthotic kit, hard, non-adhesive, shoe size, 5 full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 5 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 147.1 to 147.7		Per kit														
ITEM-147.2	New	<b>Orthotic kit, hard, non-adhesive, shoe size, 6 full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 6 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 147.1 to 147.7		Per kit														
ITEM-147.3	New	<b>Orthotic kit, hard, non-adhesive, shoe size 7, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 7 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 147.1 to 147.7		Per kit														
ITEM-147.4	New	<b>Orthotic kit, hard, non-adhesive, shoe size 8, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 8 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 147.1 to 147.7		Per kit														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATERIAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-147.5	New	<b>Orthotic kit, hard, non-adhesive, shoe size 9, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 9 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 147.1 to 147.7		Per kit														
ITEM-147.6	New	<b>Orthotic kit, hard, non-adhesive, shoe size 10, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 10 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 147.1 to 147.7		Per kit														
ITEM-147.7	New	<b>Orthotic kit, hard, non-adhesive, shoe size 11, full length</b> Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex <b>non-adhesive</b> metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: <b>hard</b> . <b>Size:</b> Shoe size 11 full length <b>Packaging requirements:</b> Per kit <b>Estimated usage:</b> New	Yes, representative sample for items 147.1 to 147.7		Per kit														
ITEM 148	<b>Burr diamond crystals</b>																		
ITEM-148.1	159098	<b>Burr diamond crystals barrel head shape</b> Burr, diamond crystals bonded directly to one piece, barrel head shape, hardened stainless steel shank <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-148.2	159095	<b>Burr diamond crystals bud head shape</b> Burr, diamond crystals bonded directly to one piece, bud head shape, hardened stainless steel shank <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-148.3	159099	<b>Burr diamond crystals cylinder-ball head</b> Burr, diamond crystals bonded directly to one piece, cylinder-ball head shape, hardened stainless steel shank <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-148.4	159096	<b>Burr diamond crystals enucleator head</b> Burr, diamond crystals bonded directly to one piece, enucleator head shape, hardened stainless steel shank <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-148.5	159093	<b>Burr diamond crystals olive head shape</b> Burr, diamond crystals bonded directly to one piece, olive head shape, hardened stainless steel shank <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCALLY MANUFACTURED
ITEM-148.6	159094	<b>Burr diamond crystals pear head shape</b> Burr, diamond crystals bonded directly to one piece, pear head shape, hardened stainless steel shank <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM-148.7	136797	<b>Burr tungsten carbide cone shape</b> Burr, tungsten carbide, cone shape with extremely fine grain size, calibrated to high isostatic pressing process, blades are accurately sharpened ensuring optimal cutting, carbide is very hard with high longevity, better sterility, hardened stainless steel shank <b>Packaging requirements:</b> Each <b>Estimated usage:</b> New item	Yes		Each														
ITEM 149	<b>Disc garnet</b>																		
ITEM-149.1	159136	<b>Disc garnet 2.25cm diameter course grid</b> Disc garnet, disposable emery disc for finishing skin and calluses with square brass centre that fits onto a mandrill stem, course grid, 2.25cm diameter, 50 units per pack <b>Dimensions:</b> 2.25cm diameter, <b>Packaging requirements:</b> 50 units per pack <b>Estimated usage:</b> 1	Yes		Per packet														
ITEM-149.2	159135	<b>Disc garnet 2.25cm diameter medium grid</b> Disc garnet, disposable emery disc for finishing skin and calluses with square brass centre that fits onto a mandrill stem, medium grid, 2.25cm diameter, 50 units per pack <b>Dimensions:</b> 2.25cm diameter, <b>Packaging requirements:</b> 50 units per pack <b>Estimated usage:</b> 1	Yes		Per packet														
ITEM 150	<b>Mandrill stem</b>																		
ITEM-150.1	159133	<b>Mandrill stem standard dental size</b> Mandrill stem, standard dental size, stainless steel shank to hold the garnet disc <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 151	<b>Silver Nitrate</b>																		
ITEM-151.1	162969	<b>Silver Nitrate stick</b> Silver nitrate sticks, 40% strength <b>Size:</b> per stick <b>Packaging requirements:</b> Pack of 10 sticks <b>Estimated usage:</b> 1	Yes		Per packet														
ITEM-151.2	137022	<b>Silver Nitrate liquid</b> 40% strength <b>Size:</b> 10ml (dropper bottle) <b>Packaging requirements:</b> Per bottle <b>Estimated usage:</b> 1	Yes		Per bottle														
ITEM-151.3	136790	<b>Silver Nitrate granules</b> Silver nitrate granules, clear and odourless crystals, very soluble in water, 40% strength <b>Size:</b> 5g granules in 50ml bottle <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 152	<b>Oil spray</b>																		
ITEM-152.1	178163	<b>Oil spray, anti rust</b> Anti rust oil spray for stainless steel instruments (pre-sterilization) t increase life span of instruments <b>Size:</b> 300ml container <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 1	Yes		Each														
ITEM 153	<b>Blade, scalpel, surgical</b>																		

ITEM NO	SAP NO	ITEM DESCRIPTION	SAMPLE REQUIRED	COMPLY YES/NO	UNIT OF MEASURE	TENDER PRICE/UNIT	STATE PACKAGING /BOX QUANTITY	BRAND	LEAD TIME (DAYS)	DIRECT MATEIRAL: % IMPORTED	DIRECT MATERIAL: % LOCAL CONTENT	%DIRECT LABOUR	%TRANSPORT	%OTHER OVERHEADS AND %PROFITS	%TOTAL	RATE OF EXCHANGE	ADJUSTMENT FREQUENCY	IMPORTED	LOCAL MANUFACTURED
ITEM-153.1	New	<b>Blade, scalpel, surgical, no. 10</b> Blade, scalpel, surgical, to fit handle no's 3, 5 and 7, with small fitting features, no. 10, must be compatible with scalpel handles for items (short no. 3, long no. 5 and long no.7), composition test for carbon steel and stainless steel, ISO900101 and ISO13485. Must be sterile, peel apart pack must peel evenly without tearing paper lining, blade size must be printed on the peel apart pack. <b>Size:</b> No. 10 blade <b>Packaging requirements:</b> Packet of 100 <b>Estimated usage:</b> 100 packets	Yes		Per packet														
ITEM-153.2	New	<b>Blade, scalpel, surgical, no. 15</b> Blade, scalpel, surgical, to fit handle no's 3, 5 and 7, with small fitting features, no. 15, must be compatible with scalpel handles for items (short no. 3, long no. 5 and long no.7), composition test for carbon steel and stainless steel, ISO900101 and ISO13485. Must be sterile, peel apart pack must peel evenly without tearing paper lining, blade size must be printed on the peel apart pack. <b>Size:</b> No. 15 blade <b>Packaging requirements:</b> Packet of 100 <b>Estimated usage:</b> 100 packets	Yes		Per packet														
ITEM 154	<b>Monofilament</b>																		
ITEM-154.1	124117	<b>Monofilament, retractable, 10g</b> Monofilament, nylon, retractable, 10g. <b>Size:</b> 10g <b>Packaging requirements:</b> Each <b>Estimated usage:</b> 100	Yes		Each														
ITEM 155	<b>Tourniquet</b>																		
ITEM-155.1	New	<b>Tourniquet, digital, small</b> Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: orange <b>Size:</b> Small, 0.6 - 0.9cm diameter <b>Packaging requirements:</b> Box of 20 <b>Estimated usage:</b> 0	Yes, representative sample for items 155.1 to 155.4		Per box														
ITEM-155.2	New	<b>Tourniquet, digital, medium</b> Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: red <b>Size:</b> Medium, ±1cm diameter <b>Packaging requirements:</b> Box of 20 <b>Estimated usage:</b> 1	Yes, representative sample for items 155.1 to 155.4		Per box														
ITEM-155.3	New	<b>Tourniquet, digital, large</b> Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: green <b>Size:</b> Large, ±1.25cm diameter <b>Packaging requirements:</b> Box of 20 <b>Estimated usage:</b> 1	Yes, representative sample for items 155.1 to 155.4		Per box														
ITEM-155.4	New	<b>Tourniquet, digital, X-Large</b> Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: blue <b>Size:</b> X-Large, ±1.5cm diameter <b>Packaging requirements:</b> Box of 20 <b>Estimated usage:</b> 1	Yes, representative sample for items 155.1 to 155.4		Per box														



## Provincial Supply Chain Management

# SAMPLE SHEET

Page 1 of 2

**SAMPLES FOR BID NO:**

PLEASE TICK WHERE SAMPLE IS SUBMITTED

[illegible]

NAME OF CONTACT PERSON

**CONTACT NO**



## Provincial Supply Chain Management

# SAMPLE SHEET

Page 2 of 2

**PLEASE TICK WHERE SAMPLE IS SUBMITTED**

[illegible]



Page 1 of 2

--

[illegible]

TEL NUMBER






## Provincial Supply Chain Management

# Sample Sheet

Page 2 of 2

**PLEASE TICK WHERE SAMPLE IS SUBMITTED**

[illegible]

 <b>GAUTENG PROVINCE</b> PROVINCIAL TREASURY REPUBLIC OF SOUTH AFRICA	<b>Provincial Supply Chain Management</b>	
	<b>PREFERENCE POINTS CLAIM FORM</b>	Page: 1 of 7

**PREFERENCE POINTS CLAIM FORM IN TERMS OF THE PREFERENTIAL PROCUREMENT REGULATIONS 2022**

This preference form must form part of all tenders invited. It contains general information and serves as a claim form for preference points for specific goals.

**NB: BEFORE COMPLETING THIS FORM, TENDERERS MUST STUDY THE GENERAL CONDITIONS, DEFINITIONS AND DIRECTIVES APPLICABLE IN RESPECT OF THE TENDER AND PREFERENTIAL PROCUREMENT REGULATIONS, 2022**

**1. GENERAL CONDITIONS**

1.1 The following preference point systems are applicable to invitations to tender:

- the 80/20 system for requirements with a Rand value of up to R50 000 000 (all applicable taxes included); and
- the 90/10 system for requirements with a Rand value above R50 000 000 (all applicable taxes included).

1.2 **To be completed by the organ of state**

[TICK APPLICABLE BOX]

	The applicable preference point system for this tender is the 90/10 preference point system.
	The applicable preference point system for this tender is the 80/20 preference point system.
	Either the 90/10 or 80/20 preference point system will be applicable in this tender. The lowest/ highest acceptable tender will be used to determine the accurate system once tenders are received.

1.3 Points for this tender (even in the case of a tender for income-generating contracts) shall be awarded for:

- (a) Price; and
- (b) Specific Goals.



# Provincial Supply Chain Management

## PREFERENCE POINTS CLAIM FORM

Page: 2 of 7

### 1.4 To be completed by the organ of state:

The maximum points for this tender are allocated as follows:

	POINTS
<b>PRICE</b>	
<b>SPECIFIC GOALS</b>	
<b>Total points for Price and SPECIFIC GOALS</b>	<b>100</b>

1.5 Failure on the part of a tenderer to submit proof or documentation required in terms of this tender to claim points for specific goals with the tender, will be interpreted to mean that preference points for specific goals are not claimed.

1.6 The organ of state reserves the right to require of a tenderer, either before a tender is adjudicated or at any time subsequently, to substantiate any claim in regard to preferences, in any manner required by the organ of state.

## 2. DEFINITIONS

- (a) **“tender”** means a written offer in the form determined by an organ of state in response to an invitation to provide goods or services through price quotations, competitive tendering process or any other method envisaged in legislation;
- (b) **“price”** means an amount of money tendered for goods or services, and includes all applicable taxes less all unconditional discounts;
- (c) **“rand value”** means the total estimated value of a contract in Rand, calculated at the time of bid invitation, and includes all applicable taxes;
- (d) **“tender for income-generating contracts”** means a written offer in the form determined by an organ of state in response to an invitation for the origination of income-generating contracts through any method envisaged in legislation that will result in a legal agreement between the organ of state and a third party that produces revenue for the organ of state, and includes, but is not limited to, leasing and disposal of assets and concession contracts, excluding direct sales and disposal of assets through public auctions; and
- (e) **“the Act”** means the Preferential Procurement Policy Framework Act, 2000 (Act No. 5 of 2000).



# Provincial Supply Chain Management

## PREFERENCE POINTS CLAIM FORM

Page: 3 of 7

### 3. FORMULAE FOR PROCUREMENT OF GOODS AND SERVICES

#### 3.1. POINTS AWARDED FOR PRICE

##### 3.1.1 THE 80/20 OR 90/10 PREFERENCE POINT SYSTEMS

A maximum of 80 or 90 points is allocated for price on the following basis:

80/20

or

90/10

$$Ps = 80 \left( 1 - \frac{Pt - P_{min}}{P_{min}} \right) \text{ or } Ps = 90 \left( 1 - \frac{Pt - P_{min}}{P_{min}} \right)$$

Where

Ps = Points scored for price of tender under consideration

Pt = Price of tender under consideration

Pmin = Price of lowest acceptable tender

#### 3.2. FORMULAE FOR DISPOSAL OR LEASING OF STATE ASSETS AND INCOME GENERATING PROCUREMENT

##### 3.2.1. POINTS AWARDED FOR PRICE

A maximum of 80 or 90 points is allocated for price on the following basis:

80/20

or

90/10

$$Ps = 80 \left( 1 + \frac{Pt - P_{max}}{P_{max}} \right) \text{ or } Ps = 90 \left( 1 + \frac{Pt - P_{max}}{P_{max}} \right)$$

Where

Ps = Points scored for price of tender under consideration

Pt = Price of tender under consideration

Pmax = Price of highest acceptable tender



# Provincial Supply Chain Management

## PREFERENCE POINTS CLAIM FORM

Page: 4 of 7

### 4. POINTS AWARDED FOR SPECIFIC GOALS

- 4.1. In terms of Regulation 4(2); 5(2); 6(2) and 7(2) of the Preferential Procurement Regulations, preference points must be awarded for specific goals stated in the tender. For the purposes of this tender the tenderer will be allocated points based on the goals stated in table 1 below as may be supported by proof/ documentation stated in the conditions of this tender:
- 4.2. In cases where organs of state intend to use Regulation 3(2) of the Regulations, which states that, if it is unclear whether the 80/20 or 90/10 preference point system applies, an organ of state must, in the tender documents, stipulate in the case of—
- (a) an invitation for tender for income-generating contracts, that either the 80/20 or 90/10 preference point system will apply and that the highest acceptable tender will be used to determine the applicable preference point system; or
  - (b) any other invitation for tender, that either the 80/20 or 90/10 preference point system will apply and that the lowest acceptable tender will be used to determine the applicable preference point system,
- then the organ of state must indicate the points allocated for specific goals for both the 90/10 and 80/20 preference point system.



# Provincial Supply Chain Management

## PREFERENCE POINTS CLAIM FORM

Page: 5 of 7

**Table 1: Specific goals for the tender and points claimed are indicated per the table below.**

**(Note to organs of state: Where either the 90/10 or 80/20 preference point system is applicable, corresponding points must also be indicated as such.)**

**Note to tenderers: The tenderer must indicate how they claim points for each preference point system.)**

The specific goals allocated points in terms of this tender	Number of points allocated (90/10 system) (To be completed by the organ of state)	Number of points allocated (80/20 system) (To be completed by the organ of state)	Number of points claimed (90/10 system) (To be completed by the tenderer)	Number of points claimed (80/20 system) (To be completed by the tenderer)



# Provincial Supply Chain Management

## PREFERENCE POINTS CLAIM FORM

Page: 6 of 7

### DECLARATION WITH REGARD TO COMPANY/FIRM

4.3. Name of company/firm

4.4. Company registration number

4.5. TYPE OF COMPANY/ FIRM

[TICK APPLICABLE BOX]

<input type="checkbox"/>	Partnership/Joint Venture / Consortium
<input type="checkbox"/>	One-person business/sole propriety
<input type="checkbox"/>	Close corporation
<input type="checkbox"/>	Public Company
<input type="checkbox"/>	Personal Liability Company
<input type="checkbox"/>	(Pty) Limited
<input type="checkbox"/>	Non-Profit Company
<input type="checkbox"/>	State Owned Company

4.6. I, the undersigned, who is duly authorised to do so on behalf of the company/firm, certify that the points claimed, based on the specific goals as advised in the tender, qualifies the company/ firm for the preference(s) shown and I acknowledge that:

- i) The information furnished is true and correct;
- ii) The preference points claimed are in accordance with the General Conditions as indicated in paragraph 1 of this form;
- iii) In the event of a contract being awarded as a result of points claimed as shown in paragraphs 1.4 and 4.2, the contractor may be required to furnish documentary proof to the satisfaction of the organ of state that the claims are correct;
- iv) If the specific goals have been claimed or obtained on a fraudulent basis or any of the conditions of contract have not been fulfilled, the organ of state may, in addition to any other remedy it may have –



# Provincial Supply Chain Management

## PREFERENCE POINTS CLAIM FORM

Page: 7 of 7

- (a) disqualify the person from the tendering process;
- (b) recover costs, losses or damages it has incurred or suffered as a result of that person's conduct;
- (c) cancel the contract and claim any damages which it has suffered as a result of having to make less favourable arrangements due to such cancellation;
- (d) recommend that the tenderer or contractor, its shareholders and directors, or only the shareholders and directors who acted on a fraudulent basis, be restricted from obtaining business from any organ of state for a period not exceeding 10 years, after the *audi alteram partem* (hear the other side) rule has been applied; and
- (e) forward the matter for criminal prosecution, if deemed necessary.

<b>SIGNATURE(S) OF TENDERER(S)</b>	
<b>SURNAME AND NAME</b>	
<b>DATE</b>	
<b>ADDRESS</b>	