



Provincial Supply Chain Management

Price Schedule - Goods

Page 1 of 3

Bidder Name

Bid No

NOTE: If the information required regarding each item cannot be inserted in the space provided on either page 1, 2 or 3 of this form, additional information may be attached to this form.

The Gauteng Provincial Treasury (GPT) is acting as Common Service Provider or buying organisation on behalf of all Gauteng Provincial Government Customer Departments / Institutions. The goods / services are therefore required by the Customer Department / Institution, as indicated on RFP 01.

CONDITIONS

1. All goods or services purchased will be subject to General Conditions of Contract.
2. All prices quoted must be inclusive of Value Added Tax (VAT) and prices must be quoted in South African currency.
3. In case of goods all delivery and/or railage charges must be included in the offered price(s) and goods must be delivered to the address indicated on RFP 01.
4. All delivery lead times submitted must be firm.
5. Firm prices, being prices deemed to be fixed prices, which are only subject to the following statutory changes, namely VAT and any levy related to customs and excise, are preferred to variable prices.
6. In case of bidding a price(s) subject to variance, please supply full details by completing RFP 07.
7. If an estimated quantity is called for, no guarantee is given that this or any quantity will be bought by GPG.

PRICING

| ITEM NO | ESTIMATED QUANTITY | DESCRIPTION | UNIT OF MEASURE |
|---------|--------------------|-------------|-----------------|
| | | | |

| | | |
|--|---|--|
| PRICE PER UNIT OF MEASURE (including VAT)* | R | (Refer to par. 16 of RFP 03 and condition no 3 on RFP 06 page 1 of 3: in case of goods, all delivery and/or railage costs must be included in the offered price*.) |
| How many days from the Closing Date of the Bid is the offer valid? | | |

DELIVERY

| | |
|--|-----|
| TIME REQUIRED FOR DELIVERY AFTER RECEIPT OF ORDER (PLEASE STATE IN HOURS, DAYS, WEEKS OR MONTHS) | |
| Is the delivery period(s) per item firm? (Also see condition no 4 on RFP 06 page 1 of 3) | Y N |

MODEL, SPECIFICATION AND GUARANTEE

| | | | |
|-----------------|---|---|---|
| Brand and model | Is the offer strictly to specification? | Y | N |
|-----------------|---|---|---|

(For specification – see RFP 05)

If not to specification, state deviation

| |
|--|
| |
| |
| |
| |

GAURANTEE



Provincial Supply Chain Management

Price Schedule - Goods

Page 2 of 3

Are the goods/services guaranteed?

Y

N

(Also see warranty requirements as required in terms of par. 15 of the General Conditions of Contract)

If a specific guarantee is required in the RFP specification, are the goods/services ACCORDINGLY guaranteed?

Y

N

(For specification – see RPF 05)

MANUFACTURING (in case of goods)

| | | | | |
|--|---|---|---|--|
| Country of manufacture | Is the bidder the manufacturer of the item? | Y | N | |
| If not, is a Letter of Supply from the relevant manufacturer(s) attached, confirming firm supply arrangement(s)? | Y | N | | |

The said manufacturer/company/supplier must confirm that it has familiarised itself with:

- o the item description;
- o specifications; and
- o bid conditions.

If the bid consist of more than one item, it should be clearly indicated in respect of which item(s) the supportive letter has been issued. Failure to submit the Letter of Supply document may invalidate the bid.

Is the bidder the accredited representative in the RSA of the manufacturer of the goods offered? (Also see RPF 01 – page 2)

Y

N

Indicate names and addresses of the factories where the goods will be manufactured/inspected?

| SITE INFORMATION | | | |
|------------------|--|------------------|--|
| Factory Name | | Physical Address | |
| Factory Name | | Physical Address | |
| Factory Name | | Physical Address | |

PHYSICAL ASSETS / EQUIPMENT (in case of goods)

What is the address in the RSA where the goods on bid, of the type offered by you may be inspected, preferably under working conditions?

What is the approximate value of spares carried in stock in the RSA for this particular make and model of equipment?

R



Provincial Supply Chain Management

Price Schedule - Goods

Page 3 of 3

At what address(es) is/are these spares held in stock?

What facilities exist for servicing / maintenance of the equipment offered?

Where (at which address) are these servicing / maintenance facilities available?

IMPORTATION (in case of goods)

NOTE: The under mentioned information must be furnished where the offered price may be affected by the fluctuation in the Rate of Exchange.

(Also ensure that RFP 07 (Non-Firm Prices per item) is completed.)

| | | | |
|--------------------------------------|----------------------------|----------------------------|--|
| Is a special import permit required? | <input type="checkbox"/> Y | <input type="checkbox"/> N | If only a general import permit is required, indicate your import permit number. |
|--------------------------------------|----------------------------|----------------------------|--|

| | |
|---|----------|
| Please state the amount in foreign currency in respect of each item, which will be remitted abroad. | R |
|---|----------|

| | | | | |
|---|--|--|--|--|
| At what stage will the amount be remitted abroad? | | | | |
|---|--|--|--|--|

| | | | | |
|---|--|--|--|--|
| State the Rate(s) of Exchange applied in the conversion of this amount to South African currency. | | | | |
|---|--|--|--|--|

| | | | |
|----------------------------|--|-------------|--|
| SIGNATURE OF BIDDER | | DATE | |
|----------------------------|--|-------------|--|

|  GAUTENG PROVINCE REPUBLIC OF SOUTH AFRICA | | | | | | | | | | | | | | | | | | | GT/GDH/002/2024 :SUPPLY AND DELIVERY OF REHABILITATION CONSUMABLES FOR A PERIOD OF THREE YEARS | | | | | | | | | | | | | | | | | | | | |
|--|--------|---|--|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|--|--|------------------------|--|----------|----------------------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGIN G/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | PROVIDE A COST BREAKDOWN AS A PERCENTAGE OF TOTAL UNIT PRICE | | IS PRODUCT IMPORTED OR | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED | | | | | | | | | | | | | | | | | | |
| A: - PHYSIOTHERAPY | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ITEM 1 Therapeutic Hand Putty | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ITEM-1.1 | 158919 | Therapeutic hand Putty Soft resistance • Used for therapeutic exercises for the hand or finger/s • Offers soft resistance for strengthening of muscles of the finger/s or hand • Non-Toxic • Non-Oily Dimensions: • 200g quantities Packaging requirements: • Packaged in re-sealable airtight durable containers Estimated usage: 0 | Yes, representative sample for items 1.1 and 1.2 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ITEM-1.2 | 177796 | Therapeutic hand Putty Soft resistance • Used for therapeutic exercises for the hand or finger/s • Offers soft resistance for strengthening of muscles of the finger/s or hand • Non-Toxic • Non-Oily Dimensions: • 500g quantities Packaging requirements: • Packaged in re-sealable airtight durable containers Estimated usage: 22 | Yes, representative sample for items 1.1 and 1.2 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ITEM-1.3 | 159010 | Therapeutic hand putty Medium resistance • Used for therapeutic exercises for the hand or finger/s • Offers medium resistance for strengthening of muscles of the finger/s or hand • Non-toxic • Non-oily Dimensions: 500g quantities Packaging requirements: Packaged in re-sealable airtight durable containers Estimated usage: 7 | Yes, representative sample for items 1.3 and 1.4 | | Each | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-1.4 | 177797 | Therapeutic hand putty Medium resistance <ul style="list-style-type: none"> Used for therapeutic exercises for the hand or finger/s Offers medium resistance for strengthening of muscles of the finger/s or hand Non-toxic Non-oily Dimensions: 200g quantities Packaging requirements: Packaged in re-sealable airtight durable containers Estimated usage: 1 | Yes, representative sample for items 1.3 and 1.4 | | Each | | | | | | | | | | | | | | |
| ITEM-1.5 | 177798 | Therapeutic hand putty Firm resistance <ul style="list-style-type: none"> Used for therapeutic exercises for the hand or finger/s Offers firm resistance for strengthening of muscles of the finger/s or hand Non-toxic Non-oily Dimensions: 500g quantities Packaging requirements: Packaged in re-sealable airtight durable containers Estimated usage: 1 | Yes, representative sample for items 1.5 and 1.6 | | Each | | | | | | | | | | | | | | |
| ITEM-1.6 | 159011 | Therapeutic hand putty Firm resistance <ul style="list-style-type: none"> Used for therapeutic exercises for the hand or finger/s Offers firm resistance for strengthening of muscles of the finger/s or hand Non-toxic Non-oily Dimensions: 200g quantities Packaging requirements: Packaged in re-sealable airtight durable containers Estimated usage: 1 | Yes, representative sample for items 1.5 and 1.6 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|--|--|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 2 | | Resistive Exercise Band | | | | | | | | | | | | | | | | | |
| ITEM-2.1 | 159012 | <p>Resistive Exercise Band – medium resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers medium resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 25m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 16</p> | Yes, representative sample for items 2.1 and 2.2 | | Per Roll | | | | | | | | | | | | | | |
| ITEM-2.2 | 177807 | <p>Resistive Exercise Band – medium resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers medium resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 50m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 2.1 and 2.2 | | Per Roll | | | | | | | | | | | | | | |
| ITEM-2.3 | 159013 | <p>Resistive Exercise Band – heavy resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers heavy resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 25m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 2.3 and 2.4 | | Per Roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-2.4 | 178716 | <p>Resistive Exercise Band – heavy resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers heavy resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 50m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 2.3 and 2.4 | | Per Roll | | | | | | | | | | | | | | |
| ITEM-2.5 | 159014 | <p>Resistive Exercise Band – extra heavy resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers extra heavy resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 25m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 18</p> | Yes, representative sample for items 2.5 and 2.6 | | Per Roll | | | | | | | | | | | | | | |
| ITEM-2.6 | 178717 | <p>Resistive Exercise Band – extra heavy resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers extra heavy resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 50m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 29</p> | Yes, representative sample for items 2.5 and 2.6 | | Per Roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-2.7 | 159016 | <p>Resistive Exercise Band – double extra heavy resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers double extra heavy resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 25m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 7</p> | Yes | | Per Roll | | | | | | | | | | | | | | |
| ITEM-2.8 | 113421 | <p>Resistive Exercise Band – light resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers light resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 25m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 2.8 and 2.9 | | Per Roll | | | | | | | | | | | | | | |
| ITEM-2.9 | 156602 | <p>Resistive Exercise Band – light resistance</p> <ul style="list-style-type: none"> • Resistive exercise band used for therapeutic strengthening of muscles • Offers light resistance to stretch • Made of elastomer • No odour <p>Dimensions: Width = 10-15cm, Length of roll = 50m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 7</p> | Yes, representative sample for items 2.8 and 2.9 | | Per Roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|---|-----------------|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 3 | | Exercise Band/Tubing Handles | | | | | | | | | | | | | | | | | |
| ITEM-3.1 | 159480 | Exercise Band/Tubing Handles • Handle for exercise band/tubing • Easy grip handle • Exercise band/tubing can be easily attached to handle Dimensions: To fit size of adult hand Packaging requirements: One Pair Estimated usage: 4 | Yes | | Pair | | | | | | | | | | | | | | |
| ITEM-3.2 | 174316 | Medium Resistant Exercise Tubing • Resistive exercise tubing used for therapeutic strengthening of muscles • Offers medium resistance to stretch Dimensions: Length of Roll: 25m Packaging requirements: Single Box Estimated usage: 1 | Yes | | Box | | | | | | | | | | | | | | |
| ITEM-3.3 | 174312 | Heavy Resistant Exercise Tubing • Resistive exercise tubing used for therapeutic strengthening of muscles • Offers heavy resistance to stretch Dimensions: Length of Roll: 25m Packaging requirements: Single Box Estimated usage: 2 | Yes | | Box | | | | | | | | | | | | | | |
| ITEM-3.4 | 174406 | Extra Heavy Resistant Exercise Tubing • Resistive exercise tubing used for therapeutic strengthening of muscles • Offers extra heavy resistance to stretch Dimensions: Length of Roll: 25m Packaging requirements: Single Box Estimated usage: 1 | Yes | | Box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|-----------|---|-----------------|---------------|-------------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 4 | Strapping | | | | | | | | | | | | | | | | | | |
| ITEM-4.1 | 159481 | <p>Strapping – underwrap</p> <ul style="list-style-type: none"> Used in conjunction with strapping Made of non-woven, breathable (porous) material which allows some stretch Adhesive-should be able to remain on the skin for up to one week Hypoallergenic Printed grid on backing sheet <p>Dimensions: Width = 5cm, Length of roll ±10m</p> <p>Packaging requirements: Single roll in a Box</p> <p>Estimated usage: 142</p> | Yes | | Per Roll in a box | | | | | | | | | | | | | | |
| ITEM-4.2 | 159482 | <p>Strapping – rigid tape</p> <ul style="list-style-type: none"> Used in conjunction with underwrap Used for therapeutic strapping of muscles and or joints Aggressively adhesive rigid strapping tape Non-stretch Made of strong rayon material, spread with zinc oxide material Porous Should be able to remain on the skin for up to one week Flesh coloured Hypoallergenic <p>Dimensions: Width 3.8cm, Length of roll ±13.7m</p> <p>Packaging requirements: Single roll</p> <p>Estimated usage: 10</p> | Yes | | Per Roll in a box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|-------------|---|--|---------------|-------------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-4.3 | 159484 | <p>Strapping – sports tape</p> <ul style="list-style-type: none"> Used for therapeutic strapping of muscles and or joints especially ankles Aggressively adhesive rigid strapping tape Non-stretch Made of strong rayon material, spread with zinc oxide material Porous Remains in place during vigorous activities Hypoallergenic Hand-tearable <p>Dimensions: Width 3.8cm, Length of roll ±13.7m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 40</p> | Yes | | Per Roll in a box | | | | | | | | | | | | | | |
| ITEM 5 | Kinesiotape | <p>Kinesiotape</p> <p>Kinesiotape, 5cm wide</p> <ul style="list-style-type: none"> Therapeutic elastic adhesive strapping tape latex free anti-allergy/hypoallergenic wave patterned adhesive design water resistant stretchable adhesive design that lifts skin to maintain flexibility, increase circulation, decrease pain, re-educate muscles Known as Kinesiotape or Kinesiology Tape <p>Dimensions: Width 5cm, Length 5m</p> <p>Packaging requirements: Single roll in a box</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 5.1 and 5.2 | | Per Roll in a box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------------------------------|--------|---|--|---------------|-------------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-5.2 | 124249 | Kinesiotape, 7.5cm wide - Therapeutic elastic adhesive strapping tape - latex free - anti-allergy/hypoallergenic - wave patterned adhesive design - water resistant - stretchable - adhesive design that lifts skin to maintain flexibility, increase circulation, decrease pain, re-educate muscles - Known as Kinesiotape or Kinesiology Tape Dimensions: Width 7.5cm, Length 5m Packaging requirements: Single roll in a box Estimated usage: 1 | Yes, representative sample for items 5.1 and 5.2 | | Per Roll in a box | | | | | | | | | | | | | | |
| ITEM 6 Dynamic tape | | | | | | | | | | | | | | | | | | | |
| ITEM-6.1 | New | Dynamic tape, 5cm • Extremely soft, highly conformable, breathable and water resistant • Allows full range of movement without limitation • Four way stretch Dimensions: 5cm x 5m Packaging requirements: Single roll in a box Estimated usage: 1 | Yes, representative sample for items 6.1 and 6.2 | | per box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-6.2 | New | Dynamic tape, 7.5cm <ul style="list-style-type: none"> Extremely soft, highly conformable, breathable and water resistant Allows full range of movement without limitation Four way stretch Dimensions: 7.5cm x 5m Packaging requirements: Single roll in a box Estimated usage: 1 | Yes, representative sample for items 6.1 and 6.2 | each | | | | | | | | | | | | | | | |
| ITEM 7 | Compression bandage short stretch | | | | | | | | | | | | | | | | | | |
| ITEM-7.1 | New | Compression bandage short stretch, 6cm <ul style="list-style-type: none"> Short stretch textile elastic Permeable, latex free, 100% cotton Flexible, high wearing comfort and reusable, can be used for management of lymphedema Dimensions: 6cm x 5m Packaging requirements: Single roll in a box Estimated usage: 1 | Yes, representative sample for items 7.1 to 7.4 | per box | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|--|---|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-7.2 | New | Compression bandage short stretch, 8cm • Short stretch textile elastic • Permeable, latex free, 100% cotton • Flexible, high wearing comfort and reusable, can be used for management of lymphedema Dimensions: 8cm x 5m Packaging requirements: Single roll in a box Estimated usage: 1 | Yes, representative sample for items 7.1 to 7.4 | | per box | | | | | | | | | | | | | | |
| ITEM-7.3 | New | Compression bandage short stretch, 10cm • Short stretch textile elastic • Permeable, latex free, 100% cotton • Flexible, high wearing comfort and reusable, can be used for management of lymphedema Dimensions: 10cm x 5m Packaging requirements: Single roll in a box Estimated usage: 1 | Yes, representative sample for items 7.1 to 7.4 | | per box | | | | | | | | | | | | | | |
| ITEM-7.4 | New | Compression bandage short stretch, 12cm • Short stretch textile elastic • Permeable, latex free, 100% cotton • Flexible, high wearing comfort and reusable, can be used for management of lymphedema Dimensions: 12cm x 5m Packaging requirements: Single roll in a box Estimated usage: 1 | Yes, representative sample for items 7.1 to 7.4 | | per box | | | | | | | | | | | | | | |
| ITEM-7.5 | New | Open cell foam bandage for padding • Bandage for padding under compression bandaging • 100% polyurethane, latex free • Air and moisture permeable. Non-sliding and tear resistant. Re-usable with same patient Dimensions: 10cm x 2.5m x 0.4cm Packaging requirements: Single roll in a box Estimated usage: 1 | Yes | | per box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|---|---|---------------|--------------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 8 | | Elastic Tubular Support Bandage | | | | | | | | | | | | | | | | | |
| ITEM-8.1 | New | <p>Elastic tubular support bandage, 4.2cm</p> <ul style="list-style-type: none"> • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <p>Dimensions: 4.2cm x 10m Packaging requirements: Single roll in a box Estimated usage: 1</p> | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |
| ITEM-8.2 | New | <p>Elastic tubular support bandage, 6.25cm</p> <ul style="list-style-type: none"> • Elasticated tubular bandage that provides firm and comfortable support, 100% • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <p>Dimensions: 6.25cm x 10m Packaging requirements: Single roll in a box Estimated usage: 1</p> | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |
| ITEM-8.3 | New | <p>Elastic tubular support bandage, 6.75cm</p> <ul style="list-style-type: none"> • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <p>Dimensions: 6.75cm x 10m Packaging requirements: Single roll in a box Estimated usage: 1</p> | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |
| ITEM-8.4 | New | <p>Elastic tubular support bandage, 7.5cm</p> <ul style="list-style-type: none"> • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <p>Dimensions: 7.5cm x 10m Packaging requirements: Single roll in a box Estimated usage: 1</p> | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |
| ITEM-8.5 | New | <p>Elastic tubular support bandage, 8.75cm</p> <ul style="list-style-type: none"> • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. <p>Dimensions: 8.75cm x 10m Packaging requirements: Single roll in a box Estimated usage: 1</p> | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|---|---|---------------|--------------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-8.6 | New | Elastic tubular support bandage, 10cm • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. Dimensions: 10cm x 10m Packaging requirements: Single roll in a box Estimated usage: 0 | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |
| ITEM-8.7 | New | Elastic tubular support bandage, 12cm • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. Dimensions: 12cm x 10m Packaging requirements: Single roll in a box Estimated usage: 0 | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |
| ITEM-8.8 | New | Elastic tubular support bandage, 14cm • Elasticated tubular bandage that provides firm and comfortable support, 100% cotton • Reduces persistent swelling • Easy to apply, re-usable with same patient, easy to wash and latex free. Dimensions: 14cm x 10m Packaging requirements: Single roll in a box Estimated usage: 0 | Yes, representative sample for items 8.1 to 8.8 | | per box / per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------|--------|---|--|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 9 | | Paraffin Wax Refill Paraffin wax refill, beads For therapeutic wax bath • Wax beads • Scent free • Paraffin wax mixed with paraffin oil Dimensions: 500g Packaging requirements: Sealed packaging Estimated usage: 7 | Yes, representative sample for items 9.1 and 9.2 | | Pkt | | | | | | | | | | | | | | |
| ITEM-9.1 | 159541 | Paraffin wax refill, beads For therapeutic wax bath • Wax beads • Scent free • Paraffin wax mixed with paraffin oil Dimensions: 500g Packaging requirements: Sealed packaging Estimated usage: 7 | Yes, representative sample for items 9.1 and 9.2 | | Pkt | | | | | | | | | | | | | | |
| ITEM-9.2 | 113522 | Paraffin wax refill, beads For therapeutic wax bath • Wax beads • Scent free • Paraffin wax mixed with paraffin oil Dimensions: 1kg Packaging requirements: Sealed packaging Estimated usage: 7 | Yes, representative sample for items 9.1 and 9.2 | | Pkt | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|---------------------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-9.3 | 159541 | Paraffin wax refill, block For therapeutic wax bath • Wax block • Scent free • Paraffin wax mixed with paraffin oil Dimensions: 500g Packaging requirements: Sealed packaging Estimated usage: 7 | Yes, representative sample for items 9.3 and 9.4 | | Pkt | | | | | | | | | | | | | | |
| ITEM-9.4 | 136424 | Paraffin wax refill, block For therapeutic wax bath • Wax Block • Scent free • Paraffin wax mixed with paraffin oil Dimensions: 1kg Packaging requirements: Sealed packaging Estimated usage: 7 | Yes, representative sample for items 9.3 and 9.4 | | Pkt | | | | | | | | | | | | | | |
| ITEM 10 Arnica Oil | | | | | | | | | | | | | | | | | | | |
| ITEM-10.1 | 159542 | Arnica Oil – 250 ml • Oil containing Arnica concentrate • Contain essential oils which promotes muscular relaxation • Easily removable Dimensions: 250 ml Packaging: Per Bottle Estimated usage: 6 | Yes, representative sample for items 10.1 to 10.3 | | Per bottle | | | | | | | | | | | | | | |
| ITEM-10.2 | 159543 | Arnica Oil – 500 ml • Massage Oil containing Arnica concentrate • Contain essential oils which promotes muscular relaxation • Easily removable Dimensions: 500 ml Packaging: Per Bottle Estimated usage: 11 | Yes, representative sample for items 10.1 to 10.3 | | Per bottle | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-10.3 | 159544 | Arnica Oil – 1L <ul style="list-style-type: none"> • Massage Oil containing Arnica concentrate • Contain essential oils which promotes muscular relaxation • Easily removable Dimensions: 1L Packaging: Per Bottle Estimated usage: 13 | Yes, representative sample for items 10.1 to 10.3 | | Per bottle | | | | | | | | | | | | | | |
| ITEM 11 | | | | | | | | | | | | | | | | | | | |
| ITEM-11.1 | 159546 | Acupuncture needles – 13 mm <ul style="list-style-type: none"> • Copper handled needles with guide tubes • Disposable: for single use only Dimensions: Length = 13mm, Gauge = 0.22mm Packaging requirements: Sterile individually packed, 100 needles per box Estimated usage: 31 | Yes | | Per box | | | | | | | | | | | | | | |
| ITEM-11.2 | 159611 | Acupuncture needles – 25mm <ul style="list-style-type: none"> • Copper handled needles with guide tubes • Disposable: for single use only Dimensions: Length = 25mm, Gauge = 0.25mm Packaging requirements: Sterile individually packed, 100 needles per box Estimated usage: 38 | Yes | | Per box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-11.3 | 162854 | Acupuncture needles – 40 mm • Copper handled needles with guide tubes • Disposable: for single use only Dimensions: Length = 40mm, Gauge = 0.30mm Packaging requirements: Sterile individually packed, 100 needles per box Estimated usage: 19 | Yes | | Per box | | | | | | | | | | | | | | |
| ITEM-11.4 | 159612 | Acupuncture needles – 50 mm • Copper handled needles with guide tubes • Disposable: for single use only Dimensions: Length = 50mm, Gauge = 0.30mm Packaging requirements: Sterile individually packed, 100 needles per box Estimated usage: 28 | Yes | | Per box | | | | | | | | | | | | | | |
| ITEM-11.5 | 136425 | Acupuncture needles – 30 mm • Copper handled needles with guide tubes • Disposable: for single use only Dimensions: Length = 0,25mm, Gauge = 0.30-35mm Packaging requirements: Sterile individually packed, 100 needles per box Estimated usage: 39 | Yes | | Per box | | | | | | | | | | | | | | |
| ITEM-11.6 | 159613 | Acupuncture needles – 75 mm • Copper handled needles with guide tubes • Disposable: for single use only Dimensions: Length = 75mm, Gauge = 0.30-35mm Packaging requirements: Sterile individually packed, 100 needles per box Estimated usage: 39 | Yes | | Per box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|-----------------|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-11.7 | 159614 | Acupuncture needles –100 mm • Copper handled needles with guide tubes • Disposable: for single use only Dimensions: Length = 100mm, Gauge = 0.30-35mm Packaging requirements: Sterile individually packed, 100 needles per box Estimated usage: 28 | Yes | | Per box | | | | | | | | | | | | | | |
| ITEM 12 | | Flutter / Broncho-vibe | | | | | | | | | | | | | | | | | |
| ITEM-12.1 | 160068 | Flutter / Broncho-vibe • A device which generates controlled oscillating positive pressure during expiration • Aids in mobilizing secretions • Hard, plastic cylindrical mouth piece around a circular core of hardened material • Contains a high density stainless steel ball • Has a perforated protective cover • Washable Dimensions: Length ± 15cm, Diameter ± 3cm Packaging requirements: Each in a sealed pack Estimated usage: 92 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 13 | | Luft Respiratory Device | | | | | | | | | | | | | | | | | |
| ITEM-13.1 | New | Luft Respiratory Device • Resistance training for the lungs • Assists with breathing • Hand held device • Re-usable with same patient Dimensions: None Packaging requirements: Each in box Estimated usage: 1 | Yes | | each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|-----------------|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 14 | | Positive Expiratory Pressure | | | | | | | | | | | | | | | | | |
| ITEM-14.1 | 137067 | Positive Expiratory Pressure (PEP) - Made of transparent hard plastic - Comes with removable mouth-piece - Small volume nebulizer - Collapsible reservoir tubing which must be able to attach to the nebulizer - It must be delivered with all accessories - It must be autoclavable Packaging: Each in a sealed plastic pack Estimated usage : 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM15 | | Respi-flow Respiratory Exerciser | | | | | | | | | | | | | | | | | |
| ITEM-15.1 | 137068 | Respi-flow Respiratory Exerciser - Must have 3 chambers - 3 balls - Disposed after a single patient use - Made of transparent plastic with removable re-usable mouth piece - Flow maximum inspiration - Help to develop & improve respiratory fitness Packaging: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 16 | | Transparent cling wrap | | | | | | | | | | | | | | | | | |
| ITEM-16.1 | 160125 | Transparent cling wrap • Transparent • Self attaching • Non breathable • Used for wrapping the body or parts thereof during dermatology treatments to facilitate uptake of preparations into the skin Dimensions: Length ± 40m, Width ± 30cm Diameter: ± 3cm Packaging requirements: Sealed packaging, Single item Estimated usage: 177 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 17 | | Rectangular sponge sleeves | | | | | | | | | | | | | | | | | |
| ITEM-17.1 | 160126 | Rectangular sponge sleeves, medium For carbon rubber electrodes of interferential machine • Durable, high quality sponge 3 Layers thick Dimensions: To fit rubber electrode of 60 x 80 mm Packaging requirements: Sets of 4 Estimated usage: 30 | Yes, representative sample for items 17.1 and 17.2 | | Per set | | | | | | | | | | | | | | |
| ITEM-17.2 | 160129 | Rectangular sponge sleeves, large For carbon rubber electrodes of interferential machine • Durable, high quality sponge 3 Layers thick Dimensions: To fit rubber electrode of 80 x 120 mm Packaging requirements: Sets of 4 Estimated usage: 1 | Yes, representative sample for items 17.1 and 17.2 | | Per set | | | | | | | | | | | | | | |
| ITEM 18 | | Electrodes | | | | | | | | | | | | | | | | | |
| ITEM-18.1 | 160131 | Reusable round/oval gel electrodes – 50mm • Solid gel with cover • Suitable for multiple applications for 10-15 days • Self adhering • Low profile with pigtail pin connection Dimensions: = ± 50mm Packaging requirements: Set of 4 Estimated usage: 28 | Yes | | Per set | | | | | | | | | | | | | | |
| ITEM-18.2 | 160428 | Reusable rectangular gel electrodes – 50 x 90mm • Solid gel with cover • Suitable for multiple applications for 10-15 days • Self adhering • Low profile with pigtail pin connection Dimensions: = 50mm x ± 90mm Packaging requirements: Set of 4 Estimated usage: 56 | Yes | | Per set | | | | | | | | | | | | | | |
| ITEM-18.3 | 160431 | Carbon electrode for multiple use – 50 x 50mm • Carbon electrodes for use with TENS machines • Suitable for multiple use Dimensions: 50 x 50mm Packaging requirements: Set of 4 Estimated usage: 27 | Yes | | Per set | | | | | | | | | | | | | | |
| ITEM 19 | | Protective eyewear | | | | | | | | | | | | | | | | | |
| ITEM-19.1 | 160434 | Protective Eyewear for Therapeutic Laser machine of ± 830 nm wavelength • Sturdy ergonomic design • Suitable for emissions of wavelength of ± 830 nm Dimensions: Suitable for most face types Packaging requirements: Single pair Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|-----------------|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-19.2 | 178075 | <p>Protective Eyewear for UV Phototherapy</p> <ul style="list-style-type: none"> • For use by patients when undergoing phototherapy • Blocks 100% UVA / UVB • Totally opaque to UV Rays <p>Dimensions: Rubber eye pieces with adjustable nose bridge for comfortable fit</p> <p>Packaging requirements: Single pair</p> <p>Estimated usage: 1</p> | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-19.3 | 178076 | <p>Protective Eyewear for UV Phototherapy</p> <ul style="list-style-type: none"> • For use by patients after UVA / UVB light treatment • Blocks 100% UVA / UVB • Totally opaque to UV rays <p>Dimensions: UV blocking wraparound glasses</p> <p>Packaging requirements: Single pair</p> <p>Estimated usage: 1</p> | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERIAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 20 | Positioner: Preformed Resting Foot Drop Splint | | | | | | | | | | | | | | | | | | |
| ITEM-20.1 | New | Positioner: Preformed Resting Foot Drop Splint - Small Left <ul style="list-style-type: none"> Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles Helps protect heel area by preventing excessive pressure and lessening friction Velcro closure system Washable/ autoclavable Dimensions: Shoe sizes 3-6 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 20.1 to 20.6 | | Each | | | | | | | | | | | | | | |
| ITEM-20.2 | New | Positioner: Preformed Resting Foot Drop Splint - Small Right <ul style="list-style-type: none"> Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles Helps protect heel area by preventing excessive pressure and lessening friction Velcro closure system Washable/ autoclavable Dimensions: Shoe sizes 3-6 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 20.1 to 20.6 | | Each | | | | | | | | | | | | | | |
| ITEM-20.3 | New | Positioner: Preformed Resting Foot Drop Splint - Medium Left <ul style="list-style-type: none"> Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles Helps protect heel area by preventing excessive pressure and lessening friction Velcro closure system Washable/ autoclavable Dimensions: Shoe sizes ± 7-9 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 20.1 to 20.6 | | Each | | | | | | | | | | | | | | |
| ITEM-20.4 | New | Positioner: Preformed Resting Foot Drop Splint - Medium Right <ul style="list-style-type: none"> Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles Helps protect heel area by preventing excessive pressure and lessening friction Velcro closure system Washable/ autoclavable Dimensions: Shoe sizes ± 7-9 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 20.1 to 20.6 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-20.5 | New | Positioner: Preformed Resting Foot Drop Splint - Large Left • Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle • Not for ambulation; used on bed-ridden clients • Used to prevent muscle and tendon shortening of the calf muscles • Helps protect heel area by preventing excessive pressure and lessening friction • Velcro closure system • Washable/ autoclavable Dimensions: Shoe sizes ± 10-12 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 20.1 to 20.6 | | Each | | | | | | | | | | | | | | |
| ITEM-20.6 | New | Positioner: Preformed Resting Foot Drop Splint - Large Right • Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle • Not for ambulation; used on bed-ridden clients • Used to prevent muscle and tendon shortening of the calf muscles • Helps protect heel area by preventing excessive pressure and lessening friction • Velcro closure system • Washable/ autoclavable Dimensions: Shoe sizes ± 10-12 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 20.1 to 20.6 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 21 | | Paediatric Positioner: Preformed Resting Foot Drop Splint | | | | | | | | | | | | | | | | | |
| ITEM-21.1 | New | Paediatric Positioner: Preformed Resting Foot Drop Splint - Small Left <ul style="list-style-type: none"> • Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle • Not for ambulation; used on bed-ridden clients • Used to prevent muscle and tendon shortening of the calf muscles • Helps protect heel area by preventing excessive pressure and lessening friction • Velcro closure system • Washable/ autoclavable Dimensions: Shoe sizes ± 4-6 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 21.1 to 21.6 | | Each | | | | | | | | | | | | | | |
| ITEM-21.2 | New | Paediatric Positioner: Preformed Resting Foot Drop Splint - Small Right <ul style="list-style-type: none"> • Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle • Not for ambulation; used on bed-ridden clients • Used to prevent muscle and tendon shortening of the calf muscles • Helps protect heel area by preventing excessive pressure and lessening friction • Velcro closure system • Washable/ autoclavable Dimensions: Shoe sizes ± 4-6 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 21.1 to 21.6 | | Each | | | | | | | | | | | | | | |
| ITEM-21.3 | New | Paediatric Positioner: Preformed Resting Foot Drop Splint - Medium Left <ul style="list-style-type: none"> • Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle • Not for ambulation; used on bed-ridden clients • Used to prevent muscle and tendon shortening of the calf muscles • Helps protect heel area by preventing excessive pressure and lessening friction • Velcro closure system • Washable/ autoclavable Dimensions: Shoe sizes 7-9 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 21.1 to 21.6 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-21.4 | New | Paediatric Positioner: Preformed Resting Foot Drop Splint - Medium Right <ul style="list-style-type: none"> Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles Helps protect heel area by preventing excessive pressure and lessening friction Velcro closure system Washable/ autoclavable Dimensions: Shoe sizes 7-9 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 21.1 to 21.6 | | Each | | | | | | | | | | | | | | |
| ITEM-21.5 | New | Paediatric Positioner: Preformed Resting Foot Drop Splint - Large Left <ul style="list-style-type: none"> Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles Helps protect heel area by preventing excessive pressure and lessening friction Velcro closure system Washable/ autoclavable Dimensions: Shoe sizes ± 10-12 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 21.1 to 21.6 | | Each | | | | | | | | | | | | | | |
| ITEM-21.6 | New | Paediatric Positioner: Preformed Resting Foot Drop Splint - Large Right <ul style="list-style-type: none"> Pre-padded, anatomically contoured, polyform splint providing 0° flexion at the ankle Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles Helps protect heel area by preventing excessive pressure and lessening friction Velcro closure system Washable/ autoclavable Dimensions: Shoe sizes ± 10-12 Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 21.1 to 21.6 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|--|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 22 Positioner: Foam Resting Foot Drop Splint | | | | | | | | | | | | | | | | | | | |
| ITEM 22.1 | New | Positioner: Foam Resting Foot drop splint - right foot <ul style="list-style-type: none"> Foam or foam coated with washable vinyl orthosis Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles and to reduce pressure on the heel Anatomically contoured right Side opening for cross ventilation and visual contact Secure fastening system Washable or autoclavable Dimensions: ± 40cm long from heel to calf, for right foot Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 22.1 and 22.2 | | Each | | | | | | | | | | | | | | |
| ITEM-22.2 | New | Positioner: Foam Resting Foot drop splint - left foot <ul style="list-style-type: none"> Foam or foam coated with washable vinyl orthosis Not for ambulation; used on bed-ridden clients Used to prevent muscle and tendon shortening of the calf muscles and to reduce pressure on the heel Anatomically contoured left Side opening for cross ventilation and visual contact Secure fastening system Washable or autoclavable Dimensions: ± 40cm long from heel to calf, for left foot Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 22.1 and 22.2 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 23 | | Positioner: Angle adjustable resting foot drop splint | | | | | | | | | | | | | | | | | |
| ITEM-23.1 | 161026 | Positioner: Angle adjustable resting foot drop splint – standard <ul style="list-style-type: none"> • Positions from 0°-45° plantar flexion. • Not for ambulation; used on bed-ridden clients • Used to prevent muscle and tendon shortening of the calf muscles and to reduce pressure on the heel • Washable • Velcro closure system Dimensions: • Calf circumference up to 36.8cm, • Can fit either foot Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 23.1 and 23.2 | | Each | | | | | | | | | | | | | | |
| ITEM-23.2 | 161027 | Positioner: Angle adjustable resting foot drop splint – extra wide <ul style="list-style-type: none"> • Positions foot from 0°-45° plantar flexion. • Not for ambulation; used on bed-ridden clients • Used to prevent muscle and tendon shortening of the calf muscles and to reduce pressure on the heel • Washable • Velcro closure system Dimensions: • For calf circumference from 36.9cm and up, • Can fit either foot Packaging requirements: Single item Estimated usage: 1 | Yes, representative sample for items 23.1 and 23.2 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 24 | | Hot pack moist heat hydro | | | | | | | | | | | | | | | | | |
| ITEM-24.1 | 172223 | Hot pack moist heat hydro, half size 13x30cm Hotpack moist head hydrocollator; half size, 13cm width; 30cm length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 24.1 to 24.7 | | Each | | | | | | | | | | | | | | |
| ITEM-24.2 | 172221 | Hot pack moist heat hydro, oversize 38x61cm Hotpack moist head hydrocollator; oversize 38cm width; 61cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 24.1 to 24.7 | | Each | | | | | | | | | | | | | | |
| ITEM-24.3 | 172236 | Hot pack moist heat hydro, spinal 25 x 46cm Hotpack moist head hydrocollator; spinal 25cm width; 46cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 24.1 to 24.7 | | Each | | | | | | | | | | | | | | |
| ITEM-24.4 | 172238 | Hot pack moist heat hydro, spinal 25 x 61cm Hotpack moist head hydrocollator; spinal 25cm width; 61cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 24.1 to 24.7 | | Each | | | | | | | | | | | | | | |
| ITEM-24.5 | 172226 | Hot pack moist heat, knee-shoulder 25x50cm Hotpack moist head hydrocollator; spinal 25cm width; 50cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 24.1 to 24.7 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-24.6 | 136433 | Hot pack moist heat hydro, cervical 60cm x 15cm Hotpack moist head hydrocollator; cervical 15cm width; 60cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. To fit the contours of the neck. Apply moist heat for a minimum of 30mins. Pack are re-usable. Packaging requirements: Each Estimated usage: 0 | Yes, representative sample for items 24.1 to 24.7 | | Each | | | | | | | | | | | | | | |
| ITEM-24.7 | 172220 | Hot pack moist heat hydro, standard 25cm x 30cm Hotpack moist head hydrocollator; standard 25cm width; 30cm in length; made of new synthetic fabric specifically to fight bacteria and its accompanying odour. Apply moist heat for a minimum of 30mins. Pack are re-usable. Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 24.1 to 24.7 | | Each | | | | | | | | | | | | | | |
| ITEM 25 | | Washable hot pack covers | | | | | | | | | | | | | | | | | |
| ITEM-25.1 | 117731 | Washable hot pack covers, standard 48x69cm Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. Size: standard 48x69cm Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 25.1 to 25.5 | | Each | | | | | | | | | | | | | | |
| ITEM-25.2 | 117733 | Washable hot pack covers, oversize 61x91cm Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. Size: oversize 61x91cm Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 25.1 to 25.5 | | Each | | | | | | | | | | | | | | |
| ITEM-25.3 | 117734 | Washable hot pack covers, cervical 64x46cm Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. Size: cervical 64x46cm Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 25.1 to 25.5 | | Each | | | | | | | | | | | | | | |
| ITEM-25.4 | 117805 | Washable hot pack covers, to fit knee/shoulder hotpack 25x51cm Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. Size: knee/shoulder 25x51cm Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 25.1 to 25.5 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-25.5 | 117806 | Washable hot pack covers, to fit half size 13x30cm Multi-layered all terry hot pack covers are soft and durable, cover are not as stiff as the foam filled version and have the feel of a plush towel. Washable hot pack covers. Size: half size 13x30cm. Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 25.1 to 25.5 | | Each | | | | | | | | | | | | | | |
| ITEM 26 | | Hemi-Arm Sling | | | | | | | | | | | | | | | | | |
| ITEM-26.1 | 136434 | Hemi-Arm Sling for Paeds - Small Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. Size: small Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 26.1 to 26.3 | | Each | | | | | | | | | | | | | | |
| ITEM-26.2 | 136435 | Hemi-Arm Sling for Paeds - Medium Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. Size: medium Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 26.1 to 26.3 | | Each | | | | | | | | | | | | | | |
| ITEM-26.3 | 136436 | Hemi-Arm Sling for Paeds - Large Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. Size: large Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 26.1 to 26.3 | | Each | | | | | | | | | | | | | | |
| ITEM-26.4 | 117813 | Hemi-Arm Sling for Adult - Small Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. Size: small Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 26.4 to 26.6 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERIAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|--|---------------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-26.5 | 117814 | Hemi-Arm Sling for Adult - Medium Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. Size: medium Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 26.4 to 26.6 | | Each | | | | | | | | | | | | | | |
| ITEM-26.6 | 117815 | Hemi-Arm Sling for Adult - Large Provides effective unilateral anti-subluxation support. Positions the head of the humerus in the glenoid fossa and scapula in external rotation and enhances normal gait pattern. Dynamic elastic straps allowing active movement; made of durable washable fabric. Size: large Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 26.4 to 26.6 | | Each | | | | | | | | | | | | | | |
| ITEM 27 Sacro-iliac Joint (SIJ) Belt | | | | | | | | | | | | | | | | | | | |
| ITEM-27.1 | 117816 | Sacro-iliac Joint (SIJ) Belt - small Provides comprehensive stabilization support at sacro-iliac joint to fit just under the anti-supra-iliac spine (ASIS); Washable and Durable; Neoprene underbelt to prevent slippage of belt; Suitable adjustable material overbelt to adjust anteriorly and posteriorly. Size: Small Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 27.1 to 27.3 | | Each | | | | | | | | | | | | | | |
| ITEM-27.2 | 117817 | Sacro-iliac Joint (SIJ) Belt - medium Provides comprehensive stabilization support at sacro-iliac joint to fit just under the anti-supra-iliac spine (ASIS); Washable and Durable; Neoprene underbelt to prevent slippage of belt; Suitable adjustable material overbelt to adjust anteriorly and posteriorly. Size: Medium Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 27.1 to 27.3 | | Each | | | | | | | | | | | | | | |
| ITEM-27.3 | 117819 | Sacro-iliac Joint (SIJ) Belt - large Provides comprehensive stabilization support at sacro-iliac joint to fit just under the anti-supra-iliac spine (ASIS); Washable and Durable; Neoprene underbelt to prevent slippage of belt; Suitable adjustable material overbelt to adjust anteriorly and posteriorly. Size: Large Packaging requirements: Each Estimated usage: 1 | Yes, representative sample for items 27.1 to 27.3 | | Each | | | | | | | | | | | | | | |
| ITEM 28 | Probes | | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|-------------------------------------|--|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-28.1 | 136437 | Intra vaginal probe Can be used with most muscle stimulators and biofeedback equipment. Single patient use. Medical grade surgical stainless steel. Unique rectangular section resist lateral movement. Countered electrode provides comfortable stimulation and accurate biofeedback. Smooth clean shape. Long integrated cables. Push fit indicator wand with extension pieces. Dimensions: length 77mm, conductive part width 14mm with a 20mm width of the probe Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-28.2 | 136438 | Intra-anal and small intra vaginal probe When inserted, the probe allows the anal canal sphincter to close around the narrow neck of the probe. Revolutionary shape for repeatable positioning and ambulatory use. Can be used with most muscle stimulators and biofeedback equipment. Single patient use. Medical grade surgical stainless steel. Unique rectangular section resist lateral movement. Countered electrode provides comfortable stimulation. Smooth clean shape. Long integrated cables. Always sits in the same position comparable EMR readings can be taken over different sessions. Allows for comfortable muscle stimulation and accurate biofeedback. Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 29 | Pelvic floor exercise system | | | | | | | | | | | | | | | | | | |
| ITEM-29.1 | 136439 | Pelvic floor exercise system Exercise system consists of two cones and a set of 4 weight or more sizes (5g,10,20g). Single patient use. Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|---------------------------------|--------|--|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| B:- OCCUPATIONAL THERAPY | | | | | | | | | | | | | | | | | | | |
| ITEM 30 | 156745 | Thermoplastic splinting material with cohesive & adhesive properties PERFORATED thermoplastic splinting material with cohesive & adhesive properties -3.2mm • Low temperature thermoplastic material used to make splints • Micro-perforated 13% • Has the following properties: - Resistance to stretch: Moderate - Drape/Conformity: Maximum - Memory/Recovery: 100% - Rigidity: minimum - Bonding: cohesive & adhesive properties - Surface finish: Fingerprint resist • Colours: Standard (fleshy pink) Dimensions: 46cm x 61cm x 3.2mm sheet Packaging requirements: Pack of 4 sheets Estimated usage: 1 | Yes | | Pack of 4 sheet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|--|--------|--|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-30.2 | 156742 | NON PERFORATED thermoplastic splinting material with cohesive & adhesive properties - 3.2mm <ul style="list-style-type: none"> Low temperature thermoplastic material used to make splints Has the following properties: <ul style="list-style-type: none"> Resistance to stretch: Moderate Drape/Conformity: Maximum Memory/Recovery: 100% Rigidity: minimum Bonding: cohesive & adhesive properties Surface finish: Fingerprint resist Colours: Standard (fleshy pink) Dimensions: 46cm x 61cm x 3.2mm sheet Packaging requirements: Pack of 4 sheets | Yes | | Pack of 4 sheet | | | | | | | | | | | | | | |
| ITEM-30.3 | 178107 | NON PERFORATED thermoplastic splinting material with cohesive & adhesive properties - 1.6mm <ul style="list-style-type: none"> Low temperature thermoplastic material used to make splints Has the following properties: <ul style="list-style-type: none"> Resistance to stretch: Moderate Drape/Conformity: Maximum Memory/Recovery: 100% Rigidity: minimum Bonding: cohesive & adhesive properties Surface finish: Fingerprint resist Colours: Standard (fleshy pink) Dimensions: 46cm x 61cm x 1.6mm sheet Packaging requirements: Pack of 4 sheets Estimated usage: 1 | Yes | | Pack of 4 sheet | | | | | | | | | | | | | | |
| ITEM-30.3 | 156744 | NON PERFORATED thermoplastic splinting material with cohesive & adhesive properties 2mm <ul style="list-style-type: none"> Low temperature thermoplastic material used to make splints Has the following properties: <ul style="list-style-type: none"> Resistance to stretch: Moderate Drape/Conformity: Maximum Memory/Recovery: 100% Rigidity: minimum Bonding: cohesive & adhesive properties Surface finish: Fingerprint resist Colours: Standard (fleshy pink) Dimensions: 46cm x 61cm x 2mm sheet Packaging requirements: Pack of 4 sheets Estimated usage: 50 | Yes | | Pack of 4 sheet | | | | | | | | | | | | | | |
| ITEM 31 Thermoplastic splinting material with temporary bonding properties | | | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-31.1 | 156746 | <p>NON-PERFORATED thermoplastic splinting material with temporary bonding properties - 3.2 mm</p> <ul style="list-style-type: none"> • Low temperature thermoplastic material used to make splints • Non-perforated • Has the following properties: <ul style="list-style-type: none"> - Resistance to stretch: Maximum - Drape/Conformity: Moderate - Memory/Recovery: 100% - Rigidity: Maximum - Bonding: Temporary bonding qualities only - Surface finish: Fingerprint resist; transparent when heated; non-stick coating • Colour: White <p>Dimensions: 46cm x 61 cm x 3,2mm sheet</p> <p>Packaging requirements: Pack of 4 sheets</p> <p>Estimated usage: 54</p> | Yes | | Sheet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-31.2 | 178108 | NON-PERFORATED thermoplastic splinting material with temporary bonding properties - 1,6 mm • Low temperature thermoplastic material used to make splints • Non-perforated • Has the following properties: - Resistance to stretch: Maximum - Drape/Conformity: Moderate - Memory/Recovery: 100% - Rigidity: Maximum - Bonding: Temporary bonding qualities only - Surface finish: Fingerprint resist; transparent when heated; non-stick coating • Colour: White Dimensions: 46cm x 61 cm x 1.6mm sheet Packaging requirements: Pack of 4 sheets Estimated usage: 1 | Yes | | Sheet | | | | | | | | | | | | | | |
| ITEM-31.3 | 178109 | NON-PERFORATED thermoplastic splinting material with temporary bonding properties - 2,4 mm • Low temperature thermoplastic material used to make splints • Non-perforated • Has the following properties: - Resistance to stretch: Maximum - Drape/Conformity: Moderate - Memory/Recovery: 100% - Rigidity: Maximum - Bonding: Temporary bonding qualities only - Surface finish: Fingerprint resist; transparent when heated; non-stick coating • Colour: White Dimensions: 46cm x 61 cm x 2.4mm sheet Packaging requirements: Pack of 4 sheets Estimated usage: 24 | Yes | | Sheet | | | | | | | | | | | | | | |
| ITEM-31.4 | 156750 | PERFORATED thermoplastic splinting material with temporary bonding properties - 2,4mm • Low temperature thermoplastic material used to make splints • 1% perforated • Has the following properties: - Resistance to stretch: Maximum - Drape/Conformity: Moderate - Memory/Recovery: 100% - Rigidity: Maximum - Bonding: Temporary bonding qualities only - Surface finish: Fingerprint resist, transparent when heated, non-stick coating • Colour: White Dimensions: 46cm x 61 cm x 2.4 mm sheet Packaging requirements: Pack of 4 sheets Estimated usage: 1 | Yes | | Sheet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 32 | | Non-Perforated thermo-plastic splinting material -Equivalent to OMEGA MAX | | | | | | | | | | | | | | | | | |
| ITEM-32.1 | New | <p>Non-Perforated thermo-plastic splinting material - 3.2mm</p> <ul style="list-style-type: none"> • Low temperature thermoplastic material used to make splints • Has the following properties: <ul style="list-style-type: none"> - Resistance to stretch: Moderate - Drape/Conformity: Moderate to maximum - Memory/Recovery: 100% - Rigidity: Moderate to maximum - Bonding: Has a powder coating to give a slightly tacky surface that holds material in place. - When cut gives smooth self sealed edges. - Working time of 3 minutes to form and harden after softening. - Similar or equivalent to OMEGA MAX material. - Surface finish: Fingerprint resist • Colour: Beige <p>Dimensions: 45cm x 60cm x 3.2mm sheet</p> <p>Packaging requirements: Pack of 4 sheets</p> <p>Estimated usage: 1</p> | Yes | each | | | | | | | | | | | | | | | |
| ITEM-32.2 | New | <p>Non-Perforated thermo-plastic splinting material - 2.4mm</p> <ul style="list-style-type: none"> • Low temperature thermoplastic material used to make splints • Has the following properties: <ul style="list-style-type: none"> - Resistance to stretch: Moderate - Drape/Conformity: Moderate to maximum - Memory/Recovery: 100% - Rigidity: Moderate to maximum - Bonding: Has a powder coating to give a slightly tacky surface that holds material in place. - When cut gives smooth self sealed edges. - Working time of 3minutes to form and harden after softening. - Similar or equivalent to OMEGA MAX material - Surface finish: Fingerprint resist • Colour: Beige <p>Dimensions: 45cm x 60cm x 2.4mm sheet</p> <p>Packaging requirements: Pack of 4 sheets</p> <p>Estimated usage: 1</p> | Yes | each | | | | | | | | | | | | | | | |
| ITEM 33 | | Black textile thermo-plastic tape | | | | | | | | | | | | | | | | | |
| ITEM-33.1 | New | <p>Black textile thermo-plastic tape - 3m, black</p> <ul style="list-style-type: none"> • Low temperature thermoplastic material used to make splints • Has the following properties: <ul style="list-style-type: none"> - Elasticity in two directions - Made of thin, breathable fabric - Self-adhesive capacity when dry heated • Colour: Black <p>Dimensions: 3cm x 3m, 1.3 mm thickness</p> <p>Packaging requirements: Per roll in a box</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 33.1 and 33.2 | per box | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-33.2 | New | Black textile thermo-plastic tape - 6m, black • Low temperature thermoplastic material used to make splints • Has the following properties: - Elasticity in two directions - Made of thin, breathable fabric - Self-adhesive capacity when dry heated • Colour: Black Dimensions: 6cm x 6m; 1.3mm thickness Packaging requirements: Per roll in a box Estimated usage: 1 | Yes, representative sample for items 33.1 and 33.2 | | Per box | | | | | | | | | | | | | | |
| ITEM-33.3 | New | Blue textile thermo-plastic tape - 3m, blue • Low temperature thermoplastic material used to make splints • Has the following properties: - Elasticity in two directions - Made of thin, breathable fabric - Self-adhesive capacity when dry heated • Colour: Blue Dimensions: 3cm x 3m; 1.3mm thickness Packaging requirements: Per roll in box Estimated usage: 1 | Yes, representative sample for items 33.3 and 33.4 | | Per box | | | | | | | | | | | | | | |
| ITEM-33.4 | New | Blue textile thermo-plastic tape - 6m, blue • Low temperature thermoplastic material used to make splints • Has the following properties: - Elasticity in two directions - Made of thin, breathable fabric - Self-adhesive capacity when dry heated • Colour: Blue Dimensions: 6cm x 6m; 1.3mm thickness Packaging requirements: Per roll in box Estimated usage: 1 | Yes, representative sample for items 33.3 and 33.4 | | Per box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|--|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 34 Wrist Brace with closure system | | | | | | | | | | | | | | | | | | | |
| ITEM-34.1 | 156756 | Wrist Brace with closure system - EXTRA SMALL, left wrist • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: To fit left wrist circumference up to 15cm Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |
| ITEM-34.2 | 113426 | Wrist Brace with closure system - EXTRA SMALL, right wrist • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: To fit right wrist circumference up to 15cm Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |
| ITEM-34.3 | 156762 | Wrist Brace with closure system - SMALL, left wrist • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: To fit left wrist circumference of 15cm to 17cm Packaging requirements: Per item Estimated usage: 120 | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |
| ITEM-34.4 | 113430 | Wrist Brace with closure system - SMALL, right wrist • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: To fit right wrist circumference of 15cm to 17cm Packaging requirements: Per item Estimated usage: 120 | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 34.5 | 156761 | <p>Wrist Brace with closure system - MEDIUM, left wrist</p> <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. <p>Dimensions: To fit left wrist circumference of 17cm to 19cm</p> <p>Packaging requirements: Per item</p> <p>Estimated usage: 20</p> | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |
| ITEM-34.6 | 113431 | <p>Wrist Brace with closure system - MEDIUM, right wrist</p> <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. <p>Dimensions: To fit right wrist circumference of 17cm to 19cm</p> <p>Packaging requirements: Per item</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |
| ITEM-34.7 | 156763 | <p>Wrist Brace with closure system - LARGE, left wrist</p> <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. <p>Dimensions: To fit left wrist circumference of 19cm to 22cm</p> <p>Packaging requirements: Per item</p> <p>Estimated usage: 120</p> | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |
| ITEM-34.8 | 113432 | <p>Wrist Brace with closure system - LARGE, right wrist</p> <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. <p>Dimensions: To fit right wrist circumference of 19cm to 22cm</p> <p>Packaging requirements: Per item</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-34.9 | 156764 | <p>Wrist Brace with closure system - EXTRA LARGE, left wrist</p> <ul style="list-style-type: none"> • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. <p>Dimensions: To fit left wrist circumference of 22cm and up</p> <p>Packaging requirements: Per item</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |
| ITEM-34.10 | 113433 | <p>Wrist Brace with closure system - EXTRA LARGE, right wrist</p> <ul style="list-style-type: none"> • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. <p>Dimensions: To fit right wrist circumference of 22cm and up</p> <p>Packaging requirements: Per item</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 34.1 to 34.10 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 35 | | Paediatric Wrist Brace with closure system | | | | | | | | | | | | | | | | | |
| ITEM-35.1 | New | Paediatric Wrist Brace with closure system - Extra Small, Left • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint • Narrow strap through 1st webspace Dimensions: wrist circumference of 12-13cm, left hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |
| ITEM-35.2 | New | Paediatric Wrist Brace with closure system - Extra Small, Right • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint • Narrow strap through 1st webspace Dimensions: wrist circumference of 12-13cm, right hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |
| ITEM-35.3 | New | Paediatric Wrist Brace with closure system - Small, Right • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint • Narrow strap through 1st webspace Dimensions: wrist circumference of 13-14cm, right hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |
| ITEM-34.4 | New | Paediatric Wrist Brace with closure system - Small, Left • Breathable; durable elastic material • Latex sealed edging • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint • Narrow strap through 1st webspace Dimensions: wrist circumference of 13-14cm, left hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERIAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-35.5 | New | Paediatric Wrist Brace with closure system - Medium, Left <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: wrist circumference of 14-15cm, left hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |
| ITEM-35.6 | New | Paediatric Wrist Brace with closure system - Medium, Right <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: wrist circumference of 14-15cm, right hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |
| ITEM-35.7 | New | Paediatric Wrist Brace with closure system - Large, Right <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: wrist circumference of 15-17cm, right hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |
| ITEM-35.8 | New | Paediatric Wrist Brace with closure system - Large, Left <ul style="list-style-type: none"> Breathable; durable elastic material Latex sealed edging Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Narrow strap through 1st webspace. Dimensions: wrist circumference of 15-17cm, left hand Packaging requirements: Per item Estimated usage: 10 | Yes, representative sample for items 35.1 to 35.8 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 36 | | Wrist and Thumb Splint | | | | | | | | | | | | | | | | | |
| ITEM-36.1 | New | Wrist and Thumb Splint - Extra Small Right <ul style="list-style-type: none"> Used to prevent thumb MP joint motion, Made of breathable material with D-Ring or hook and loop straps Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit right hand circumference of 13cm to 14cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |
| ITEM-36.2 | New | Wrist and Thumb Splint - Extra Small Left <ul style="list-style-type: none"> Used to prevent thumb MP joint motion, Made of breathable material with D-Ring or hook and loop straps Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit left hand circumference of 13cm to 14cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |
| ITEM-36.3 | New | Wrist and Thumb Splint - Small Right <ul style="list-style-type: none"> Used to prevent thumb MP joint motion, Made of breathable material with D-Ring or hook and loop straps Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit right hand circumference of 14cm to 16cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-36.4 | New | Wrist and Thumb Splint - Small Left • Used to prevent thumb MP joint motion, • Made of breathable material with D-Ring or hook and loop straps • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit left hand circumference of 14cm to 16cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |
| ITEM-36.5 | New | Wrist and Thumb Splint - Medium Right • Used to prevent thumb MP joint motion, • Made of breathable material with D-Ring or hook and loop straps • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint . Must include thumb support Dimensions: To fit right hand circumference of 16cm to 19cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |
| ITEM-36.6 | New | Wrist and Thumb Splint - Medium Left • Used to prevent thumb MP joint motion, • Made of breathable material with D-Ring or hook and loop straps • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit left hand circumference of 16cm to 19cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-36.7 | New | Wrist and Thumb Splint - Large Right <ul style="list-style-type: none"> Used to prevent thumb MP joint motion, Made of breathable material with D-Ring or hook and loop straps Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit right hand circumference of 19cm to 20cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |
| ITEM-36.8 | New | Wrist and Thumb Splint - Large Left <ul style="list-style-type: none"> Used to prevent thumb MP joint motion, Made of breathable material with D-Ring or hook and loop straps Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit left hand circumference of 19cm to 20cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |
| ITEM-36.9 | New | Wrist and Thumb Splint - Extra Large Right <ul style="list-style-type: none"> Used to prevent thumb MP joint motion, Made of breathable material with D-Ring or hook and loop straps Malleable aluminium strip covered with leather or vinyl Removable volar metal stay Palmar bar should not extend past palmar crease Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support Dimensions: To fit right hand circumference of 20cm to 22cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-36.10 | New | Wrist and Thumb Splint - Extra Large Left • Used to prevent thumb MP joint motion, • Made of breathable material with D-Ring or hook and loop straps • Malleable aluminium strip covered with leather or vinyl • Removable volar metal stay • Palmar bar should not extend past palmar crease • Length of gutter should be a minimum of half of forearm length as appropriate for size of splint. Must include thumb support. Dimensions: To fit left hand circumference of 20cm to 22cm Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 36.1 to 36.10 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 37 | | Pressure garment material | | | | | | | | | | | | | | | | | |
| ITEM-37.1 | 156777 | Pressure garment material - 50m, beige • Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation • Highly durable • Is resistant to stretch - has 100% memory • Has extended shelf life • Colours: Beige Dimensions: Width of material: 180cm Packaging requirements: 50 metre roll Estimated usage: 13 | Yes, representative sample for items 37.1 to 37.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-37.2 | New | Pressure garment material - 25m, beige • Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation • Highly durable • Is resistant to stretch - has 100% memory • Has extended shelf life • Colours: Beige Dimensions: Width of material: 180cm Packaging requirements: 25 metre roll Estimated usage: 1 | Yes, representative sample for items 37.1 to 37.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-37.3 | New | Pressure garment material - 10m, beige • Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation • Highly durable • Is resistant to stretch - has 100% memory • Has extended shelf life • Colours: Beige Dimensions: Width of material: 180cm Packaging requirements: 10 metre roll Estimated usage: 1 | Yes, representative sample for items 37.1 to 37.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-37.4 | 156776 | Pressure garment material - 50m, brown • Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation • Highly durable • Is resistant to stretch - has 100% memory • Has extended shelf life • Colours: Brown Dimensions: Width of material: 180cm Packaging requirements: 50 metre roll Estimated usage: 13 | Yes, representative sample for items 37.1 to 37.6 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-37.5 | New | <p>Pressure garment material - 25m, brown</p> <ul style="list-style-type: none"> • Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation • Highly durable • Is resistant to stretch - has 100% memory • Has extended shelf life • Colours: Brown <p>Dimensions: Width of material: 180cm Packaging requirements: 25 metre roll Estimated usage: 1</p> | Yes, representative sample for items 37.1 to 37.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-37.6 | New | <p>Pressure garment material - 10m, brown</p> <ul style="list-style-type: none"> • Woven, multidirectional stretch material specifically designed for applying pressure to skin to reduce scar formation • Highly durable • Is resistant to stretch - has 100% memory • Has extended shelf life • Colours: Brown <p>Dimensions: Width of material: 180cm Packaging requirements: 10 metre roll Estimated usage: 1</p> | Yes, representative sample for items 37.1 to 37.6 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 38 | | Zipper | | | | | | | | | | | | | | | | | |
| ITEM-38.1 | New | Zipper Beige 30cm long • Can be used to fabricate pressure garments • Highly durable • Plastic zipper • Must have top and bottom stop • Colours: Beige Dimensions: • 30cm long Packaging requirements: • 10 in a box Estimated usage: 1 | Yes, representative sample for items 38.1 to 38.6 | | Each | | | | | | | | | | | | | | |
| ITEM-38.2 | New | Zipper Beige 20cm Long • Can be used to fabricate pressure garments • Highly durable • Plastic zipper • Must have top and bottom stop • Colours: Beige Dimensions: 20cm long Packaging requirements: 10 in a box Estimated usage: 1 | Yes, representative sample for items 38.1 to 38.6 | | Each | | | | | | | | | | | | | | |
| ITEM-38.3 | New | Zipper Beige 15cm Long • Can be used to fabricate pressure garments • Highly durable • Plastic zipper • Must have top and bottom stop • Colours: Beige Dimensions: 15cm long Packaging requirements: 10 in a box Estimated usage: 1 | Yes, representative sample for items 38.1 to 38.6 | | Each | | | | | | | | | | | | | | |
| ITEM-38.4 | New | Zipper Brown 30cm Long • Can be used to fabricate pressure garments • Highly durable • Plastic zipper • Must have top and bottom stop • Colours: Brown Dimensions: 30cm long Packaging requirements: 10 in a box Estimated usage: 1 | Yes, representative sample for items 38.1 to 38.6 | | Each | | | | | | | | | | | | | | |
| ITEM-38.5 | New | Zipper Brown 20cm Long • Can be used to fabricate pressure garments • Highly durable • Plastic zipper • Must have top and bottom stop • Colours: Brown Dimensions: 20cm long Packaging requirements: 10 in a box Estimated usage: 1 | Yes, representative sample for items 38.1 to 38.6 | | Each | | | | | | | | | | | | | | |
| ITEM-38.6 | New | Zipper Brown 15cm Long • Can be used to fabricate pressure garments • Highly durable • Plastic zipper • Must have top and bottom stop • Colours: Brown Dimensions: 15cm long Packaging requirements: 10 in a box Estimated usage: 1 | Yes, representative sample for items 38.1 to 38.6 | | Each | | | | | | | | | | | | | | |
| ITEM 39 | | Neoprene | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-39.1 | 156779 | Neoprene 3.2mm thick Non-Perforated • Nylon lining on one side • Pile material on the other over a rubber core, • Has a degree of elasticity resulting in comfortable and secure strapping. • Adhesive to male velcro • Latex free Dimensions: 46cm x 61cm x 3.2mm sheet Packaging requirements: Single sheet Estimated usage: 29 | Yes, representative sample for items 39.1 and 39.2 | | Per sheet | | | | | | | | | | | | | | |
| ITEM-39.2 | 156781 | Neoprene 1.5mm thick Non-Perforated • Nylon lining on one side • Pile material on the other over a rubber core, • Has a degree of elasticity resulting in comfortable and secure strapping. • Adhesive to male velcro • Latex free Dimensions: 46cm x 61 cm x 1.5mm sheet Packaging requirements: Single sheet Estimated usage: 61 | Yes, representative sample for items 39.1 and 39.2 | | Per sheet | | | | | | | | | | | | | | |
| ITEM-39.3 | New | Neoprene 3.2mm thick Perforated • Nylon lining on one side • Pile material on the other over a rubber core, • Has a degree of elasticity resulting in comfortable and secure strapping. • Adhesive to male velcro • Latex free Dimensions: 46cm x 61cm x 3.2mm sheet Packaging requirements: Single sheet Estimated usage: 29 | Yes, representative sample for items 39.3 and 39.4 | | Per sheet | | | | | | | | | | | | | | |
| ITEM-39.4 | New | Neoprene 1.5mm thick Perforated • Nylon lining on one side • Pile material on the other over a rubber core, • Has a degree of elasticity resulting in comfortable and secure strapping. • Adhesive to male velcro • Latex free Dimensions: 46cm x 61 cm x 1.5mm sheet Packaging requirements: Single sheet Estimated usage: 61 | Yes, representative sample for items 39.3 and 39.4 | | Per sheet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 40 | | Neoprene like wrap on a ROLL | | | | | | | | | | | | | | | | | |
| ITEM-40.1 | 156782 | <p>Neoprene like wrap on a ROLL; 7.6cm wide</p> <ul style="list-style-type: none"> • Piled fabric on one side • Woven blend of nylon & lycra outer layer • Elastomeric foam inner layer. • Medium conformability • Min to Moderate elasticity • Adhesive to male velcro • Mostly linear stretch • Non fraying <p>Dimensions: 7.6cm x 4.5m roll x ±2mm thick</p> <p>Packaging requirements: Per roll</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 40.1 to 40.4 | Per roll | | | | | | | | | | | | | | | |
| ITEM-40.2 | 156787 | <p>Neoprene like wrap on a ROLL; 15cm wide</p> <ul style="list-style-type: none"> • Piled fabric on one side • Woven blend of nylon & lycra outer layer • Elastomeric foam inner layer. • Medium conformability • Min to Moderate elasticity • Adhesive to male velcro • Mostly linear stretch • Non fraying <p>Dimensions: 15cm x 4.5m roll x ±2mm thick</p> <p>Packaging requirements: Per roll</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 40.1 to 40.4 | Per roll | | | | | | | | | | | | | | | |
| ITEM-40.3 | 156788 | <p>Neoprene like wrap on a ROLL; 6.4cm wide</p> <ul style="list-style-type: none"> • Piled fabric on one side • Woven blend of nylon & lycra outer layer • Elastomeric foam inner layer. • Moderate elasticity • Multi-directional stretch • Adhesive to male velcro • Non fraying • Medium conformability <p>Dimensions: 6.4cm x 4.5m roll x 1.5mm</p> <p>Packaging requirements: Per roll</p> <p>Estimated usage: 2</p> | Yes, representative sample for items 40.1 to 40.4 | Per roll | | | | | | | | | | | | | | | |
| ITEM-40.4 | 156789 | <p>Neoprene like wrap on a ROLL; 15cm wide</p> <ul style="list-style-type: none"> • Piled fabric on one side • Woven blend of nylon & lycra outer layer • Elastomeric foam inner layer. • Moderate elasticity • Multi-directional stretch • Adhesive to male velcro • Non fraying • Medium conformability <p>Dimensions: 15cm x 4.5m roll x 1.5mm</p> <p>Packaging requirements: Per roll</p> <p>Estimated usage: 37</p> | Yes, representative sample for items 40.1 to 40.4 | Per roll | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|---------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-40.5 | 156415 | Neoprene-like-wrap in a SHEET form - Black • Piled fabric on one side • Woven blend of nylon & lycra outer layer • Elastomeric foam inner layer. • High conformability • Moderate elasticity • Multi-directional stretch • Adhesive to male velcro • Non fraying • Colour: Black Dimensions: ±43.5cm x ±56.4cm x ±1.5mm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 40.5 and 40.6 | | Sheet | | | | | | | | | | | | | | |
| ITEM-40.6 | 156429 | Neoprene-like-wrap in a SHEET form - Dark blue • Piled fabric on one side • Woven blend of nylon & lycra outer layer • Elastomeric foam inner layer. • High conformability • Moderate elasticity • Multi-directional stretch • Adhesive to male velcro • Non fraying • Colour: Dark blue Dimensions: ±43.5cm x ±56.4cm x ±1.5mm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 40.5 and 40.6 | | Sheet | | | | | | | | | | | | | | |
| ITEM 41 | Padding | Padding for thermoplastic splinting - Closed cell, non perforated, self adhesive, PVC • Closed cell padding • Made of plasticized PVC with adhesive backing • Non perforated • Conformable • Cleanable Dimensions: • 100cm x 51cm x 3,2mm Packaging requirements: •4 per box Estimated usage: 1 | Yes | | Sheet | | | | | | | | | | | | | | |
| ITEM-41.1 | 156442 | Pressure relief padding • Insulates against shock and vibration • Use to cushion and protects hand splints • Conforms to all body contours • Effectively eliminating abrasion and shearing • Latex free, cleanable, 2mm thick gel pad with self adhesive backing Dimensions: 20cm x 91cm Packaging requirements: Per sheet Estimated usage: 1 | Yes | | per sheet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|---|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-41.3 | 156443 | Padding material for thermoplastic splinting – adhesive, towelling • Self adhesive • Towelling material • Non-stretch • Washable Dimensions: 50cm x 100cm x 3mm Packaging requirements: Per sheet Estimated usage: 25 | Yes | | Sheet | | | | | | | | | | | | | | |
| ITEM 42 Leather suede (for making wheelchair gloves) | | | | | | | | | | | | | | | | | | | |
| ITEM-42.1 | 174261 | Leather suede (for making wheelchair gloves) • suede leather • Colour: Black Dimensions: ± 2 mm thick Packaging requirements: Per square metre Estimated usage: 1 | Yes, representative sample for items 42.1 and 42.2 | | M ² | | | | | | | | | | | | | | |
| ITEM-42.2 | 178113 | Leather suede (for making wheelchair gloves) • suede leather • Colour: Tan Dimensions: ± 2 mm thick Packaging requirements: Per square metre Estimated usage: 1 | Yes, representative sample for items 42.1 and 42.2 | | M ² | | | | | | | | | | | | | | |
| ITEM 43 Towelling | | | | | | | | | | | | | | | | | | | |
| ITEM-43.1 | 173609 | Towelling • Piled on both sides of material • White Dimensions: Minimum 150cm wide Packaging requirements: 10 metre Estimated usage: 1 | Yes, representative sample for items 43.1 and 43.2 | | Per meter | | | | | | | | | | | | | | |
| ITEM-43.2 | 156497 | Towelling • Piled on both sides of material • Dark colours Dimensions: Minimum 150cm wide Packaging requirements: 10 metre Estimated usage: 1 | Yes, representative sample for items 43.1 and 43.2 | | Per meter | | | | | | | | | | | | | | |
| ITEM 44 Calico | | | | | | | | | | | | | | | | | | | |
| ITEM-44.1 | 156498 | Calico White • Natural cotton • Pre-shrunk • Colour: White Dimensions: Minimum 115 cm wide Packaging requirements: 10 metre Estimated usage: 1 | Yes, representative sample for items 44.1 and 44.2 | | Per meter | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|----------------------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-44.2 | 173610 | Calico Cream • Natural cotton • Pre-shrunk • Colour: Cream Dimensions: Minimum 115 cm wide Packaging requirements: 10 metre Estimated usage: 1 | Yes, representative sample for items 44.1 and 44.2 | | Per meter | | | | | | | | | | | | | | |
| ITEM 45 Foam | | | | | | | | | | | | | | | | | | | |
| ITEM-45.1 | 156499 | Foam, medium density sheet, 5cm thick • Open cell foam • Medium density Non flammable Dimensions: Minimum 180cm x 91cm x 5cm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 45.1 and 45.2 | | Sheet | | | | | | | | | | | | | | |
| ITEM-45.2 | 161592 | Foam, medium density sheet, 10cm thick • Open cell foam • Medium density Non flammable Dimensions: Minimum 180cm x 91cm x 10cm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 45.1 and 45.2 | | Sheet | | | | | | | | | | | | | | |
| ITEM-45.3 | 161593 | Foam, high density sheets, 5cm thick • Open cell foam • High density Non flammable Dimensions: Minimum 180cm x 91cm x 5cm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 45.3 and 45.4 | | Sheet | | | | | | | | | | | | | | |
| ITEM-45.4 | 156500 | Foam, high density sheets, 10cm thick • Open cell foam • High density Non flammable Dimensions: Minimum 180cm x 91cm x 10cm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 45.3 and 45.4 | | Sheet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|--------------------------------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 46 Memory Foam | | | | | | | | | | | | | | | | | | | |
| ITEM-46.1 | New | Memory Foam, medium density sheets, 2cm thick • 16-80kg/m3 , Visco elastic polymer • Medium density Non flammable Dimensions: Minimum 100cm x 180cm x 2cm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 46.1 and 46.2 | | Per sheet | | | | | | | | | | | | | | |
| ITEM-46.2 | New | Memory Foam, medium density sheets, 5cm thick • 16-80kg/m3, Visco elastic polymer • Medium density Non flammable Dimensions: Minimum 100cm x 180cm x 5cm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 46.1 and 46.2 | | Per sheet | | | | | | | | | | | | | | |
| ITEM 47 Wood, SA pine | | | | | | | | | | | | | | | | | | | |
| ITEM-47.1 | 156506 | Wood, SA Pine, 22mm thickness • South African grown pine • Solid • Planed • To be cut in meter square blocks Dimensions: •22 mm thick Packaging requirements: • Per plank; cut as specified Estimated usage: 1 | Yes | | M ² | | | | | | | | | | | | | | |
| ITEM 48 Wood, Superwood | | | | | | | | | | | | | | | | | | | |
| ITEM-48.1 | 156513 | Wood, Superwood, 3mm thick • Smooth on both sides • Durable • To be cut in meter square blocks Dimensions: .Thickness = 3mm Packaging requirements: Per plank; cut as specified Estimated usage: 1 | Yes, representative sample for items 48.1 to 48.7 | | M ² | | | | | | | | | | | | | | |
| ITEM-48.2 | 156509 | Wood, Superwood, 6mm thick • Smooth on both sides • Durable • To be cut in meter square blocks Packaging requirements: Per plank; cut as specified Estimated usage: 1 | Yes, representative sample for items 48.1 to 48.7 | | M ² | | | | | | | | | | | | | | |
| ITEM-48.3 | 156510 | Wood, Superwood, 9mm thick • Smooth on both sides • Durable • To be cut in meter square blocks Dimensions: .Thickness = 9mm Packaging requirements: Per plank; cut as specified Estimated usage: 1 | Yes, representative sample for items 48.1 to 48.7 | | M ² | | | | | | | | | | | | | | |
| ITEM-48.4 | 156511 | Wood, Superwood, 12mm thick • Smooth on both sides • Durable • To be cut in meter square blocks Dimensions: .Thickness = 12mm Packaging requirements: Per plank; cut as specified Estimated usage: 1 | Yes, representative sample for items 48.1 to 48.7 | | M ² | | | | | | | | | | | | | | |
| ITEM-48.5 | 156515 | Wood, Superwood, 18mm thick • Smooth on both sides • Durable • To be cut in meter square blocks Dimensions: .Thickness =18mm Packaging requirements: Per plank; cut as specified Estimated usage: 1 | Yes, representative sample for items 48.1 to 48.7 | | M ² | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-48.6 | 156517 | Wood, Superwood, 22mm thick • Smooth on both sides • Durable • To be cut in meter square blocks Dimensions: .Thickness = 22mm Packaging requirements: Per plank; cut as specified Estimated usage: 1 | Yes, representative sample for items 48.1 to 48.7 | | M ² | | | | | | | | | | | | | | |
| ITEM-48.7 | 156522 | Wood, Superwood, 32mm thick • Smooth on both sides • Durable • To be cut in meter square blocks Dimensions: .Thickness = 32mm Packaging requirements: Per plank; cut as specified Estimated usage: 1 | Yes, representative sample for items 48.1 to 48.7 | | M ² | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 49 | | Stretch loop, non-adhesive (equivalent to velcro) | | | | | | | | | | | | | | | | | |
| ITEM-49.1 | 156601 | Stretch loop, 2.5cm, non-adhesive (equivalent to velcro) • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: white Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 2 | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-49.2 | 156602 | Stretch loop, 2.5cm, non-adhesive (equivalent to velcro) • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: beige Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 2 | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-49.3 | 173617 | Stretch loop, 2.5cm, non-adhesive (equivalent to velcro) • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: Dark Brown Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-49.4 | New | Stretch loop, 2.5cm, non-adhesive (equivalent to velcro) • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: Black Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-49.5 | 156603 | Stretch loop, 5cm, non-adhesive (equivalent to velcro) • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: white Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 15 | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|---|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-49.6 | 156604 | <p>Stretch loop, 5cm, non-adhesive (equivalent to velcro)</p> <ul style="list-style-type: none"> • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: Beige <p>Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 1</p> | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-49.7 | 173618 | <p>Stretch loop, 5cm, non-adhesive (equivalent to velcro)</p> <ul style="list-style-type: none"> • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: Dark Brown <p>Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 1</p> | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-49.8 | New | <p>Stretch loop, 5cm, non-adhesive (equivalent to velcro)</p> <ul style="list-style-type: none"> • Stretch loop • Elastic loop provides even tension • Made of rubber with soft edges • Moderate resistance to stretch for securing splints • Snug fit and secure closure • Mates with all Velcro type hook strapping • Colours: Black <p>Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 1</p> | Yes, representative sample for items 49.1 to 49.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM 50 | Soft loop, non-adhesive (equivalent to velcro) | | | | | | | | | | | | | | | | | | |
| ITEM-50.1 | 156523 | <p>Soft loop, 2.5cm, non-adhesive (equivalent to velcro)</p> <p>Loop – 2.5cm wide</p> <ul style="list-style-type: none"> • Soft • Non-abrasive loop Velcro • Highly durable • Maximum bonding capability • Colours: white <p>Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 39</p> | Yes, representative sample for items 50.1 to 50.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-50.2 | 156525 | <p>Soft loop, 5cm, non-adhesive (equivalent to velcro)</p> <p>Loop – 5cm wide</p> <ul style="list-style-type: none"> • Soft • Non-abrasive loop Velcro • Highly durable • Maximum bonding capability • Colours: white <p>Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 50</p> | Yes, representative sample for items 50.1 to 50.6 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-50.3 | 156524 | Soft loop, 2.5cm, non-adhesive (equivalent to velcro) Loop – 2.5cm wide • Soft • Non-abrasive loop Velcro • Highly durable • Maximum bonding capability • Colours: beige Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 9 | Yes, representative sample for items 50.1 to 50.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-50.4 | 156526 | Soft loop 5cm non-adhesive (equivalent to velcro) Loop – 5cm wide • Soft • Non-abrasive loop Velcro • Highly durable • Maximum bonding capability • Colours: beige Dimensions: . 25m x 5cm Packaging requirements: . 25m reel Estimated usage: 9 | Yes, representative sample for items 50.1 to 50.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-50.5 | 174266 | Soft loop, 2.5cm. non-adhesive (equivalent to velcro) Loop – 2.5cm wide • Soft • Non-abrasive loop Velcro • Highly durable • Maximum bonding capability • Colours: Dark Brown Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 50.1 to 50.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-50.6 | 174269 | Soft loop, 5cm, non-adhesive (equivalent to velcro) Loop – 5cm wide • Soft • Non-abrasive loop Velcro • Highly durable • Maximum bonding capability • Colours: Dark Brown Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 50.1 to 50.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM 51 | Loop, non-adhesive (equivalent to velcro) | | | | | | | | | | | | | | | | | | |
| ITEM-51.1 | 164587 | Loop, non-adhesive, 2.5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: White Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.2 | 164588 | Loop, non-adhesive, 5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: White Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-51.3 | 118877 | Loop, non-adhesive, 2cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: White Dimensions: 50m x 2cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.4 | 118265 | Loop, non-adhesive, 2.5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: Beige Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.5 | 118264 | Loop, non-adhesive, 5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: Beige Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.6 | 173614 | Loop, non-adhesive, 2cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: Brown Dimensions: 50m x 2cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.7 | 171336 | Loop, non-adhesive, 2.5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: Brown Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.8 | 171326 | Loop, non-adhesive, 5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: Brown Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.9 | 164567 | Loop, non-adhesive, 2.5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: Black Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM-51.10 | 164569 | Loop, non-adhesive, 5cm (equivalent to velcro) Loop, non-adhesive, nylon Colours: Black Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 90 rolls over range of sizes | Yes, representative sample for items 51.1 to 51.10 | | Per roll | | | | | | | | | | | | | | |
| ITEM 52 | | Hook, non-adhesive (equivalent to velcro) | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-52.1 | 164561 | Hook, non-adhesive, 2.5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: White Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.2 | 156527 | Hook, non-adhesive, 2.5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: White Dimensions: 25m x 2.5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.3 | 164563 | Hook, non-adhesive, 5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: White Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.4 | 156529 | Hook, non-adhesive, 5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: White Dimensions: 25m x 5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.5 | 156528 | Hook, non-adhesive, 2.5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: Beige Dimensions: 25m x 2.5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.6 | 156530 | Hook, non-adhesive, 5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: Beige Dimensions: 25m x 5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.7 | 171299 | Hook, non-adhesive, 2.5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: Brown Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.8 | 171311 | Hook, non-adhesive, 5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: Brown Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-52.9 | 173612 | Hook, non-adhesive, 2.5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: Dark Brown Dimensions: 25m x 2.5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.10 | 173613 | Hook, non-adhesive, 5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: Dark Brown Dimensions: 25m x 5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.11 | 164557 | Hook, non-adhesive, 2.5cm (equivalent to velcro) Hook, non-adhesive, nylon Colours: Black Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |
| ITEM-52.12 | 164560 | Hook, non-adhesive, 5cm (equivalent to velcro) Hook non adhesive nylon Colours: Black Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 52.1 to 52.12 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 53 | | Self adhesive loop (equivalent to velcro) | | | | | | | | | | | | | | | | | |
| ITEM-53.1 | 173615 | Self adhesive loop 2.5cm (equivalent to velcro) Loop – 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: white • Highly durable Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 8 | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-53.2 | 171273 | Self adhesive loop 5cm (equivalent to velcro) Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: white • Highly durable Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 45 rolls over range of sizes | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-53.3 | 118878 | Self adhesive loop 2.5cm (equivalent to velcro) Loop - 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: beige • Highly durable Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 45 rolls over range of sizes | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-53.4 | 118879 | Self adhesive loop 5cm (equivalent to velcro) Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: beige • Highly durable Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 45 rolls over range of sizes | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-53.5 | 171336 | Self adhesive loop 2.5cm (equivalent to velcro) Loop - 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: brown • Highly durable Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 45 rolls over range of sizes | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-53.6 | 171279 | Self adhesive loop 5 cm (equivalent to velcro) Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: brown • Highly durable Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 45 rolls over range of sizes | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-53.7 | 171284 | Self adhesive loop 2.5 cm (equivalent to velcro) Loop - 2.5cm self adhesive, nylon • Maximum bonding capability • Colours: black • Highly durable Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated requirements: 45 rolls over range of sizes | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-53.8 | 171289 | <p>Self adhesive loop 5 cm (equivalent to velcro) Loop - 5cm self adhesive, nylon • Maximum bonding capability • Colours: black • Highly durable Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 45 rolls over range of sizes</p> | Yes, representative sample for items 53.1 to 53.8 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 54 | | Hook self adhesive (equivalent to velcro) | | | | | | | | | | | | | | | | | |
| ITEM-54.1 | 156531 | Hook self adhesive 2cm (equivalent to velcro) Hook – 2cm self adhesive • Maximum bonding capability • Colours: white • Highly durable Dimensions: 25m x 2cm Packaging requirements: • 25m reel Estimated usage: 0 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.2 | 156533 | Hook self adhesive 2,5cm (equivalent to velcro) Hook – 2,5cm self adhesive • Maximum bonding capability • Colours: white • Highly durable Dimensions: 25m x 2,5cm Packaging requirements: • 25m reel Estimated usage: 42 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.3 | 171270 | Hook self adhesive 2,5cm (equivalent to velcro) Hook - 2,5cm self adhesive Maximum bonding capability • Colours: white • Highly durable Dimensions: 50m x 2,5cm Packaging requirements: Per roll Estimated quantity: 45 rolls over range of sizes | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.4 | 156535 | Hook self adhesive 5cm (equivalent to velcro) Hook – 5cm self adhesive • Maximum bonding capability • Colours: white • Highly durable Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 11 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.5 | 164583 | Hook self adhesive 5cm (equivalent to velcro) Hook - 5cm self adhesive. • Maximum bonding capability • Colours: white • Highly durable Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 10 rolls over range of sizes | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.6 | 156532 | Hook self adhesive 2cm (equivalent to velcro) Hook – 2cm self adhesive • Maximum bonding capability • Colours: beige • Highly durable Dimensions: 25m x 2cm Packaging requirements: • 25m reel Estimated usage: 1 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.7 | 156534 | Hook self adhesive 2,5cm (equivalent to velcro) Hook – 2,5cm self adhesive • Maximum bonding capability • Colours: beige • Highly durable Dimensions: 25m x 2,5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-54.8 | 156536 | Hook self adhesive 5cm equivalent to velcro Hook – 5cm self adhesive • Maximum bonding capability • Colours: beige • Highly durable Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 9 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.9 | 173614 | Hook self adhesive 2cm (equivalent to velcro) Hook – 2cm self adhesive • Maximum bonding capability • Colours: Dark Brown • Highly durable Dimensions: 25m x 2cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.10 | 174272 | Hook self-adhesive 2,5cm (equivalent to velcro) Hook – 2,5cm self adhesive • Maximum bonding capability • Colours: Dark Brown • Highly durable Dimensions: 25m x 2,5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.11 | 173616 | Hook self adhesive 5cm (equivalent to velcro) Hook – 5cm self adhesive • Maximum bonding capability • Colours: Dark Brown • Highly durable Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.12 | 164584 | Hook self adhesive 2,5cm (equivalent to velcro) Hook - 2,5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable Dimensions: 50m x 2,5cm Packaging requirements: Per roll Estimated quantity: 10 rolls over range of sizes | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.13 | 164586 | Hook self adhesive 5cm (equivalent to velcro) Hook - 5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated requirements: 10 rolls over range of sizes | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM-54.14 | 164557 | Hook self adhesive 2,5cm (equivalent to velcro) Hook - 2,5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable Dimensions: 50m x 2,5cm Packaging requirements: Per roll Estimated requirements: 10 rolls over range of sizes | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-54.15 | 171268 | Hook self adhesive 5cm (equivalent to velcro) Hook - 5cm self adhesive. • Maximum bonding capability • Colours: Brown • Highly durable Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 10 rolls over range of sizes | Yes, representative sample for items 54.1 to 54.15 | | Per roll | | | | | | | | | | | | | | |
| ITEM 55 | | Hook maximum bonding (equivalent to velcro) | | | | | | | | | | | | | | | | | |
| ITEM-55.1 | 156527 | Hook, maximum bonding, 2.5cm (equivalent to velcro) Hook – 2.5cm wide • Highly durable • Maximum bonding capability • Colours: white Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 2 | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-55.2 | 164563 | Hook, maximum bonding, 5cm (equivalent to velcro) Hook - 5cm wide • Highly durable • Maximum bonding capability • Colours: white Dimensions: 50m x 5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-55.3 | 156529 | Hook, maximum bonding, 5cm (equivalent to velcro) Hook – 5cm wide • Highly durable • Maximum bonding capability • Colours: white Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 3 | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-55.4 | 156528 | Hook, maximum bonding, 2.5cm (equivalent to velcro) Hook – 2.5cm wide • Highly durable • Maximum bonding capability • Colours: beige Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-55.5 | 156530 | Hook, maximum bonding, 5cm (equivalent to velcro) Hook – 5cm wide • Highly durable • Maximum bonding capability • Colours: beige Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-55.6 | 171299 | Hook, maximum bonding, 2.5cm (equivalent to velcro) Hook - 2.5cm wide • Highly durable • Maximum bonding capability • Colours: brown Dimensions: 50m x 2.5cm Packaging requirements: Per roll Estimated quantity: 25 rolls over range of sizes | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|---|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-55.7 | 173612 | Hook, maximum bonding, 2.5cm (equivalent to velcro) Hook – 2.5cm wide • Highly durable • Maximum bonding capability • Colours: Dark Brown Dimensions: 25m x 2.5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM-55.8 | 173613 | Hook, maximum bonding, 5cm (equivalent to velcro) Hook – 5cm wide • Highly durable • Maximum bonding capability • Colours: Dark Brown Dimensions: 25m x 5cm Packaging requirements: 25m roll Estimated usage: 1 | Yes, representative sample for items 55.1 to 55.8 | | Per roll | | | | | | | | | | | | | | |
| ITEM 56 | Self gripping (equivalent to velcro) | | | | | | | | | | | | | | | | | | |
| ITEM-56.1 | New | Self gripping 25mm (equivalent to velcro) Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. Colour: Tan/beige Dimensions: 25m x 25mm Packaging requirements: 25m roll Estimated usage: New | Yes, representative sample for items 56.1 to 56.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-56.2 | New | Self gripping 50mm (equivalent to velcro) Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. Colour: Tan/beige Dimensions: 25m x 50mm Packaging requirements: 25m roll Estimated usage: New | Yes, representative sample for items 56.1 to 56.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-56.3 | New | Self gripping 25mm (equivalent to velcro) Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. Colour: Black Dimensions: 25m x 25mm Packaging requirements: 25m roll Estimated usage: New | Yes, representative sample for items 56.1 to 56.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-56.4 | New | Self gripping 50mm (equivalent to velcro) Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. Colour: Black Dimensions: 25m x 50mm Packaging requirements: 25m roll Estimated usage: New | Yes, representative sample for items 56.1 to 56.6 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERIAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-56.5 | New | Self gripping 25mm (equivalent to velcro) Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. Colour: White Dimensions: 25m x 25mm Packaging requirements: 25m roll Estimated usage: New | Yes, representative sample for items 56.1 to 56.6 | | Per roll | | | | | | | | | | | | | | |
| ITEM-56.6 | New | Self gripping 50mm (equivalent to velcro) Self gripping, with integrated hook and loops that adheres to itself. Highly durable. Maximum bonding capability. Colour: White Dimensions: 25m x 50mm Packaging requirements: 25m roll Estimated usage: New | Yes, representative sample for items 56.1 to 56.6 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 57 | | Silicon Gel Sheets | | | | | | | | | | | | | | | | | |
| ITEM-57.1 | 118267 | Silicon Gel Sheets Reduce the appearance of scars; Re-usable; Self- adhesive; Dimensions: 12 x 15cm Packaging Requirements: Pack of 10 sheets Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 58 | | Silicon | | | | | | | | | | | | | | | | | |
| ITEM-58.1 | New | Silicon Gel 15ml Reduce the appearance of scars; improves pliability and mobility Self drying silicon topical gel that is absorbable by the skin and can be applied to microporous tape Ingredients: Silicon gel, centella, asiatica Dimensions: 15ml tube Packaging Requirements: Pack of 5 Estimated usage: 1 | Yes, representative sample for items 58.1 and 58.2 | | Per pack | | | | | | | | | | | | | | |
| ITEM-58.2 | New | Silicon Gel 50ml Reduce the appearance of scars; improves pliability and mobility Self drying silicon topical gel that is absorbable by the skin and can be applied to microporous tape Ingredients: Silicon gel, centella, asiatica Dimensions: 50ml tube Packaging Requirements: Each Estimated usage: 1 | Yes, representative sample for items 58.1 and 58.2 | | Each | | | | | | | | | | | | | | |
| ITEM-58.3 | New | Silicon Spray Treatment of scars to improve pliability, motility and appearance. Self drying silicon spray Ingredients: Silicon dioxide, polysiloxanes. Dimensions: 100ml Aerosol Can Packaging Requirements: Box of 6 Estimated usage: 1 | Yes | | Per box | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|----------------------------|--|--|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 59 | Squeeze Ball | | | | | | | | | | | | | | | | | | |
| ITEM-59.1 | 118263 | Squeeze Ball Hand Exerciser 6 cm Palm size ball; Returns to original diameter; Dimensions: 6cm diameter Packaging Requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 60 | Pen and Pencil holder | | | | | | | | | | | | | | | | | | |
| ITEM-60.1 | 118194 | Pen and Pencil holder Soft PVC moulding that fits onto a normal sized pen/pencil to make them easier to grip and control writing. Dimensions: Internal Diameter 8mm. Packaging Requirements: Pack of 3 Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 61 | Long-Handled sponge | | | | | | | | | | | | | | | | | | |
| ITEM-61.1 | 118193 | Long-Handled sponge Mounted on plastic mouldable handle; Dimensions: Length 57.2cm. Packaging Requirements: Each Estimated usage: 1 | Yes, representative sample for items 61.1 and 61.2 | | Each | | | | | | | | | | | | | | |
| ITEM-61.2 | New | Long-Handled sponge Mounted on plastic mouldable handle; Dimensions: Length of more than 60cm Packaging Requirements: Each Estimated usage: 1 | Yes, representative sample for items 61.1 and 61.2 | | Each | | | | | | | | | | | | | | |
| ITEM 62 | Sock Aid | | | | | | | | | | | | | | | | | | |
| ITEM-62.1 | 118192 | Sock Aid Cord attached assists in pulling up sock over foot; Cord style; Dimensions: Width 10cm; Cord style 25.5cm. Packaging Requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 63 | Button hook and Zip puller | | | | | | | | | | | | | | | | | | |
| ITEM-63.1 | 118191 | Button hook and Zip puller Assist in tying buttons and pulling up zips with one hand. Packaging Requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|-----------------|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 64 | | Plate guard | | | | | | | | | | | | | | | | | |
| ITEM-64.1 | 118741 | Plate guard Fits around ordinary plate to assist with one handed eating; Microwave and Dishwasher safe. Dimensions: 19 to 25cm Packaging Requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 65 | | Finger coil spring | | | | | | | | | | | | | | | | | |
| ITEM-65.1 | 118190 | Finger coil spring stainless steel spring; Gauge 0,7mm; Dimensions: Coil Diameter of 5mm. Packaging Requirements: Pack of 10 Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 66 | | Knuckle bender coil spring | | | | | | | | | | | | | | | | | |
| ITEM-66.1 | 118188 | Knuckle bender coil spring Stainless steel spring; Designed to assist intrinsic muscle weakness; Gauge 1mm. Torque force at full finger extension: ca.1000g.cm Dimensions: Coil Diameter 1cm. Packaging Requirements: Pack of 10 Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 67 | | Finger extension assist | | | | | | | | | | | | | | | | | |
| ITEM-67.1 | 118187 | Finger extension assist Stainless steel spring; Assist in extending fingers. Gauge 1mm. At total finger flexion the spring is tightened up to a torque force of : ca.2000g.cm Dimensions: Coil Diameter 4mm. Packaging Requirements: Pack of 10 Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 68 | | Wrist extension coil spring | | | | | | | | | | | | | | | | | |
| ITEM-68.1 | 118180 | Wrist extension coil spring Stainless steel spring; Two toes - left and right. Gauge 1,5 mm. Torque force when tightened over an angle of 130 degrees: 5000g.cm Diameter: 15mm. Packaging Requirements: Pack of 10 (5 left and 5 right) Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 69 | | Plastic shoehorn | | | | | | | | | | | | | | | | | |
| ITEM 69.1 | 124872 | Plastic shoehorn To assist in putting on shoes. Dimensions: Length 43cm; Weight 97g Packaging Requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 70 | | Webbing | | | | | | | | | | | | | | | | | |
| ITEM-70.1 | 156599 | Webbing 2.5cm • Cotton • Non-elastic Dimensions: 2.5cm wide on a roll Packaging requirements: Single roll Estimated usage: 6 | Yes, representative sample for items 70.1 and 70.2 | | Per roll | | | | | | | | | | | | | | |
| ITEM-70.2 | 156600 | Webbing 5cm • Cotton • Non-elastic Dimensions: 5cm wide on a roll Packaging requirements: Single roll Estimated usage: 6 | Yes, representative sample for items 70.1 and 70.2 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------------------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| item 71 | Soft Cast | | | | | | | | | | | | | | | | | | |
| ITEM-71.1 | 178631 | Soft Cast • Semi-rigid casting tape Dimensions: 5cm width Packaging requirements: Box of 10 Estimated usage: 1 | Yes, representative sample for items 71.1 to 71.3 | | Per box | | | | | | | | | | | | | | |
| ITEM-71.2 | 113434 | Soft Cast • Semi-rigid casting tape Dimensions: 7,5cm width Packaging requirements: Box of 10 Estimated usage: 1 | Yes, representative sample for items 71.1 to 71.3 | | Per box | | | | | | | | | | | | | | |
| ITEM-71.3 | 117664 | Soft Cast • Semi-rigid casting tape Dimensions: 10cm width Packaging requirements: Box of 10 Estimated usage: 1 | Yes, representative sample for items 71.1 to 71.3 | | Per box | | | | | | | | | | | | | | |
| item 72 | Elastic thread | | | | | | | | | | | | | | | | | | |
| ITEM-72.1 | 178638 | Elastic thread Stretch ability at least 200%; +/- 200g. Dimensions: Spool of 10m Packaging Requirements: 10 meter Estimated usage: 1 | Yes | | Per spool | | | | | | | | | | | | | | |
| ITEM-72.2 | 178639 | Elastic thread Stretch ability at least 200%; +/- 400g. Dimensions: Spool of 10m Packaging Requirements: 10 meter Estimated usage: 1 | Yes | | Per spool | | | | | | | | | | | | | | |
| ITEM-72.3 | 117665 | Elastic thread Stretch ability at least 200%; +/- 600g. Dimensions: Spool of 10m Packaging Requirements: 10 meter Estimated usage: 1 | Yes | | Per spool | | | | | | | | | | | | | | |
| item 73 | Thermoplastic tube | | | | | | | | | | | | | | | | | | |
| ITEM-73.1 | 178644 | Thermoplastic tube bender • plastic bender without stress on the thumbs • Intended for bending thermoplastic orifices Packaging requirements: Single item Estimated usage: 1 | Yes | | Per item | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|---|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-73.2 | 178659 | Thermoplastic tube <ul style="list-style-type: none"> Thermoplastic tube intended for making dynamic hand orthoses Dimensions: Length 33cm, external diameter 6.5mm, internal diameter 3.3mm Packaging requirements: Packet of 10 tubes with two wires Estimated usage: 1 | Yes | | Per packet | | | | | | | | | | | | | | |
| ITEM 74 | Outrigger Wire Copper Coated steel | | | | | | | | | | | | | | | | | | |
| ITEM-74.1 | 113437 | Outrigger Wire Copper Coated steel, 1.6mm Welding Rod 1.6mm Bends easily with pliers or with bender Length: 12 x 91cm Packaging requirements: Pack of 10 Estimated usage: 1 | Yes, representative sample for items 74.1 to 74.3 | | Per pack | | | | | | | | | | | | | | |
| ITEM-74.2 | 113723 | Outrigger Wire Copper Coated steel, 2.4mm Welding Rod 2.4mm Bends easily with pliers or with bender Length: 12 x 91cm Packaging requirements: Pack of 10 Estimated usage: 1 | Yes, representative sample for items 74.1 to 74.3 | | Per pack | | | | | | | | | | | | | | |
| ITEM-74.3 | 113724 | Outrigger Wire Copper Coated steel, 3.2mm Welding Rod 3.2mm Bends easily with pliers or with bender Length: 12 x 91cm Packaging requirements: Pack of 10 Estimated usage: 1 | Yes, representative sample for items 74.1 to 74.3 | | Per pack | | | | | | | | | | | | | | |
| ITEM 75 | Tap turners | | | | | | | | | | | | | | | | | | |
| ITEM-75.1 | 156618 | Tap turners for cross type taps <ul style="list-style-type: none"> Strong moulded tap turner used to open a cross shaped tap head with leverage Hooks over the top of a '+' type tap, to make it easy to turn. Supplied as a pair for hot and cold taps. Dimensions: Length ±140mm, Weight = 65g Packaging requirements: Per pair Estimated usage: 54 | Yes | | Per pair | | | | | | | | | | | | | | |
| ITEM-75.2 | 156619 | Tap turners for crystal type taps <ul style="list-style-type: none"> Highly durable tap turner which clamps around the top of cylindrical tap tops and provides the leverage to turn them easily Barrel of the handle is spring loaded to allow the jaws to clamp down on the tap, which is gripped by soft non-slip rubber linings in the jaws. Jaws can also be locked in position. Comfortable and easy to grip handle. For use with taps up to 60 mm in diameters. Dimensions: Closed length = 190mm, Weight = 112g Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 76 | Bath bench | | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|-----------------|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-76.1 | 156628 | <p>Bath bench</p> <ul style="list-style-type: none"> • To fit over top of bath to facilitate transfer into bath / enable client to sit above bath to wash with hand shower. • Must be adjustable to fit any bath and fit / secure safely to bath. <p>Dimensions: To fit any bath Packaging requirements: Single unit Estimated usage: 1</p> | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 77 | | Bath seat | | | | | | | | | | | | | | | | | |
| ITEM-77.1 | 156629 | <p>Bath seat</p> <ul style="list-style-type: none"> • Seat / bench which fits firmly into bath (seat $\pm\frac{1}{2}$ way between bath rim and bottom of bath). • Adjustable to fit firmly into contour of various standard baths to enable client to sit safely above the bottom of the bath or transfer safely into the bottom of the bath. • Non-slip handles. <p>Dimensions: To fit standard bath Packaging requirements: Single unit Estimated usage: 4</p> | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 78 | | Rubber band post | | | | | | | | | | | | | | | | | |
| ITEM-78.1 | 156605 | Rubber band post • Used to anchor rubber bands Dimensions: 6.4mm x 13mm Packaging requirements: Pack of 10 Estimated usage: 10 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 79 | | Closed cell foam tubing for assistive devices | | | | | | | | | | | | | | | | | |
| ITEM-79.1 | 156606 | Closed cell foam tubing for assistive devices, Size 1 • Closed cell foam tubing to build up tool & utensil handles to give greater control to people who lack coordination, strength or gripping ability • Slip resistant outer layer • Dishwasher safe • Can be cut to the required length • Available in bright colours Dimensions: • Length = 300mm, • Inner diameter = 6mm, • Outer diameter 22mm Packaging requirements: Pack of 6 Estimated usage: 159 | Yes, representative sample for items 79.1 to 79.3 | | Pack of 6 | | | | | | | | | | | | | | |
| ITEM-79.2 | 156607 | Closed cell foam tubing for assistive devices, Size 2 • Closed cell foam tubing to build up tool & utensil handles to give greater control to people who lack coordination, strength or gripping ability • Slip resistant outer layer • Dishwasher safe • Can be cut to the required length • Available in bright colours Dimensions: • Length = 300mm, • Inner diameter = 10mm, • Outer diameter 28mm Packaging requirements: Pack of 6 Estimated usage: 150 | Yes, representative sample for items 79.1 to 79.3 | | Pack of 6 | | | | | | | | | | | | | | |
| ITEM-79.3 | 156608 | Closed cell foam tubing for assistive devices, Size 3 • Closed cell foam tubing to build up tool & utensil handles to give greater control to people who lack coordination, strength or gripping ability • Slip resistant outer layer • Dishwasher safe • Can be cut to the required length • Available in bright colours Dimensions: • Length = 300mm, • Inner diameter = 16mm, • Outer diameter 28mm Packaging requirements: Pack of 6 Estimated usage: 100 | Yes, representative sample for items 79.1 to 79.3 | | Pack of 6 | | | | | | | | | | | | | | |
| ITEM 80 | | Non slip rubber sheeting | | | | | | | | | | | | | | | | | |
| ITEM-80.1 | 156610 | Non slip rubber sheeting • Non woven fabric, coated with non-slip properties. • Suitable for moulding or wrapping around an object. Dimensions: 40cm x 9m x 1mm Packaging requirements: 9 meter roll Estimated usage: 3 | Yes | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|----------------------------------|---|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 81 | Adapted knife with dagger handle | | | | | | | | | | | | | | | | | | |
| ITEM-81.1 | 156615 | Adapted knife with dagger handle with serrated blade • High quality stainless steel blade at ± 90° angle to handle to keep wrist in stress-free position. • Serrated edge- ideal for slicing bread. • Dishwasher safe. Dimensions: • Blade length = 254mm, • Weight 125g Packaging requirements: Each <i>Estimated usage: 1</i> | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-81.2 | 156616 | Adapted knife with dagger handle with scalloped blade • Stainless steel blade at ± 90° angle to handle to keep wrist in stress-free position. • Permanently sharp blade • Blade is ideal for cutting meat and vegetables • Dishwasher safe. Dimensions: • Blade length = 203mm, • Weight 155g Packaging requirements: Each <i>Estimated usage: 1</i> | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 82 | Long handled reacher | | | | | | | | | | | | | | | | | | |
| ITEM-82.1 | 156617 | Long handled reacher • Ultra light weight aluminium shaft • Assist those with restricted reach or mobility problems • Single moving jaw operated with a trigger action from the handle Dimensions: • Weight ±180 grams, • Length ± 600 mm Packaging requirements: Each <i>Estimated usage: 0</i> | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 83 | | Raised toilet seat | | | | | | | | | | | | | | | | | |
| ITEM-83.1 | 156625 | <p>Raised toilet seat</p> <ul style="list-style-type: none"> Light weight, plastic raised toilet seat with non-slip clips / brackets to attach firmly to / detach from most toilets easily. Has a front cut-out to facilitate personal hygiene and a slight slope from back to front making sitting & rising easier. Height-adjustable: can be raised 50mm, 100mm and 150mm above toilet rim Moulded white polyethylene/ or similar compound that is chemically resilient, and easy to clean. <p>Dimensions: • To fit standard toilet, • Height from toilet rim 50mm, 100mm and 150mm.</p> <p>Packaging requirements: Each</p> <p>Estimated usage: 71</p> | Yes | Each | | | | | | | | | | | | | | | |
| ITEM 84 | | Nylon D rings | | | | | | | | | | | | | | | | | |
| ITEM-84.1 | 162856 | <p>Nylon D rings, 25mm, White</p> <ul style="list-style-type: none"> Nylon D-rings- rectangular White <p>Dimensions: 25 mm</p> <p>Packaging requirements: 50 per pack</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 84.1 and 84.2 | Per packet | | | | | | | | | | | | | | | |
| ITEM-84.2 | 162857 | <p>Nylon D rings, 50mm, White</p> <ul style="list-style-type: none"> Nylon D-rings- rectangular White <p>Dimensions: 50mm</p> <p>Packaging requirements: 25 per pack</p> <p>Estimated usage: 1</p> | Yes, representative sample for items 84.1 and 84.2 | Per packet | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERIAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 85 | | PATIENT TRANSFER ROLLBOARD | | | | | | | | | | | | | | | | | |
| ITEM-85.1 | 161149 | PATIENT TRANSFER ROLLBOARD, 170cm A patient transfer roller-slide to provide a painless and comfortable transfer for patients with special needs. <input type="checkbox"/> Shape retaining sliding board covered with polyvinylchloride and PU-covered nylon <input type="checkbox"/> Able to bridge a gap of approx. 20cm, without bending, bearing a patient weighing up to 100kg <input type="checkbox"/> Allows for a patient to be transferred between to points even if there is a height difference <input type="checkbox"/> X-ray translucent <input type="checkbox"/> Foldable and lightweight <input type="checkbox"/> Removable machine washable cover, made of Ultra Low Friction material <input type="checkbox"/> Able to withstand surface disinfectants with 70% Ethanol, 1% Virkon or 45% Isopropyl alcohol Dimensions: Approx. 170cm x 50cm Packaging requirements: Single unit Estimated usage: 10 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-85.2 | 161150 | PATIENT TRANSFER ROLLBOARD, 85cm A patient transfer roller-slide to provide a painless and comfortable transfer for patients with special needs. <input type="checkbox"/> Shape retaining sliding board covered with polyvinylchloride and PU-covered nylon <input type="checkbox"/> Able to bridge a gap of approx. 20cm, without bending, bearing a patient weighing up to 100kg <input type="checkbox"/> Allows for a patient to be transferred between to points even if there is a height difference <input type="checkbox"/> X-ray translucent <input type="checkbox"/> Foldable and lightweight <input type="checkbox"/> Removable machine washable cover, made of Ultra Low Friction material <input type="checkbox"/> Able to withstand surface disinfectants with 70% Ethanol, 1% Virkon or 45% Isopropyl alcohol Dimensions: Approx. 85cm x 50cm Packaging requirements: Single unit Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 86 | | PATIENT TRANSFER SHEET / ROLL-OVER SLIDE TYPE | | | | | | | | | | | | | | | | | |
| ITEM-86.1 | 161174 | PATIENT TRANSFER SHEET / ROLL-OVER SLIDE TYPE, 50cm x 60cm A patient transfer sheet to provide a painless and comfortable lateral transfer for patients, also to turn the patient. <input type="checkbox"/> A soft one piece padded roll bag which has flame retardant properties <input type="checkbox"/> Cotton outer material and inner Ultra low friction material <input type="checkbox"/> Has a polyurethane coated Nylon protective cover <input type="checkbox"/> Machine washable and or autoclavable <input type="checkbox"/> Able to withstand surface disinfectants with 45% Isopropyl alcohol Dimensions: Approx. 50cm x 60cm Packaging requirements: Single unit Estimated usage: 0 | Yes, representative sample for items 86.1 and 86.2 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|--|--------|--|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-86.2 | 161173 | PATIENT TRANSFER SHEET / ROLL-OVER SLIDE TYPE, 190cm x 60cm A patient transfer sheet to provide a painless and comfortable lateral transfer for patients, also to turn the patient. <input type="checkbox"/> A soft one piece padded roll bag which has flame retardant properties <input type="checkbox"/> Cotton outer material and inner Ultra low friction material <input type="checkbox"/> Has a polyurethane coated Nylon protective cover <input type="checkbox"/> Machine washable and or autoclavable <input type="checkbox"/> Able to withstand surface disinfectants with 45% Isopropyl alcohol Dimensions: Approx. 190 x 60cm Packaging requirements: Single unit Estimated usage: 1 | Yes, representative sample for items 86.1 and 86.2 | | Each | | | | | | | | | | | | | | |
| ITEM 87 DISPOSABLE COVER FOR TRANSFER SHEET/ ROLL-OVER SLIDE TYPE | | | | | | | | | | | | | | | | | | | |
| ITEM-87.1 | 161175 | DISPOSABLE COVER FOR TRANSFER SHEET/ ROLL-OVER SLIDE TYPE, 190cm x 60cm <input type="checkbox"/> Disposable covers shall be available for use with the transfer sheet Dimensions: Approx. 190 x 60cm Packaging requirements: Roll of 100 Estimated usage: 1 | Yes, representative sample for items 87.1 and 87.2 | | Per roll | | | | | | | | | | | | | | |
| ITEM-87.2 | 161176 | DISPOSABLE COVER FOR TRANSFER SHEET/ ROLL-OVER SLIDE TYPE, 50cm x 60cm <input type="checkbox"/> Disposable covers shall be available for use with transfer sheet Dimension: 50cm X 60cm Packaging requirements: Roll of 100 Estimated usage: 1 | Yes, representative sample for items 87.1 and 87.2 | | Per roll | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 88 | | PATIENT TRANSFER SUPPORT AND MOVING AID | | | | | | | | | | | | | | | | | |
| ITEM-88.1 | 161159 | PATIENT TRANSFER SUPPORT AND MOVING AID A patient support and moving aid, assisting and supporting when raising and sitting the patient down, transfer to and from bed. <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface for patient comfort and safety <input type="checkbox"/> Two pairs of handles on either end, giving a choice of grip Dimensions: Approx. 50x 20cm Packaging requirements: Single unit Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 89 | | PATIENT SUPPORT AND MOVING BELT | | | | | | | | | | | | | | | | | |
| ITEM-89.1 | 161164 | PATIENT SUPPORT AND MOVING BELT, Small Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles <input type="checkbox"/> Seven handles Dimensions: To accommodate waist sizes of 60 – 100cm Packaging requirements: Single unit Estimated usage: 1 | Yes, representative sample for items 89.1 to 89.4 | | Each | | | | | | | | | | | | | | |
| ITEM-89.2 | 161165 | PATIENT SUPPORT AND MOVING BELT, Medium Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles Nine handles Dimensions: To accommodate waist sizes of 72-120 cm Packaging requirements: Single unit Estimated usage: 55 | Yes, representative sample for items 89.1 to 89.4 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-89.3 | 161166 | PATIENT SUPPORT AND MOVING BELT, Large Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles Eleven handles Dimensions: To accommodate waist sizes of 100 – 160cm Packaging requirements: Single unit Estimated usage: 5 | Yes, representative sample for items 89.1 to 89.4 | | Each | | | | | | | | | | | | | | |
| ITEM-89.4 | 161167 | PATIENT SUPPORT AND MOVING BELT, X-Large Transfer aid to assist with standing, walking & transfers <input type="checkbox"/> Manufactured from polyester <input type="checkbox"/> Padded inner surface, for patient comfort and safety <input type="checkbox"/> Has both vertical and horizontal handles, offering a variety of grips <input type="checkbox"/> Adjustable strap and snap buckle, enabling a proper fit and tightening of the belt <input type="checkbox"/> Padded vertical and horizontal handles Thirteen handles Dimensions: To accommodate waist sizes of 115 – 175cm Packaging requirements: Single unit Estimated usage: 1 | Yes, representative sample for items 89.1 to 89.4 | | Each | Each | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|---------|---------------------|--|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 90 | Patient glide board | Patient glide board <input type="checkbox"/> For seating transfers <input type="checkbox"/> Manufactured from a Polyethylene and reinforced with non-rigid PVC <input type="checkbox"/> Flexible, allowing contouring to take place from different heights, during transfers from one surface to another <input type="checkbox"/> Elongated break allowing a side section to be bent upwards approx. 60 degree, to protect the patient during transfers <input type="checkbox"/> Anti-slip strips underneath <input type="checkbox"/> Handles at either end, for ease of use <input type="checkbox"/> Suitable for use during transfers, between bed and chair/wheelchair, between car and wheelchair and between wheelchair and toilet <input type="checkbox"/> Able to bridge gaps of approx. 10cm <input type="checkbox"/> Has no weight limit providing there is a stable surface, during transfer <input type="checkbox"/> Can be washed down, cleaned with alcohol or a surface disinfectant Dimensions: 75cm x 33cm Packaging requirements: Each Estimated usage: 1 | ITEM-90.1 | 161172 | Yes | Each | | | | | | | | | | | | | |
| ITEM 91 | Non-slip mat | Non-slip mat <input type="checkbox"/> Used on floor to prevent patients from slipping <input type="checkbox"/> Made of polyethylene Dimensions: Approx. 45 x 38cm Packaging requirements: Each Estimated usage: 0 | ITEM-91.1 | New | Yes | Each | | | | | | | | | | | | | |
| ITEM 92 | Transfer board | Transfer board, wooden <input type="checkbox"/> Tapered ends, light weight, smooth finish for easy sliding over board <input type="checkbox"/> Varnished Dimensions: 68cm x 30cm x 0.5cm Packaging requirements: Each Estimated usage: 0 | ITEM-92.1 | New | Yes | Each | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 93 | | Conventional vacuum erection device | | | | | | | | | | | | | | | | | |
| ITEM-93.1 | New | Conventional vacuum erection device □ Kit must include pump, quick seal tension rings of different sizes, service oil, personal lubrication Dimensions: Standard Packaging requirements: Each Estimated usage: 0 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 94 | | Wheelchair gloves (fingerless glove) | | | | | | | | | | | | | | | | | |
| ITEM-94.1 | New | Wheelchair gloves, heavy duty, extra small □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Extra Small Packaging requirements: per pair Estimated usage: 1 | Yes, representative sample for items 94.1 to 94.5 | | Per pair | | | | | | | | | | | | | | |
| ITEM-94.2 | New | Wheelchair gloves, heavy duty, small □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Small Packaging requirements: per pair Estimated usage: 1 | Yes, representative sample for items 94.1 to 94.5 | | Per pair | | | | | | | | | | | | | | |
| ITEM-94.3 | New | Wheelchair gloves, heavy duty, medium □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Medium Packaging requirements: per pair Estimated usage: 1 | Yes, representative sample for items 94.1 to 94.5 | | Per pair | | | | | | | | | | | | | | |
| ITEM-94.4 | New | Wheelchair gloves, heavy duty, large □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: Large Packaging requirements: per pair Estimated usage: 1 | Yes, representative sample for items 94.1 to 94.5 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-94.5 | New | Wheelchair gloves, heavy duty, X-Large <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove. Webbing covered with leather. Maximum impact protection. Dimensions: X-Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 94.1 to 94.5 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 95 | | Wheelchair gloves (fully covered) | | | | | | | | | | | | | | | | | |
| ITEM-95.1 | New | Wheelchair gloves, heavy duty, extra small □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Extra Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 95.1 to 95.5 | | Per pair | | | | | | | | | | | | | | |
| ITEM-95.2 | New | Wheelchair gloves, heavy duty, small □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 95.1 to 95.5 | | Per pair | | | | | | | | | | | | | | |
| ITEM-95.3 | New | Wheelchair gloves, heavy duty, medium □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 95.1 to 95.5 | | Per pair | | | | | | | | | | | | | | |
| ITEM-95.4 | New | Wheelchair gloves, heavy duty, large □ Durable foam padded leather on wrist and palm □ Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 95.1 to 95.5 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-----------|------------------------------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-95.5 | New | Wheelchair gloves, heavy duty, X-Large <input type="checkbox"/> Durable foam padded leather on wrist and palm <input type="checkbox"/> Breathable spandex on the back of the hand. Visco - elastic polymer pad over the palm and between thumb and index finger. Fingerless glove with thumb fully covered. Webbing covered with leather. Maximum impact protection Dimensions: X-Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 95.1 to 95.5 | | Per pair | | | | | | | | | | | | | | |
| ITEM 96 | Universal cuff | | | | | | | | | | | | | | | | | | |
| ITEM-96.1 | New | Universal Cuff utensil pocket <input type="checkbox"/> Made out of webbing material <input type="checkbox"/> Used to help patients with limited hand function to grasp eating utensils Has hook and loop closure on a D-Ring Includes utensil pocket Dimensions: can be adjusted to fit palms measuring 18cm to 23cm in circumference Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-96.2 | New | Universal Cuff cylindrical foam <input type="checkbox"/> Made out of webbing material <input type="checkbox"/> Used to help patients with limited hand function to grasp eating utensils Has hook and loop closure on a D-Ring Includes cylindrical foam to hold utensils Dimensions: can be adjusted to fit palms measuring 18cm to 23cm in circumference Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 97 | Contact adhesive glue | | | | | | | | | | | | | | | | | | |
| ITEM-97.1 | New | Contact adhesive glue <input type="checkbox"/> Non - staining, quick drying, versatile on different material <input type="checkbox"/> Similar or equal to Genkem Dimensions: 500ml tin Packaging requirements: Each Estimated usage: 0 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 98 | Patient turn table | | | | | | | | | | | | | | | | | | |
| ITEM-98.1 | 161674 | PATIENT TURN TABLE Diameter. 38cm & 16mm thick <input type="checkbox"/> Used on floor to assist with turning patients that have difficulty in turning, or moving their feet. <input type="checkbox"/> Two discs manufactured from Polyethylene, fixed together in the centre on a axle. Bottom disc has non slip surface which prevents the patient from moving or slipping on it, Top disc rotates on axle, while bottom disc stays static. <input type="checkbox"/> Has carrying handle, <input type="checkbox"/> Easy to clean Dimensions: Diameter of 38cm, Thickness 12 mm with handle Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 99 | | Cup specialized feeding | | | | | | | | | | | | | | | | | |
| ITEM-99.1 | 107265 | <p>Cup specialized feeding transparent/semi</p> <p>Cup, specialized feeding, should be transparent or semi-transparent with no sharp edge and the surface should be smooth, made of high grade polypropylene and be Bisphenol A (toxin) free., must be able to withstand 125°C temp, must be calibrated in 5ml and must hold 50-60ml, must be autoclavable, the inner and outer surface of the cup must be reachable by hand and/or finger in order to clean, "Sinapi" equal or better.</p> <p>Packaging requirements: Per cup</p> <p>Estimated usage: 4,500 per year</p> | Yes | | Per cup | | | | | | | | | | | | | | |
| ITEM 100 | | Cup feeding system | | | | | | | | | | | | | | | | | |
| ITEM-100.1 | 107270 | <p>Cup feeding system</p> <p>Cup feeding system, consists of main holder, lid for the main holder and the feeding cup, main holder and feeding cup should be from virgin high-grade polypropylene material, be able to withstand a temperature of 120°C, the main holder must hold 200ml and have volume calibrations marked every 25ml, the feeding cup must hold 50 - 60ml, the feeding cup must be able to be stored inside the main holder with the lid in place, have no sharp edges and the surface texture should be smooth, the lid must fit watertight on the main holder, the inner and outer surface of the different components must be reachable by hand and/or finger in order to clean, the unit must be individually packaged and sealed in a plastic bag with an illustrated leaflet printed, "Sinapi" equal or better.</p> <p>Packaging requirements: Per system</p> <p>Estimated usage: 6,000 per year</p> | Yes | | Per system | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 101 | | ITEM 101 Item 101: Breastmilk storage container, disposable | | | | | | | | | | | | | | | | | |
| ITEM-101.1 | New | Breastmilk storage container, disposable, 120-125ml Breastmilk storage container, disposable, to store expressed breast milk. Manufactured from high density food grade polyethylene. BPA (Bisphenol-A) and DEHP (Di-Ethylhexyl Phthalate) free. Include foil cap. Include snap on cap made of polypropylene or polyethylene (water tight seal). Clean, filtered air used during blow moulding. Withstand pasteurisation up to 90°C. Size: 120-125ml capacity Packaging requirements: Each Estimated usage: 5,000 per year | Yes, representative sample for items 101.1 and 101.2 | | Each | | | | | | | | | | | | | | |
| ITEM-101.2 | New | Breastmilk storage container, disposable, 250ml Breastmilk storage container, disposable, to store expressed breast milk. Manufactured from high density food grade polyethylene. BPA (Bisphenol-A) and DEHP (Di-Ethylhexyl Phthalate) free. Include foil cap. Include snap on cap made of polypropylene or polyethylene (water tight seal). Clean, filtered air used during blow moulding. Withstand pasteurisation up to 90°C. Size: 250ml capacity Packaging requirements: Each Estimated usage: 10,000 per year | Yes, representative sample for items 101.1 and 101.2 | | Each | | | | | | | | | | | | | | |
| ITEM 102 | | C:- PODIATRY | | | | | | | | | | | | | | | | | |
| ITEM-102.1 | 159280 | ITEM 102 Item 102: Foam-o-felt Foam-o-felt, 45x22.5cmx5mm Foam-o-felt, mixture of 70% wool and 30% viscose (good for shock absorption when natural fibrofatty pad has worn out), 45cm minimum length, 22.5cm width, 5mm thickness. Dimensions: 45cm x 22.5cm x 5mm Packaging requirements: Pack of 4 Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|------------------------|--|--|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 103 | Poron | | | | | | | | | | | | | | | | | | |
| ITEM-103.1 | 159282 | Poron, 31.5x22.5cmx3.2mm Poron, soft supporting material, micro porous polyurethane with shock-absorbing properties, medical grade, 3.2mm thickness, 31.5cm length, 22.5cm width, per sheet, pack of 4 Dimensions: 31.5cm x 22.5cm x 3.2mm Packaging requirements: Pack of 4 Estimated usage: 72 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 104 | PPT ultralux open cell | | | | | | | | | | | | | | | | | | |
| ITEM-104-1 | 159287 | PPT ultralux open cell, 152x30cmx1.6mm PPT ultralux, durable open cell material that doesn't bottom out, 152cm length, 30cm width, 1.6mm thickness, good shock absorber, ideal covering material Dimensions: 152cm x 30cm x 1.6mm Packaging requirements: Pack of 2 Estimated usage: 1 | Yes | | Per pack | | | | | | | | | | | | | | |
| ITEM 105 | Cortex cork material | | | | | | | | | | | | | | | | | | |
| ITEM-105.1 | 159150 | Cortex cork material, 1mx600x1.33mm Cortex, cork material, 1m length per sheet, 600mm width, 1.33mm thickness, used as base for corrective foot devices Dimensions: 1m x 600mm x 1.33mm Packaging requirements: Per sheet Estimated usage: 1 | Yes, representative sample for items 105.1 and 105.2 | | Per sheet | | | | | | | | | | | | | | |
| ITEM-105.2 | 159151 | Cortex cork material, 1mx600x1.5mm Cortex, cork material, 1m length per sheet, 600mm width, 1.5mm thickness, used as base for corrective foot devices Dimensions: 1m x 600mm x 1.5mm Packaging requirements: Per sheet Estimated usage: 50 | Yes, representative sample for items 105.1 and 105.2 | | Per sheet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 106 | | Neoprene material | | | | | | | | | | | | | | | | | |
| ITEM-106.1 | 178709 | Neoprene material, 1mx3mm Neoprene material 1m length, 3mm thick Dimensions: 1m x 0,5cm Packaging requirements: Per sheet Estimated usage: 1 | Yes | | Per sheet | | | | | | | | | | | | | | |
| ITEM 107 | | Polyethylene, mould plastic | | | | | | | | | | | | | | | | | |
| ITEM-107.1 | 159152 | Polyethylene, mould plastic, 1-2mm thick Polyethylene, heat mouldable plastic, 1-2mm thickness, 1.2x1m length Dimensions: 1-2mm Packaging requirements: Per sheet Estimated usage: 1 | Yes | | Per sheet | | | | | | | | | | | | | | |
| ITEM 108 | | Polypropylene, heat mould plastic | | | | | | | | | | | | | | | | | |
| ITEM-108.1 | 159153 | Polypropylene, heat mould plastic, 1.2x1mx2mm Polypropylene, heat mouldable plastic, 2mm thickness and above, 1.2x1m length Dimensions: 1.2 x 1m x 2mm Packaging requirements: Per sheet Estimated usage: 1 | Yes | | Per sheet | | | | | | | | | | | | | | |
| ITEM 109 | | Adhesive contact liquid | | | | | | | | | | | | | | | | | |
| ITEM-109.1 | 136798 | Adhesive contact liquid 1L Adhesive contact, used for laminating all components of corrective foot devices, liquid, 1L Size: 1L Packaging requirements: Per 1L Estimated usage: 1 | Yes, representative sample for items 109.1 and 109.2 | | Per L | | | | | | | | | | | | | | |
| ITEM-109.2 | 159162 | Adhesive contact liquid 5L Adhesive contact, used for laminating all components of corrective foot devices, liquid, 5L Size: 5L Packaging requirements: Per 5L Estimated usage: 10 | Yes, representative sample for items 109.1 and 109.2 | | Per L | | | | | | | | | | | | | | |
| ITEM 110 | | Remover solvent | | | | | | | | | | | | | | | | | |
| ITEM-110.1 | 159163 | Remover solvent to dissolve adhesive, 5L Remover solvent, chemical to dissolve contact adhesive, 5L Size: 5L Packaging requirements: Per 5L Estimated usage: 10 | Yes | | Per L | | | | | | | | | | | | | | |
| ITEM 111 | | Fleecy web, cotton | | | | | | | | | | | | | | | | | |
| ITEM-111.1 | 159283 | Fleecy web, cotton, 40 x 22.5cm Fleecy web, 100% cotton padding with a raised fleecy surface, flexible application, 40cm length, 22.5cm width, pack of 4 Dimensions: 40 x 22.5cm Packaging requirements: Pack of 4 Estimated usage: 1 | Yes | | per pack | | | | | | | | | | | | | | |
| ITEM 112 | | Felt, semi-compressed | | | | | | | | | | | | | | | | | |
| ITEM-112.1 | 114120 | Felt, semi-compressed, 3mm Felt, semi-compressed, 100% pure wool, used as a highly versatile padding which can be easily cut and shaped, providing excellent pressure absorption. Size: 3mm Packaging requirements: Pack of 4 Estimated usage: 1 | Yes, representative sample for items 112.1 and 112.2 | | per pack | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|---|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-112.2 | 114121 | Felt, semi-compressed, 5mm Felt, semi-compressed, 100% pure wool, used as a highly versatile padding which can be easily cut and shaped, providing excellent pressure absorption. Size: 5mm Packaging requirements: Pack of 4 Estimated usage: 1 | Yes, representative sample for items 112.1 and 112.2 | | Per pack | | | | | | | | | | | | | | |
| ITEM 113 Tubifoam sponge tubular | | | | | | | | | | | | | | | | | | | |
| ITEM-113.1 | 159284 | Tubifoam, sponge tubular, 18mmx25cm Tubifoam, sponge, tubular, 18mm width, 25cm length, requires no strapping or tape, protects toes from friction, per tube (pack of 10 tubes) Dimensions: 18mm x 25cm Packaging requirements: Pack of 10 tubes Estimated usage: 1 | Yes, representative sample for items 113.1 to 113.3 | | Per pack | | | | | | | | | | | | | | |
| ITEM 113.2 Tubifoam, sponge tubular, 21mmx25cm | | | | | | | | | | | | | | | | | | | |
| ITEM-113.2 | 159285 | Tubifoam, sponge, tubular 21mm width, 25cm length, requires no strapping or tape, protects toes from friction, per tube (pack of 10 tubes) Dimensions: 21mm x 25cm Packaging requirements: Pack of 10 tubes Estimated usage: 1 | Yes, representative sample for items 113.1 to 113.3 | | Per pack | | | | | | | | | | | | | | |
| ITEM-113.3 | 159286 | Tubifoam, sponge tubular, 25mmx25cm Tubifoam, sponge, tubular, 25mm width, 25cm length, requires no strapping or tape, protects toes from friction, per tube (pack of 10 tubes) Dimensions: 25mm x 25cm Packaging requirements: Pack of 10 tubes Estimated usage: 1 | Yes, representative sample for items 113.1 to 113.3 | | Per pack | | | | | | | | | | | | | | |
| ITEM 114 Heel lifts, PPT pad | | | | | | | | | | | | | | | | | | | |
| ITEM-114.1 | 136776 | Heel lifts, PPT pad, small Heel lifts, PPT pad indicated for heel spurs, plantar fasciitis, very soft density to provide relief to sensitive areas. Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 114 to 114.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-114.2 | 136777 | Heel lifts, PPT pad, medium Heel lifts, PPT pad indicated for heel spurs, plantar fasciitis, very soft density to provide relief to sensitive areas. Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 114 to 114.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-114.3 | 136778 | Heel lifts, PPT pad, large Heel lifts, PPT pad indicated for heel spurs, plantar fasciitis, very soft density to provide relief to sensitive areas. Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 114 to 114.3 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|---|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 115 Heel raise | | | | | | | | | | | | | | | | | | | |
| ITEM-115.1 | 159288 | Heel raise, small Heel raise, indicated for leg length difference and heel pain, polyurethane, density 0.25-0.28 grams/cm ³ , hardness of 15-20° shore A, covering material is micro fibre material. Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 115.1 to 115.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-115.2 | 159289 | Heel raise, medium Heel raise, indicated for leg length difference and heel pain, polyurethane, density 0.25-0.28 grams/cm ³ , hardness of 15-20° shore A, covering material is micro fibre material. Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 115.1 to 115.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-115.3 | 159300 | Heel raise, large Heel raise, indicated for leg length difference and heel pain, polyurethane, density 0.25-0.28 grams/cm ³ , hardness of 15-20° shore A, covering material is micro fibre material. Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 115.1 to 115.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 116 Heel pad, horse shoe | | | | | | | | | | | | | | | | | | | |
| ITEM-116.1 | 114147 | Heel pad, horse shoe, small Heel pad, horse shoe, indicated for off loading, latex sponge rubber, 8-10mm, with density of 0.25 - 0.28 grams/cm ³ , hardness of 15-20° shore A. Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 116.1 to 116.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-116.2 | 114148 | Heel pad, horse shoe, medium Heel pad, horse shoe, indicated for off loading, latex sponge rubber, 8-10mm, with density of 0.25 - 0.28 grams/cm ³ , hardness of 15-20° shore A. Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 116.1 to 116.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-116.3 | 114149 | Heel pad, horse shoe, large Heel pad, horse shoe, indicated for off loading, latex sponge rubber, 8-10mm, with density of 0.25 - 0.28 grams/cm ³ , hardness of 15-20° shore A. Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 116.1 to 116.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 117 Heel pad horse shoe with extractable central core | | | | | | | | | | | | | | | | | | | |
| ITEM-117.1 | 114150 | Heel pad horse shoe with extractable central core Heel pad, horse shoe, with extractable central core to relieved calcaneal spurs, indicated for heel spurs and plantar fasciitis and off loading, latex sponge rubber, 5-10mm, with density of 0.25 - 0.28 grams/cm ³ , hardness of 15-20° shore A. Size: One size fits all Packaging requirements: Per pair Estimated usage: 1 | Yes | | Per pair | | | | | | | | | | | | | | |
| ITEM 118 Cobra pad with valgus pad extension | | | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------------------------------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-118.1 | 114153 | Cobra pad with valgus pad extension, small Cobra pad with valgus pad extension, sponge rubber, density of 0.25 - 0.28 grams/cm ³ , hardness of 15-20° shore A. Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 118.1 to 118.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-118.2 | 114152 | Cobra pad with valgus pad extension, medium Cobra pad with valgus pad extension, sponge rubber, density of 0.25 - 0.28 grams/cm ³ , hardness of 15-20° shore A. Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 118.1 to 118.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-118.3 | 114151 | Cobra pad with valgus pad extension, large Cobra pad with valgus pad extension, sponge rubber, density of 0.25 - 0.28 grams/cm ³ , hardness of 15-20° shore A. Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 118.1 to 118.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 119 | Metatarsal bar, PPT pad | | | | | | | | | | | | | | | | | | |
| ITEM-119.1 | 159301 | Metatarsal bar, PPT pad, small Metatarsal bar, PPT pad, small, allowing for extra protection under metatarsal heads (prevents growth of hard callouses and reduces pain and irritation). Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 119.1 to 119.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-119.2 | 159302 | Metatarsal bar, PPT pad, medium Metatarsal bar, PPT pad, medium, allowing for extra protection under metatarsal heads (prevent growth of hard areas and callouses and reduces pain and irritation). Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 119.1 to 119.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-119.3 | 159303 | Metatarsal bar, PPT pad, large Metatarsal bar, PPT pad, large, allowing for extra protection under metatarsal heads (prevents growth of hard areas and callouses and reduces pain and irritation) Size: Large Packaging requirements: Per pair Estimated usage: 50 | Yes, representative sample for items 119.1 to 119.3 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|---|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 120 | Metatarsal bar, synthetic latex pad with cut-out | | | | | | | | | | | | | | | | | | |
| ITEM-120.1 | 114122 | Metatarsal bar, synthetic latex pad with cut-out, small Metatarsal bar, synthetic latex pad, small, with a cut-out on the first and fifth metatarsal, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm ³ Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 120.1 to 120.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-120.2 | 114123 | Metatarsal bar, synthetic latex pad with cut-out, medium Metatarsal bar, synthetic latex pad, medium, with a cut-out on the first and fifth metatarsal, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm ³ Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 120.1 to 120.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-120.3 | 114124 | Metatarsal bar, synthetic latex pad with cut-out, large Metatarsal bar, synthetic latex pad, large, with a cut-out on the first and fifth metatarsal, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm ³ Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 120.1 to 120.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 121 | Metatarsal bar, synthetic latex pad, straight | | | | | | | | | | | | | | | | | | |
| ITEM-121.1 | 114125 | Metatarsal bar, synthetic latex pad, straight, small Metatarsal bar, synthetic latex pad, small, straight across metatarsal heads, with no cut-out on the first and fifth metatarsal, med pad made of synthetic latex foam, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm ³ Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 121.1 and 121.2 | | Per pair | | | | | | | | | | | | | | |
| ITEM-121.2 | 114126 | Metatarsal bar, synthetic latex pad, straight, large Metatarsal bar, synthetic latex pad, large, straight across metatarsal heads, with no cut-out on the first and fifth metatarsal, med pad made of synthetic latex foam, hardness 15-20° shore A, density 0.25 / 0.28 grams/cm ³ Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 121.1 and 121.2 | | Per pair | | | | | | | | | | | | | | |
| ITEM 122 | Metatarsal pad with arch support, synthetic latex | | | | | | | | | | | | | | | | | | |
| ITEM-122.1 | 114142 | Metatarsal pad with arch support, synthetic latex, small Metatarsal pad with arch support, made of synthetic latex, with hardness of 15-20° shore A, density of 0.25gram/cm ³ Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 122.1 to 122.3 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-122.2 | 114143 | Metatarsal pad with arch support, synthetic latex, medium Metatarsal pad with arch support, made of synthetic latex, with hardness of 15-20° shore A, density of 0.25gram/cm ³ Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 122.1 to 122.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-122.3 | 114144 | Metatarsal pad with arch support, synthetic latex, large Metatarsal pad with arch support, made of synthetic latex, with hardness of 15-20° shore A, density of 0.25gram/cm ³ Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 122.1 to 122.3 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 123 | | Pad, metatarsal, foam rubber | | | | | | | | | | | | | | | | | |
| ITEM-123.1 | 159159 | Pad, metatarsal, foam rubber, small Pad, metatarsal, soft foam rubber, made of foam rubber for its supportive and cushioning properties, small Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 123.1 to 123.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-123.2 | 159160 | Pad, metatarsal foam, rubber, medium Pad, metatarsal, soft foam rubber, made of foam rubber for its supportive and cushioning properties, small Size: Medium Packaging requirements: Per pair Estimated usage: 6 | Yes, representative sample for items 123.1 to 123.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-123.3 | 159161 | Pad, metatarsal foam, rubber, large Pad, metatarsal, soft foam rubber, made of foam rubber for its supportive and cushioning properties, small Size: Large Packaging requirements: Per pair Estimated usage: 6 | Yes, representative sample for items 123.1 to 123.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 124 | | Pads, valgus arch support, EVA | | | | | | | | | | | | | | | | | |
| ITEM-124.1 | 159156 | Pads, valgus arch support, EVA, small Pads valgus, arch support, made of EVA of soft density sponge rubber, small Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 124.1 to 124.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-124.2 | 159157 | Pads, valgus arch support, EVA, medium Pads valgus, arch support, made of EVA of soft density sponge rubber, small Size: Medium Packaging requirements: Per pair Estimated usage: 6 | Yes, representative sample for items 124.1 to 124.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-124.3 | 159158 | Pads, valgus arch support, EVA, large Pads valgus, arch support, made of EVA of soft density sponge rubber, small Size: Large Packaging requirements: Per pair Estimated usage: 6 | Yes, representative sample for items 124.1 to 124.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 125 | | Pads, valgus arch support, synthetic latex, hard | | | | | | | | | | | | | | | | | |
| ITEM-125.1 | 114139 | Pads, valgus arch support, synthetic latex, hard, small Pads, valgus arch support, made of synthetic latex foam, hard, density of 30-35° shore, hardness 1 grams/cm ³ Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 125.1 to 125.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-125.2 | 114140 | Pads, valgus arch support, synthetic latex, hard, medium Pads, valgus arch support, made of synthetic latex foam, hard, density of 30-35° shore, hardness 1 grams/cm ³ Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 125.1 to 125.3 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-125.3 | 114138 | Pads, valgus arch support, synthetic latex, hard, large Pads, valgus arch support, made of synthetic latex foam, hard, density of 30-35° shore, hardness 1 grams/cm ³ Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 125.1 to 125.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-126 | | Pads, valgus arch support, synthetic latex, soft | | | | | | | | | | | | | | | | | |
| ITEM-126.1 | 114136 | Pads, valgus arch support, synthetic latex, soft, small Pads, valgus arch support, made of synthetic latex foam, soft, density between 15-20° shore Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 126.1 to 126.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-126.2 | 114137 | Pads, valgus arch support, synthetic latex, soft, medium Pads, valgus arch support, made of synthetic latex foam, soft, density between 15-20° shore Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 126.1 to 126.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-126.3 | 114141 | Pads, valgus arch support, synthetic latex, soft, large Pads, valgus arch support, made of synthetic latex foam, soft, density between 15-20° shore Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 126.1 to 126.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-127 | | Digi-tube, ribbed elastic, strip | | | | | | | | | | | | | | | | | |
| ITEM-127.1 | 159305 | Digi-tube, ribbed, elastic strip, small Digi-tube, ribbed elastic material, tubular strip, fully lined with silicone, small, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) Size: Small Packaging requirements: Each Estimated usage: 30 | Yes, representative sample for items 127.1 to 127.3 | | Each | | | | | | | | | | | | | | |
| ITEM-127.2 | 159307 | Digi-tube, ribbed, elastic strip, medium Digi-tube, ribbed elastic material, tubular strip, fully lined with silicone, medium, durable, washable which can be cut into different lengths. Size: Medium Packaging requirements: Each Estimated usage: 20 | Yes, representative sample for items 127.1 to 127.3 | | Each | | | | | | | | | | | | | | |
| ITEM-127.3 | 159308 | Digi-tube, ribbed elastic, strip, large Digi-tube, ribbed elastic material, tubular strip, fully lined with silicone, large, durable, washable which can be cut into different lengths. Size: Large Packaging requirements: Each Estimated usage: 8 | Yes, representative sample for items 127.1 to 127.3 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|----------------------------|---|-----------------|--|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 128 | Digi-tube, pure gel, strip | Digi-tube, pure gel, strip, small Digi-tube, pure gel, tubular strip, small, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) Size: Small Packaging requirements: Each Estimated usage: New item | 136779 | Yes, representative sample for items 128.1 to 128.3 | Each | | | | | | | | | | | | | | |
| ITEM-128.1 | | | | | | | | | | | | | | | | | | | |
| ITEM-128.2 | 136780 | Digi-tube, pure gel, strip, medium Digi-tube, pure gel, tubular strip, medium, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) Size: Medium Packaging requirements: Each Estimated usage: New item | | Yes, representative sample for items 128.1 to 128.3 | Each | | | | | | | | | | | | | | |
| ITEM-128.3 | 136781 | Digi-tube, pure gel, strip, large Digi-tube, pure gel, tubular strip, large, durable, washable which can be cut into different lengths, (relief pressure and friction on the toes) Size: Large Packaging requirements: Each Estimated usage: New item | | Yes, representative sample for items 128.1 to 128.3 | Each | | | | | | | | | | | | | | |
| ITEM 129 | Toe loop | Toe loop, narrow Toe loop, narrow loop, foam-lined material with velcro strap, ideal for repositioning the toes, strap with loops secures two or more toes together to stabilise deformed toes Size: Narrow loop, 1.5 x 7.6cm Packaging requirements: Each Estimated usage: 1 | 114157 | Yes, representative sample for items 129.1 and 129.2 | Each | | | | | | | | | | | | | | |
| ITEM-129.1 | | | | | | | | | | | | | | | | | | | |
| ITEM-129.2 | 114158 | Toe loop, wider Toe loop, wider loop, foam-lined material with velcro strap, ideal for repositioning the toes, strap with loops secures two or more toes together to stabilise deformed toes Size: Wider loop, 4 x 10cm Packaging requirements: Each Estimated usage: 1 | | Yes, representative sample for items 129.1 and 129.2 | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 130 | | Toe bar | | | | | | | | | | | | | | | | | |
| ITEM-130.1 | 114154 | Toe bar, small Toe bar, with adjustable elastic band that fits around toes, made of chamois leather, recommended for hammer toes. Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 130.1 to 130.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-130.2 | 114155 | Toe bar, medium Toe bar, with adjustable elastic band that fits around toes, made of chamois leather, recommended for hammer toes. Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 130.1 to 130.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-130.3 | 114156 | Toe bar, large Toe bar, with adjustable elastic band that fits around toes, made of chamois leather, recommended for hammer toes. Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 130.1 to 130.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 131 | | CAP digital, pure gel | | | | | | | | | | | | | | | | | |
| ITEM-131.1 | 159720 | CAP digital, pure gel, small Cap digital, pure gel, small, eases pressure and friction caused by blisters, corns and it is closed in front, transparent, durable and washable, single Size: Small Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 131.1 to 131.3 | | Each | | | | | | | | | | | | | | |
| ITEM-131.2 | 159719 | CAP digital, pure gel, medium Cap digital, pure gel, medium, eases pressure and friction caused by blisters, corns and it is closed in front, transparent, durable and washable, single Size: Medium Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 131.1 to 131.3 | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|---|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-131.3 | 159722 | CAP digital, pure gel, large Cap digital, pure gel, large, eases pressure and friction caused by blisters, corns and it is closed in front, transparent, durable and washable, single Size: Large Packaging requirements: Per item Estimated usage: 1 | Yes, representative sample for items 131.1 to 131.3 | | Each | | | | | | | | | | | | | | |
| ITEM 132 Spreader toe, pure gel interdigital wedge | | | | | | | | | | | | | | | | | | | |
| ITEM-132.1 | 159797 | Spreader toe, pure gel interdigital wedge, small Spreader toe, pure gel interdigital wedge, small, for deformities and overlapping toes, protection in helomas and interdigital pressure, helps with re-alignment of toes Size: Small Packaging item: Per pack of 12 Estimated usage: 1 | Yes, representative sample for items 132.1 to 132.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-132.2 | 159798 | Spreader toe, pure gel interdigital wedge, medium Spreader toe, pure gel interdigital wedge, medium, for deformities and overlapping toes, protection in helomas and interdigital pressure, helps with re-alignment of toes. Size: Medium Packaging requirements: Per pack of 12 Estimated usage: 1 | Yes, representative sample for items 132.1 to 132.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-132.3 | 159800 | Spreader toe, pure gel interdigital wedge, large Spreader toe, pure gel interdigital wedge, large, for deformities and overlapping toes, protection in heloma and interdigital, pressure, helps with re-alignment of toes. Size: Large Packaging requirements: Per pack of 12 Estimated usage: 1 | Yes, representative sample for items 132.1 to 132.3 | | Per packet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING /BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|-------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 133 | | Spreader, hammer, made of soft polymer gel, big toe spreader | | | | | | | | | | | | | | | | | |
| ITEM-133.1 | 137019 | Spreader, hammer, made of soft polymer gel, big toe spreader, small Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Small Packaging requirements: Pack of 12 units Estimated usage: 6 | Yes, representative sample for items 133.1 to 133.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-133.2 | 137020 | Spreader, hammer, made of soft polymer gel, big toe spreader, medium Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Medium Packaging requirements: Pack of 12 units Estimated usage: 6 | Yes, representative sample for items 133.1 to 133.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-133.3 | 170844 | Spreader, hammer, made of soft polymer gel, big toe spreader, large Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Large Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 133.1 to 133.3 | | Per packet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 134 | | Spreader, hammer, made of soft polymer gel, advanced toe spreader | | | | | | | | | | | | | | | | | |
| ITEM-134.1 | 170848 | Spreader, hammer, made of soft polymer gel, advanced toe spreader, small Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Small Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 134.1 to 134.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-134.2 | 170849 | Spreader, hammer, made of soft polymer gel, advanced toe spreader, medium Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Medium Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 134.1 to 134.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-134.3 | 170860 | Spreader, hammer, made of soft polymer gel, advanced toe spreader, large Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Large Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 134.1 to 134.3 | | Per packet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 135 | | Spreader, hammer, made of soft polymer gel, universal toe spreader | | | | | | | | | | | | | | | | | |
| ITEM-135.1 | 170845 | Spreader, hammer, made of soft polymer gel, universal toe spreader, small Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Small Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 135.1 to 135.4 | | Per packet | | | | | | | | | | | | | | |
| ITEM-135.2 | 170846 | Spreader, hammer, made of soft polymer gel, universal toe spreader, medium Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Medium Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 135.1 to 135.4 | | Per packet | | | | | | | | | | | | | | |
| ITEM-135.3 | 170847 | Spreader, hammer, made of soft polymer gel, universal toe spreader, large Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: Large Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 135.1 to 135.4 | | Per packet | | | | | | | | | | | | | | |
| ITEM-135.4 | 170839 | Spreader, hammer, made of soft polymer gel, universal toe spreader, X-Large Stretch, for easily spread crooked, overlapping and rubbing toes, to reduce friction related injuries such as blisters, corns and calluses, ±2.5cm length, 12 units Size: X-Large Packaging requirements: Pack of 12 units Estimated usage: 1 | Yes, representative sample for items 135.1 to 135.4 | | Per packet | | | | | | | | | | | | | | |
| ITEM 136 | | Pad, hammer, polymer gel | | | | | | | | | | | | | | | | | |
| ITEM-136.1 | 137021 | Pad, hammer, polymer gel Hammer toes crest pad, made of polymer gel, soft, good stretch, relieves forefoot and toe pain caused by hammer, claw and mallet toes, toe crest provides maximum comfort and reduces pressure from the toes and metatarsals Packaging requirements: Pack of 10 units Estimated usage: 1 | Yes | | Per packet | | | | | | | | | | | | | | |
| ITEM 137 | | Separator toe, gel mould fitting | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-137.1 | 159803 | Separator toe, gel mould fitting, small Separator toe, gel mould fitting, small, to reduce interdigital pressure between toes, protection from inter-digital pressure, helps with re-alignment of toes. Size: Small Packaging requirements: Per pack of 12 Estimated usage: 1 | Yes, representative sample for items 137.1 to 137.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-137.2 | 178702 | Separator toe, gel mould fitting, medium Separator toe, gel mould fitting, medium, to reduce interdigital pressure between toes, protection from inter-digital pressure, helps with re-alignment of toes. Size: Medium Packaging requirements: Per pack of 12 Estimated usage: 1 | Yes, representative sample for items 137.1 to 137.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-137.3 | 159805 | Separator toe, gel mould fitting, large Separator toe, gel mould fitting, large, to reduce interdigital pressure between toes, protection from inter-digital pressure, helps with re-alignment of toes. Size: Large Packaging requirements: Per pack of 12 Estimated usage: 1 | Yes, representative sample for items 137.1 to 137.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM 138 | 159712 | Shield bunion, foam shell, gel disc, small Shield bunion, toes, soft foam shell for pressure absorption and gel disc for pressure positioning, small, protects the bunion joint. Reduces pain in the bone protuberance, thin silicone for better patient compliance. Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 138.1 to 138.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-138.2 | 159713 | Shield bunion, foam shell, gel disc, medium Shield bunion, toes, soft foam shell for pressure absorption and gel disc for pressure positioning, medium, protects the bunion joint. Reduces pain in bone protuberance, thin silicone for patient compliance. Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 138.1 to 138.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-138.3 | 159714 | Shield bunion, foam shell, gel disc, large Shield bunion, toes, soft foam shell for pressure absorption and gel disc for pressure positioning, large, protects the bunion joint. Reduces pain in bone protuberance, thin silicone for better patient compliance. Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 138.1 to 138.3 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|-----------------------------------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 139 | Shield, bunion, polymer gel, soft | Shield, bunion, polymer gel, soft, small Shield, bunion, toes, made of soft polymer gel, good stretch, in 2 levels of density, protector is slim, yet very effective and can be used with any shoe, protectors help guard a painful bunion from rubbing and pressure and for protection with loop around big toe. Dimensions: Length ± 2.5cm Size: Small Packaging requirements: Pack of 6 units Estimated usage: 1 | Yes, representative sample for items 139.1 to 139.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-139.1 | 170827 | Shield, bunion, polymer gel, soft, medium Shield, bunion, toes, made of soft polymer gel, good stretch, in 2 levels of density, protector is slim, yet very effective and can be used with any shoe, protectors help guard a painful bunion from rubbing and pressure and for protection with loop around big toe. Dimensions: Length ± 2.5cm Size: Medium Packaging requirements: Pack of 6 units Estimated usage: 1 | Yes, representative sample for items 139.1 to 139.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-139.2 | 137066 | Shield, bunion, polymer gel, soft, medium Shield, bunion, toes, made of soft polymer gel, good stretch, in 2 levels of density, protector is slim, yet very effective and can be used with any shoe, protectors help guard a painful bunion from rubbing and pressure and for protection with loop around big toe. Dimensions: Length ± 2.5cm Size: Medium Packaging requirements: Pack of 6 units Estimated usage: 1 | Yes, representative sample for items 139.1 to 139.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-139.3 | 170829 | Shield, bunion, polymer gel, soft, large Shield, bunion, toes, made of <u>soft</u> polymer gel, good stretch, in 2 levels of density, protector is slim, yet very effective and can be used with any shoe, protectors help guard a painful bunion from rubbing and pressure and for protection with loop around big toe. Dimensions: Length ± 2.5cm Size: Large Packaging requirements: Pack of 6 units Estimated usage: 1 | Yes, representative sample for items 139.1 to 139.3 | | Per packet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------------------------------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 140 | Shield bunion, silicone | | | | | | | | | | | | | | | | | | |
| ITEM-140.1 | 178468 | Shield bunion, silicone, small Shield bunion toes layer of silicone gel placed to protect prominent big toe joint, covered with elastic material to keep silicone gel in place, reduce pain in bone protuberance to relief pressure Size: Small Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 140.1 to 140.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-140.2 | 178470 | Shield bunion, silicone, medium Shield bunion toes layer of silicone gel placed to protect prominent big toe joint, covered with elastic material to keep silicone gel in place, reduce pain in bone protuberance to relief pressure Size: Medium Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 140.1 to 140.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM-140.3 | 178472 | Shield bunion, silicone, large Shield bunion toes layer of silicone gel placed to protect prominent big toe joint, covered with elastic material to keep silicone gel in place, reduce pain in bone protuberance to relief pressure Size: Large Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 140.1 to 140.3 | | Per pair | | | | | | | | | | | | | | |
| ITEM 141 | Splint, Bunion / Hallux Valgus | | | | | | | | | | | | | | | | | | |
| ITEM-141.1 | 114145 | Splint, Bunion / Hallux Valgus Splint, for correction of bunion / hallux valgus, keeps first radius of foot in correct position, thin enough to wear in shoes, non-slip, breathable foam, holds the correction even during weight-bearing. Size: One size fits all Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|----------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 142 | ITEM 142 | Arch lift plantar fasciitis | | | | | | | | | | | | | | | | | |
| ITEM-142.1 | 114146 | Arch lift plantar fasciitis Arch lift, for plantar fasciitis, foam lined material with velcro straps, alleviates symptoms of plantar fasciitis and heel spurs, provides upward compression to the arch, supporting and stabilising the plantar fascia ligaments and includes a built-in pad for additional support to medial arch, helping to prevent excessive pronation. Size: One size fits all Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 143 | ITEM 143 | CAP digital, polymer gel, soft | | | | | | | | | | | | | | | | | |
| ITEM-143.1 | 137063 | CAP digital, polymer gel, soft, small CAP digital textile fabric inside coated with polymer gel, <u>soft</u> , stretchable, protects toes against rubbing and pressure, transparent, durable, washable, due to its softness, elasticity and moistening qualities (medical oil + vaseline). It works to cushion apical areas, as digital amputation cover, in-growing nail care, as keloid scar prevention. Size: Small Packaging requirements: Packet of 12 units Estimated usage: 1 | Yes, representative sample for items 143.1 to 143.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-143.2 | 137064 | CAP digital, polymer gel, soft, medium CAP digital textile fabric inside coated with polymer gel, <u>soft</u> , stretchable, protects toes against rubbing and pressure, transparent, durable, washable, due to its softness, elasticity and moistening qualities (medical oil + vaseline). It works to cushion apical areas, as digital amputation cover, in-growing nail care, as keloid scar prevention. Size: Medium Packaging requirements: Packet of 12 units Estimated usage: 1 | Yes, representative sample for items 143.1 to 143.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM-143.3 | 137065 | CAP digital, polymer gel, soft, large CAP digital textile fabric inside coated with polymer gel, <u>soft</u> , stretchable, protects toes against rubbing and pressure, transparent, durable, washable, due to its softness, elasticity and moistening qualities (medical oil + vaseline). It works to cushion apical areas, as digital amputation cover, in-growing nail care, as keloid scar prevention. Size: Large Packaging requirements: Packet of 12 units Estimated usage: 1 | Yes, representative sample for items 143.1 to 143.3 | | Per packet | | | | | | | | | | | | | | |
| ITEM 144 | ITEM 144 | Heel Spur, polymer gel, soft | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|---|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-144.1 | 170811 | Heel Spur, polymer gel, soft, small Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. Size: Small Packaging requirements: Per pair <u>Estimated usage:</u> 1 | Yes, representative sample for items 144.1 to 144.4 | | Per pair | | | | | | | | | | | | | | |
| ITEM-144.2 | 137058 | Heel Spur, polymer gel, soft, medium Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. Size: Medium Packaging requirements: Per pair <u>Estimated usage:</u> 1 | Yes, representative sample for items 144.1 to 144.4 | | Per pair | | | | | | | | | | | | | | |
| ITEM-144.3 | 137059 | Heel Spur, polymer gel, soft, large Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. Size: Large Packaging requirements: Per pair <u>Estimated usage:</u> 1 | Yes, representative sample for items 144.1 to 144.4 | | Per pair | | | | | | | | | | | | | | |
| ITEM-144.4 | 174308 | Heel Spur, Polymer gel, soft, X-Large Heel spur, made of polymer gel, <u>soft</u> , good stretch, transparent, washable, consists of a pair of shock absorbing heel cushions with removable plugs to fit specific comfort needs. Cushioning plugs must be <u>extra soft</u> , provides relief for general or specific heel pain. Size: X-Large Packaging requirements: Per pair <u>Estimated usage:</u> 1 | Yes, representative sample for items 144.1 to 144.4 | | Per pair | | | | | | | | | | | | | | |
| ITEM 145 | Rings, digital, textile, polymer gel | | | | | | | | | | | | | | | | | | |
| ITEM-145.1 | 170822 | Rings, digital, textile, polymer gel, small Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of soft polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. Length: ± 2.5cm Size: Small Packaging requirements: Pack of 10 rings <u>Estimated usage:</u> 1 | Yes, representative sample for items 145.1 to 145.4 | | Per packet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-145.2 | 170823 | <p>Rings, digital, textile, polymer gel, medium Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of <u>soft</u> polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. Length: ± 2.5cm Size: Medium Packaging requirements: Pack of 10 rings <u>Estimated usage:</u> 1</p> | Yes, representative sample for items 145.1 to 145.4 | | Per packet | | | | | | | | | | | | | | |
| ITEM-145.3 | 170825 | <p>Rings, digital, textile, polymer gel, large Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of <u>soft</u> polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. Length: ± 2.5cm Size: Large Packaging requirements: Pack of 10 rings <u>Estimated usage:</u> 1</p> | Yes, representative sample for items 145.1 to 145.4 | | Per packet | | | | | | | | | | | | | | |
| ITEM-145.4 | 170826 | <p>Rings, digital, textile, polymer gel, X-Large Rings, digital textile made of polymer gel, good stretch, washable discs / tubes of <u>soft</u> polymer gel attached to flexible textile fabric, provides specific protection while adjusting to the toes for comfort and durable use. Length: ± 2.5cm Size: X-Large Packaging requirements: Pack of 10 rings <u>Estimated usage:</u> 1</p> | Yes, representative sample for items 145.1 to 145.4 | | Per packet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 146 | 159806 | Silicone putty kit, soft mouldable for relief Silicone putty kit, mouldable for pressure-relief and correcting digital deformities, normal density, shore density of between 12-14, 250g base with 250g catalyst, mixed 1:1; it is soft, quick setting, self-mixing, pliable and washable. Size: Base with catalyst = 500g product Packaging requirements: Per kit Estimated usage: 1 | Yes | | Per kit | | | | | | | | | | | | | | |
| ITEM-146.1 | 159806 | Silicone putty kit, hard, mouldable for relief Silicone putty kit, mouldable for pressure-relief and correcting digital deformities, normal density, shore density of 20-25, 250g base with 250g catalyst, mixed 1:1; it is hard, quick setting, self-mixing, pliable and washable. Size: Base with catalyst = 500g product Packaging requirements: Per kit Estimated usage: 1 | Yes | | Per kit | | | | | | | | | | | | | | |
| ITEM 147 | 114127 | Silicone putty kit, hard, mouldable for relief Silicone putty kit, mouldable for pressure-relief and correcting digital deformities, normal density, shore density of 20-25, 250g base with 250g catalyst, mixed 1:1; it is hard, quick setting, self-mixing, pliable and washable. Size: Base with catalyst = 500g product Packaging requirements: Per kit Estimated usage: 1 | Yes | | Per kit | | | | | | | | | | | | | | |
| ITEM 148 | 159695 | Insole, shoe, thin moulded EVA Insole, shoe size 4, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 4 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-148.1 | 159697 | Insole, shoe size 5, thin moulded EVA Insole, shoe size 5, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 5 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-148.2 | 136782 | Insole, shoe size 6, thin moulded EVA Insole, shoe size 6, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 6 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-148.3 | 136783 | Insole, shoe size 7, thin moulded EVA Insole, shoe size 7, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 7 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------------------------------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-148.5 | 136784 | Insole, shoe size 8, thin moulded EVA Insole, shoe size 8, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 8 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-148.6 | 159700 | Insole, shoe size 9, thin moulded EVA Insole, shoe size 9, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 9 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-148.7 | 136786 | Insole, shoe size 10, thin moulded EVA Insole, shoe size 10, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 10 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-148.8 | 159702 | Insole, shoe size 11, thin moulded EVA Insole, shoe size 11, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 11 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-148.9 | 159703 | Insole, shoe size 12, thin moulded EVA Insole, shoe size 12, thin moulded EVA (Ethyl Vinyl Acetate) innersoles with a strengthened deep heel cup covered with quick dry moisture absorbing material. Size: Shoe size 12 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 148.1 to 148.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM 149 | Insole, shoe, thin cork | | | | | | | | | | | | | | | | | | |
| ITEM-149.1 | 136787 | Insole, shoe size 4, thin cork Insole, shoe size 4, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 4 Packaging requirements: Per pair Estimated usage: New item | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-149.2 | 136788 | Insole, shoe size 5, thin cork Insole, shoe size 5, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 5 Packaging requirements: Per pair Estimated usage: New item | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-149.3 | 136789 | Insole, shoe size 6, thin cork Insole, shoe size 6, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 6 Packaging requirements: Per pair Estimated usage: New item | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-149.4 | 136790 | Insole, shoe size 7, thin cork Insole, shoe size 7, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 7 Packaging requirements: Per pair Estimated usage: New item | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-149.5 | 136791 | Insole, shoe size 8, thin cork Insole, shoe size 8, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 8 Packaging requirements: Per pair Estimated usage: New item | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-149.6 | 136792 | Insole, shoe size 9, thin cork Insole, shoe size 9, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 9 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-149.7 | 136793 | Insole, shoe size 10, thin cork Insole, shoe size 10, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 10 Packaging requirements: Per pair Estimated usage: New item | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-149.8 | 136794 | Insole, shoe size 11, thin cork Insole, shoe size 11, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 11 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| ITEM-149.9 | 136795 | Insole, shoe size 12, thin cork Insole, shoe size 12, thin cork, pre-cut, 1-2mm thick. Size: Shoe size 12 Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 149.1 to 149.9 | | Per pair | | | | | | | | | | | | | | |
| item 150 | | Innersole, polyurethane | | | | | | | | | | | | | | | | | |
| ITEM-150.1 | 114128 | Innersole, polyurethane, shoe size 4, full length Innersole, polyurethane foam, shoe size 4, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 4 full length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.1 to 150.7 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.2 | 114130 | Innersole, polyurethane, shoe size 5, full length Innersole, polyurethane foam, shoe size 5, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 5 full length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.1 to 150.7 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-150.3 | 159696 | Innersole, polyurethane, shoe size 6, full length Innersole, polyurethane foam, shoe size 6, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 6 full length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.1 to 150.7 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.4 | 159698 | Innersole, polyurethane, shoe size 7, full length Innersole, polyurethane foam, shoe size 7, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 7 full length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.1 to 150.7 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.5 | 159699 | Innersole, polyurethane, shoe size 8, full length Innersole, polyurethane foam, shoe size 8, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 8 full length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.1 to 150.7 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.6 | 114135 | Innersole, polyurethane, shoe size 9, full length Innersole, polyurethane foam, shoe size 9, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 9 full length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.1 to 150.7 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.7 | 159701 | Innersole, polyurethane, shoe size 10, full length Innersole, polyurethane foam, shoe size 10, full length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is black micro fibre material. Size: Shoe size 10 full length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.1 to 150.7 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-150.8 | 114129 | Innersole, polyurethane, shoe size 4, ¾ length Innersole, polyurethane foam, shoe size 4, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 4 ¾ length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.8 to 150.12 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.9 | 114131 | Innersole, polyurethane, shoe size 5, ¾ length Innersole, polyurethane foam, shoe size 5, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 5 ¾ length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.8 to 150.12 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|-------------|--------|---|--|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-150.10 | 114132 | Innersole, polyurethane, shoe size 6, ¾ length Innersole, polyurethane foam, shoe size 6, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 6 ¾ length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.8 to 150.12 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.11 | 114133 | Innersole, polyurethane, shoe size 7, ¾ length Innersole, polyurethane foam, shoe size 7, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 7 ¾ length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.8 to 150.12 | | Per pair | | | | | | | | | | | | | | |
| ITEM-150.12 | 114134 | Innersole, polyurethane, shoe size 8, ¾ length Innersole, polyurethane foam, shoe size 8, ¾ length, with a density of 15-25° shore A, hardness of 0.27 grams/cm ³ , must be breathable and absorbs sweat, polyurethane foam is well known for comfort, its elasticity and flexibility, covering material is micro fibre material. Size: Shoe size 8 ¾ length Packaging requirements: Per pair Estimated usage: 1 | Yes, representative sample for items 150.8 to 150.12 | | Per pair | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 151 | | Orthotic kit, medium, adhesive | | | | | | | | | | | | | | | | | |
| ITEM-151.1 | New | Orthotic kit, medium, adhesive, shoe size 5, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 5 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 151.1 to 151.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-151.2 | New | Orthotic kit, medium, adhesive, shoe size 6, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 6 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 151.1 to 151.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-151.3 | New | Orthotic kit, medium, adhesive, shoe size 7, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 7 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 151.1 to 151.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-151.4 | New | Orthotic kit, medium, adhesive, shoe size 8, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 8 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 151.1 to 151.7 | | Per kit | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-151.5 | New | Orthotic kit, medium, adhesive, shoe size 9, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 9 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 151.1 to 151.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-151.6 | New | Orthotic kit, medium, adhesive, shoe size 10, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 10 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 151.1 to 151.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-151.7 | New | Orthotic kit, medium, adhesive, shoe size 11, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 11 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 151.1 to 151.7 | | Per kit | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|------------------------------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 152 | Orthotic kit, hard, adhesive | | | | | | | | | | | | | | | | | | |
| ITEM-152.1 | New | Orthotic kit, hard, adhesive, shoe size 5, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 5 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 152.1 to 152.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-152.2 | New | Orthotic kit, hard, adhesive, shoe size 6, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 6 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 152.1 to 152.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-152.3 | New | Orthotic kit, hard, adhesive, shoe size 7, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 7 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 152.1 to 152.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-152.4 | New | Orthotic kit, hard, adhesive, shoe size 8, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 8 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 152.1 to 152.7 | | Per kit | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|---|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-152.5 | New | Orthotic kit, hard, adhesive, shoe size 9, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 9 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 152.1 to 152.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-152.6 | New | Orthotic kit, hard, adhesive, shoe size 10, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 10 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 152.1 to 152.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-152.7 | New | Orthotic kit, hard, adhesive, shoe size 11, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex adhesive wedge for left and right foot (pair); latex adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 11 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 152.1 to 152.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM 153 | Orthotic kit, medium, non-adhesive | | | | | | | | | | | | | | | | | | |
| ITEM-153.1 | New | Orthotic kit, medium, non-adhesive, shoe size 5, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 5 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 153.1 to 153.7 | | Per kit | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | % DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|-----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-153.2 | New | Orthotic kit, medium, non-adhesive, shoe size 6, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 6 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 153.1 to 153.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-153.3 | New | Orthotic kit, medium, non-adhesive, shoe size 7, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 7 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 153.1 to 153.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-153.4 | New | Orthotic kit, medium, non-adhesive, shoe size 8, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 8 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 153.1 to 153.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-153.5 | New | Orthotic kit, medium, non-adhesive, shoe size 9, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 9 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 153.1 to 153.7 | | Per kit | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-153.6 | New | Orthotic kit, medium, non-adhesive, shoe size, 10 full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 10 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 153.1 to 153.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-153.7 | New | Orthotic kit, medium, non-adhesive, shoe size 11, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: medium . Size: Shoe size 11 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 153.1 to 153.7 | | Per kit | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-154 | | Orthotic kit, hard, non-adhesive | | | | | | | | | | | | | | | | | |
| ITEM-154.1 | New | Orthotic kit, hard, non-adhesive, shoe size, 5 full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 5 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 154.1 to 154.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-154.2 | New | Orthotic kit, hard, non-adhesive, shoe size, 6 full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 6 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 154.1 to 154.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-154.3 | New | Orthotic kit, hard, non-adhesive, shoe size 7, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 7 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 154.1 to 154.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-154.4 | New | Orthotic kit, hard, non-adhesive, shoe size 8, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 8 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 154.1 to 154.7 | | Per kit | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|------------------------------|---|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-154.5 | New | Orthotic kit, hard, non-adhesive, shoe size 9, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 9 full length Packaging requirements: Per kit | Yes, representative sample for items 154.1 to 154.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-154.6 | New | Orthotic kit, hard, non-adhesive, shoe size 10, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 10 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 154.1 to 154.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM-154.7 | New | Orthotic kit, hard, non-adhesive, shoe size 11, full length Kit must contain the following components: moulded resin neutral orthotic base with breathable, conformable cover for left and right foot (pair); latex non-adhesive metatarsal pad with the same lift for all 5 metatarsal heads for left and right foot (pair); latex non-adhesive wedge for left and right foot (pair); latex non-adhesive cobra/horse shoe pad with retractable core for left and right foot (pair). Degree of latex hardness: hard . Size: Shoe size 11 full length Packaging requirements: Per kit Estimated usage: New | Yes, representative sample for items 154.1 to 154.7 | | Per kit | | | | | | | | | | | | | | |
| ITEM 155 | Burr diamond crystals | | | | | | | | | | | | | | | | | | |
| ITEM-155.1 | 159098 | Burr diamond crystals barrel head shape Burr, diamond crystals bonded directly to one piece, barrel head shape, hardened stainless steel shank Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-155.2 | 159095 | Burr diamond crystals bud head shape Burr, diamond crystals bonded directly to one piece, bud head shape, hardened stainless steel shank Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|-----------------------|---|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-155.3 | 159099 | Burr diamond crystals cylinder-ball head Burr, diamond crystals bonded directly to one piece, cylinder-ball head shape, hardened stainless steel shank Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-155.4 | 159096 | Burr diamond crystals enucleator head Burr, diamond crystals bonded directly to one piece, enucleator head shape, hardened stainless steel shank Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-155.5 | 159093 | Burr diamond crystals olive head shape Burr, diamond crystals bonded directly to one piece, olive head shape, hardened stainless steel shank Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-155.6 | 159094 | Burr diamond crystals pear head shape Burr, diamond crystals bonded directly to one piece, pear head shape, hardened stainless steel shank Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-155.7 | 136797 | Burr tungsten carbide cone shape Burr, tungsten carbide, cone shape with extremely fine grain size, calibrated to high isostatic pressing process, blades are accurately sharpened ensuring optimal cutting, carbide is very hard with high longevity, better sterility, hardened stainless steel shank Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 156 | Disc garnet | | | | | | | | | | | | | | | | | | |
| ITEM-156.1 | 159136 | Disc garnet 2.25cm diameter course grid Disc garnet, disposable emery disc for finishing skin and calluses with square brass centre that fits onto a mandrill stem, course grid, 2.25cm diameter, 50 units per pack Dimensions: 2.25cm diameter, Packaging requirements: 50 units per pack Estimated usage: 1 | Yes | | Per packet | | | | | | | | | | | | | | |
| ITEM-156.2 | 159135 | Disc garnet 2.25cm diameter medium grid Disc garnet, disposable emery disc for finishing skin and calluses with square brass centre that fits onto a mandrill stem, medium grid, 2.25cm diameter, 50 units per pack Dimensions: 2.25cm diameter, Packaging requirements: 50 units per pack Estimated usage: 1 | Yes | | Per packet | | | | | | | | | | | | | | |
| ITEM 157 | Mandrill stem | | | | | | | | | | | | | | | | | | |
| ITEM-157.1 | 159133 | Mandrill stem standard dental size Mandrill stem, standard dental size, stainless steel shank to hold the garnet disc Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 158 | Silver Nitrate | | | | | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATEIRAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND %PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|-----------------------------|----------------------------------|----------------|-------------|--------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM-158.1 | 162969 | Silver Nitrate stick Silver nitrate sticks, 40% strength Size: per stick Packaging requirements: Pack of 10 sticks Estimated usage: 1 | Yes | | Per packet | | | | | | | | | | | | | | |
| ITEM-158.2 | 137022 | Silver Nitrate liquid 40% strength Size: 10ml (dropper bottle) Packaging requirements: Per bottle Estimated usage: 1 | Yes | | Per bottle | | | | | | | | | | | | | | |
| ITEM-158.3 | 136790 | Silver Nitrate granules Silver nitrate granules, clear and odourless crystals, very soluble in water, 40% strength Size: 5g granules in 50ml bottle Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|---|-----------------|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 159 | 178163 | Oil spray, anti rust Anti rust oil spray for stainless steel instruments (pre-sterilization) to increase life span of instruments Size: 300ml container Packaging requirements: Each Estimated usage: 1 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM-159.1 | | | | | | | | | | | | | | | | | | | |
| ITEM 160 | 178163 | Blade, scalpel, surgical Blade, scalpel, surgical, no. 10 Blade, scalpel, surgical, to fit handle no's 3, 5 and 7, with small fitting features, no. 10, must be compatible with scalpel handles for items (short no. 3, long no. 5 and long no.7), composition test for carbon steel and stainless steel, ISO900101 and ISO13485. Must be sterile, peel apart pack must peel evenly without tearing paper lining, blade size must be printed on the peel apart pack. Size: No. 10 blade Packaging requirements: Packet of 100 Estimated usage: 100 packets | Yes | | Per packet | | | | | | | | | | | | | | |
| ITEM-160.1 | New | | | | | | | | | | | | | | | | | | |
| ITEM-160.2 | New | Blade, scalpel, surgical, no. 15 Blade, scalpel, surgical, to fit handle no's 3, 5 and 7, with small fitting features, no. 15, must be compatible with scalpel handles for items (short no. 3, long no. 5 and long no.7), composition test for carbon steel and stainless steel, ISO900101 and ISO13485. Must be sterile, peel apart pack must peel evenly without tearing paper lining, blade size must be printed on the peel apart pack. Size: No. 15 blade Packaging requirements: Packet of 100 Estimated usage: 100 packets | Yes | | Per packet | | | | | | | | | | | | | | |

| ITEM NO | SAP NO | ITEM DESCRIPTION | SAMPLE REQUIRED | COMPLY YES/NO | UNIT OF MEASURE | TENDER PRICE/UNIT | STATE PACKAGING/BOX QUANTITY | BRAND | LEAD TIME (DAYS) | DIRECT MATERAL: % IMPORTED | DIRECT MATERIAL: % LOCAL CONTENT | %DIRECT LABOUR | % TRANSPORT | % OTHER OVERHEADS AND % PROFITS | %TOTAL | RATE OF EXCHANGE | ADJUSTMENT FREQUENCY | IMPORTED | LOCALLY MANUFACTURED |
|------------|--------|--|---|---------------|-----------------|-------------------|------------------------------|-------|------------------|----------------------------|----------------------------------|----------------|-------------|---------------------------------|--------|------------------|----------------------|----------|----------------------|
| ITEM 161 | | Monofilament | | | | | | | | | | | | | | | | | |
| ITEM-161.1 | 124117 | Monofilament, retractable, 10g Monofilament, nylon, retractable, 10g. Size: 10g Packaging requirements: Each Estimated usage: 100 | Yes | | Each | | | | | | | | | | | | | | |
| ITEM 162 | | Tourniquet | | | | | | | | | | | | | | | | | |
| ITEM-162.1 | New | Tourniquet, digital, small Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: orange Size: Small, 0.6 - 0.9cm diameter Packaging requirements: Box of 20 Estimated usage: 0 | Yes, representative sample for items 162.1 to 162.4 | | Per box | | | | | | | | | | | | | | |
| ITEM-162.2 | New | Tourniquet, digital, medium Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: red Size: Medium, ±1cm diameter Packaging requirements: Box of 20 Estimated usage: 1 | Yes, representative sample for items 162.1 to 162.4 | | Per box | | | | | | | | | | | | | | |
| ITEM-162.3 | New | Tourniquet, digital, large Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: green Size: Large, ±1.25cm diameter Packaging requirements: Box of 20 Estimated usage: 1 | Yes, representative sample for items 162.1 to 162.4 | | Per box | | | | | | | | | | | | | | |
| ITEM-162.4 | New | Tourniquet, digital, X-Large Tourniquet, digital, latex free, circular incorporating a sever gap for ease of removal. For podiatry use to achieve safe and even pressure to stop blood flow in digits. One per sterile peel pouch for single use. Colour: blue Size: X-Large, ±1.5cm diameter Packaging requirements: Box of 20 Estimated usage: 1 | Yes, representative sample for items 162.1 to 162.4 | | Per box | | | | | | | | | | | | | | |



GAUTENG PROVINCE
PROVINCIAL TREASURY
REPUBLIC OF SOUTH AFRICA

Provincial Supply Chain Management

Sample Sheet

Page 1 of 2

SAMPLES FOR RFP

PLEASE TICK WHERE SAMPLE IS SUBMITTED

NAME OF CONTACT PERSON

TEI NUMBER



GAUTENG PROVINCE

PROVINCIAL TREASURY
REPUBLIC OF SOUTH AFRICA

Provincial Supply Chain Management

Sample Sheet

Page 2 of 2

PLEASE TICK WHERE SAMPLE IS SUBMITTED



GAUTENG PROVINCE
PROVINCIAL TREASURY
REPUBLIC OF SOUTH AFRICA

Provincial Supply Chain Management

Sample Sheet

Page 1 of 2

SAMPLES FOR RFP

PLEASE TICK WHERE SAMPLE IS SUBMITTED

NAME OF CONTACT PERSON

TEI NUMBER



GAUTENG PROVINCE

PROVINCIAL TREASURY
REPUBLIC OF SOUTH AFRICA

Provincial Supply Chain Management

Sample Sheet

Page 2 of 2

PLEASE TICK WHERE SAMPLE IS SUBMITTED



GAUTENG PROVINCE

PROVINCIAL TREASURY
REPUBLIC OF SOUTH AFRICA

Provincial Supply Chain Management

Sample Sheet

Page 1 of 2

SAMPLES FOR RFP

PLEASE TICK WHERE SAMPLE IS SUBMITTED

NAME OF CONTACT PERSON

TEI NUMBER



GAUTENG PROVINCE

PROVINCIAL TREASURY
REPUBLIC OF SOUTH AFRICA

Provincial Supply Chain Management

Sample Sheet

Page 2 of 2

PLEASE TICK WHERE SAMPLE IS SUBMITTED



Provincial Supply Chain Management

Non-firm Prices

Page 1 of 5

NOTE: If the information required regarding each item cannot be inserted in the space provided on either page 1, 2, 3, 4 or 5 of this form, additional information may be provided on a photocopy of this form.

The *Gauteng Provincial Treasury (GPT) is acting as Common Service Provider or buying organisation on behalf of all Gauteng Provincial Government Customer Departments / Institutions. The goods / services are therefore required by the Customer Department / Institution, as indicated on RFP 01.

"Firm" prices are deemed to be fixed prices, which are only subject to the following statutory changes, namely VAT and any levy related to customs and excise. "Non-firm" prices are all prices other than "firm" prices and which are either subject to Rate of Exchange claims and/or subject to published indices or documentary evidence.

In case of a contract being awarded as a result of this RFP, all applications for price adjustments must be referred to the GPT, Private Bag X112, Marshalltown. 2107. Price adjustments shall not be backdated and shall be applicable from date of approval by the GPT. The successful supplier shall submit an application, in writing and supported by documentary proof, to the GPT thirty (30) days before the effective date of any price adjustment, whether based on RoE (par. 1), formula(e) (par. 1 and/or 2) or documentary proof (par. 3).

NOTE: Non-firm prices may be loaded when calculating comparative prices.

1. CLAIMS SUBJECT TO RATE OF EXCHANGE (RoE)

In case of bid prices being subject to Rate of Exchange (RoE) variations, please specify the amount, as well as the applicable rate of exchange at the time of bidding applicable to each item.

Price adjustments, based on changes in the rate of exchange, shall only be considered if the change, based on the previous rate, had moved by more than 5%.

1.1 COMPILING THE RoE FORMULA - Cost breakdown structure

A cost breakdown in terms of local & imported components must be submitted for each item as follows:

| Imported portion | % | Local portion | % | Total | 100% |
|------------------|---|---------------|---|-------|------|
|------------------|---|---------------|---|-------|------|

Whatever breakdown chosen in terms of local and imported components cannot be changed over the duration of any contract awarded based on this RFP.

1.2 AMOUNT(S) SUBJECT TO RoE VARIANCE - IMPORTED PORTION OF RoE FORMULA

Amount(s) subject to exchange rate variations – please specify for each item number individually.

| | |
|--------------|--|
| Item N° | |
| ROE Variance | |

(as per item number on RFP 06 Price Schedule)

1.3 BASE RoE - IMPORTED PORTION OF RoE FORMULA

Exchange rate applicable as at the time of bidding – please specify every currency applicable and indicate which items are affected by which exchange rate

| | |
|---------------|--|
| Item N° | |
| EXCHANGE RATE | |

(as per item number on RFP 06 Price Schedule)



Provincial Supply Chain Management

Non-firm Prices

Page 2 of 5

2. CLAIMS SUBJECT TO PUBLISHED INDICES

For all items to which non-firm prices apply, bidders must submit a cost breakdown, which indicates the proportion of the prices subject to adjustment.

| Item № |
|--------|
|--------|

(as per item number on RFP 06 Price Schedule)

Where price adjustments based on an escalation formula is chosen as the method of price adjustment, variations in the actual costs for whatever reason are for the account of the successful bidder. This implies that:

- formula-based adjustments (as per this method in par 2) and
- proven cost increases (as per the method in par 3) cannot both be entertained at the same time for the same cost element.

2.1 METHOD 1: FORMULA BASED ADJUSTMENTS

2.1.1 The escalation (or price adjustment) formula shall provide for a minimum fixed element of at least 15% of the price, which shall not be subject to adjustment. The following standard formula structure is prescribed for calculating escalation / price adjustments:

IN CASES OF PERIOD CONTRACTS, NON FIRM PRICES WILL BE ADJUSTED (LOADED) WITH THE ASSESSED CONTRACT PRICE ADJUSTMENTS IMPLICIT IN NON FIRM PRICES WHEN CALCULATING THE COMPARATIVE PRICES

IN THIS CATEGORY PRICE ESCALATIONS WILL ONLY BE CONSIDERED IN TERMS OF THE FOLLOWING FORMULA:

$$Pa = (1 - V)Pt \left(D1 \frac{R1t}{R1o} + D2 \frac{R2t}{R2o} + D3 \frac{R3t}{R3o} + D4 \frac{R4t}{R4o} \right) + VPt$$

Where:

| | | |
|---------------|---|--|
| Pa | = | The new escalated price to be calculated. |
| (1-V)Pt | = | 85% of the original bid price. Note that Pt must always be the original bid price and not an escalated price. |
| D1, D2.. | = | Each factor of the bid price eg. labour, transport, clothing, footwear, etc. The total of the various factors D1, D2...etc. must add up to 100%. |
| R1t, R2t..... | = | Index figure obtained from new index (depends on the number of factors used). |
| R1o, R2o | = | Index figure at time of bidding. |
| VPt | = | 15% of the original bid price. This portion of the bid price remains firm i.e. it is not subject to any price escalations. |

The following index/indices must be used to calculate your bid price:

| Index | Dated | Index | Dated | Index | Dated |
|-------|-------|-------|-------|-------|-------|
| | | | | | |
| | | | | | |



Provincial Supply Chain Management

Non-firm Prices

Page 3 of 5

2.1.2 Cost breakdown structure

A percentage (%) cost breakdown for each item must be submitted according to the following example – actual figures to be added to the table (2nd column) in par. 2.1.3 below:

– the various elements specified in this paragraph will form the series (D1, D2 ... up to Dn) in the formula to be indicated in par 2.1.1 above.

EXAMPLE OF COST BREAKDOWN STRUCTURE

| | |
|-----------------|-------|
| Direct labour | 25%** |
| Direct material | 30%** |
| Transport | 25%** |
| Other overheads | 20%** |
| <i>Total</i> | 100% |

FURNISH A BREAKDOWN OF YOUR PRICE IN TERMS OF ABOVE-MENTIONED FORMULA. THE TOTAL OF THE VARIOUS FACTORS MUST ADD UP TO 100%.

| FACTOR (D1, D2 etc. eg. Labour, transport etc.) | PERCENTAGE OF BID PRICE |
|--|-------------------------|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

Whatever breakdown of index series are chosen cannot be changed over the duration of any contract awarded based on this RFP.

2.2 APPLICATION OF THE PRICE ADJUSTMENT FORMULA

2.2.1. Fixed price period

GPG suggests an initial fixed period of at least six (6) months from the effective date of any agreement, which may be awarded as a result of this RFP. Specify from what date the first price adjustment would be applicable (as per item number on RFP 06 Price Schedule)

2.2.2 Frequency of price adjustments after fixed price period

GPG suggests quarterly adjustments, after the initial fixed price period. Longer periods than three (3) months (quarterly) between adjustments will be considered even more favourably.

(as per item number on RFP 06 Price Schedule)



Provincial Supply Chain Management

Non-firm Prices

Page 4 of 5

3. METHOD 2: ADJUSTMENTS BASED ON DOCUMENTARY PROOF

Price adjustments shall also be considered based on independent trust worthy documentary proof of cost, price variances or payments (forward cover cost).

PGP would however prefer claims based on published indices (as per the method in par 2 above), as a result of the fact that the confidentiality regarding source and price of goods can be maintained with a formula. Providing substantiating or supporting documentation (as per this method in par 3) often reveals information that may be considered confidential by the bidder.

Please again note that

- formula-based adjustments (as per the method in par 2) and
- proven cost increases (as per this method in par 3) cannot both be entertained at the same time for the same cost element.

3.1. Amount(s) subject to change based on documentary proof –

Please specify for each item individually.

| | |
|----------------|--|
| Item N° | |
| Amount | |

(as per item number on RFP 06 Price Schedule)

3.2.1. Fixed price period

PGP suggests an initial fixed period of at least six (6) months from the effective date of any agreement, which may be awarded as a result of this RFP. Specify from what date the first price adjustment would be applicable.

(as per item number on RFP 06 Price Schedule)

3.2.2 Frequency of price adjustments after fixed price period

PGP suggests quarterly adjustments, after the initial fixed price period. Longer periods than three (3) months (quarterly) between adjustments will be considered even more favourably.

(as per item number on RFP 06 Price Schedule)



Provincial Supply Chain Management

Non-firm Prices

Page 5 of 5

4 PRICES SUBJECT TO RATE OF EXCHANGE VARIATIONS

4.1. Please furnish full particulars of your financial institution, state the currencies used in the conversion of the prices of the items to South African currency, which portion of the price is subject to rate of exchange variations and the amounts remitted abroad.

| PARTICULARS OF FINANCIAL INSTITUTION | ITEM NO | PRICE | CURRENCY | RATE | PORTION OF PRICE SUBJECT TO ROE | AMOUNT IN FOREIGN CURRENCY REMITTED ABROAD |
|--------------------------------------|---------|-------|----------|------|---------------------------------|--|
| | | | | ZAR= | | |
| | | | | ZAR= | | |
| | | | | ZAR= | | |
| | | | | ZAR= | | |
| | | | | ZAR= | | |
| | | | | ZAR= | | |

4.2 Adjustments for rate of exchange variations during the contract period will be calculated by using the average monthly exchange rates as issued by your commercial bank for the periods indicated hereunder: (Proof from bank required)

| AVERAGE MONTHLY EXCHANGE RATES FOR THE PERIOD: | DATE DOCUMENTATION MUST BE SUBMITTED TO THIS OFFICE | DATE FROM WHICH NEW CALCULATED PRICES WILL BECOME EFFECTIVE | DATE UNTIL WHICH NEW CALCULATED PRICE WILL BE EFFECTIVE |
|--|---|---|---|
| | | | |
| | | | |
| | | | |

Name(s) & Signature(s) of Bidder(s)

CAPACITY

DATE:



Provincial Supply Chain Management

PREFERENCE POINTS CLAIM FORM

Page: 1 of 7

PREFERENCE POINTS CLAIM FORM IN TERMS OF THE PREFERENTIAL PROCUREMENT REGULATIONS 2022

This preference form must form part of all tenders invited. It contains general information and serves as a claim form for preference points for specific goals.

NB: BEFORE COMPLETING THIS FORM, TENDERERS MUST STUDY THE GENERAL CONDITIONS, DEFINITIONS AND DIRECTIVES APPLICABLE IN RESPECT OF THE TENDER AND PREFERENTIAL PROCUREMENT REGULATIONS, 2022

1. GENERAL CONDITIONS

1.1 The following preference point systems are applicable to invitations to tender:

- the 80/20 system for requirements with a Rand value of up to R50 000 000 (all applicable taxes included); and
- the 90/10 system for requirements with a Rand value above R50 000 000 (all applicable taxes included).

1.2 **To be completed by the organ of state**

[TICK APPLICABLE BOX]

| | |
|--|--|
| | The applicable preference point system for this tender is the 90/10 preference point system. |
| | The applicable preference point system for this tender is the 80/20 preference point system. |
| | Either the 90/10 or 80/20 preference point system will be applicable in this tender. The lowest/ highest acceptable tender will be used to determine the accurate system once tenders are received. |

1.3 Points for this tender (even in the case of a tender for income-generating contracts) shall be awarded for:

- (a) Price; and
- (b) Specific Goals.



Provincial Supply Chain Management

PREFERENCE POINTS CLAIM FORM

Page: 2 of 7

1.4 To be completed by the organ of state:

The maximum points for this tender are allocated as follows:

| | POINTS |
|--|---------------|
| PRICE | |
| SPECIFIC GOALS | |
| Total points for Price and SPECIFIC GOALS | 100 |

1.5 Failure on the part of a tenderer to submit proof or documentation required in terms of this tender to claim points for specific goals with the tender, will be interpreted to mean that preference points for specific goals are not claimed.

1.6 The organ of state reserves the right to require of a tenderer, either before a tender is adjudicated or at any time subsequently, to substantiate any claim in regard to preferences, in any manner required by the organ of state.

2. DEFINITIONS

- (a) “**tender**” means a written offer in the form determined by an organ of state in response to an invitation to provide goods or services through price quotations, competitive tendering process or any other method envisaged in legislation;
- (b) “**price**” means an amount of money tendered for goods or services, and includes all applicable taxes less all unconditional discounts;
- (c) “**rand value**” means the total estimated value of a contract in Rand, calculated at the time of bid invitation, and includes all applicable taxes;
- (d) “**tender for income-generating contracts**” means a written offer in the form determined by an organ of state in response to an invitation for the origination of income-generating contracts through any method envisaged in legislation that will result in a legal agreement between the organ of state and a third party that produces revenue for the organ of state, and includes, but is not limited to, leasing and disposal of assets and concession contracts, excluding direct sales and disposal of assets through public auctions; and
- (e) “**the Act**” means the Preferential Procurement Policy Framework Act, 2000 (Act No. 5 of 2000).



Provincial Supply Chain Management

PREFERENCE POINTS CLAIM FORM

Page: 3 of 7

3. FORMULAE FOR PROCUREMENT OF GOODS AND SERVICES

3.1. POINTS AWARDED FOR PRICE

3.1.1 THE 80/20 OR 90/10 PREFERENCE POINT SYSTEMS

A maximum of 80 or 90 points is allocated for price on the following basis:

80/20

or

90/10

$$Ps = 80 \left(1 - \frac{Pt - P_{min}}{P_{min}} \right) \text{ or } Ps = 90 \left(1 - \frac{Pt - P_{min}}{P_{min}} \right)$$

Where

Ps = Points scored for price of tender under consideration

Pt = Price of tender under consideration

Pmin = Price of lowest acceptable tender

3.2. FORMULAE FOR DISPOSAL OR LEASING OF STATE ASSETS AND INCOME GENERATING PROCUREMENT

3.2.1. POINTS AWARDED FOR PRICE

A maximum of 80 or 90 points is allocated for price on the following basis:

80/20

or

90/10

$$Ps = 80 \left(1 + \frac{Pt - P_{max}}{P_{max}} \right) \text{ or } Ps = 90 \left(1 + \frac{Pt - P_{max}}{P_{max}} \right)$$

Where

Ps = Points scored for price of tender under consideration

Pt = Price of tender under consideration

Pmax = Price of highest acceptable tender



Provincial Supply Chain Management

PREFERENCE POINTS CLAIM FORM

Page: 4 of 7

4. POINTS AWARDED FOR SPECIFIC GOALS

- 4.1. In terms of Regulation 4(2); 5(2); 6(2) and 7(2) of the Preferential Procurement Regulations, preference points must be awarded for specific goals stated in the tender. For the purposes of this tender the tenderer will be allocated points based on the goals stated in table 1 below as may be supported by proof/ documentation stated in the conditions of this tender:
- 4.2. In cases where organs of state intend to use Regulation 3(2) of the Regulations, which states that, if it is unclear whether the 80/20 or 90/10 preference point system applies, an organ of state must, in the tender documents, stipulate in the case of—
 - (a) an invitation for tender for income-generating contracts, that either the 80/20 or 90/10 preference point system will apply and that the highest acceptable tender will be used to determine the applicable preference point system; or
 - (b) any other invitation for tender, that either the 80/20 or 90/10 preference point system will apply and that the lowest acceptable tender will be used to determine the applicable preference point system,

then the organ of state must indicate the points allocated for specific goals for both the 90/10 and 80/20 preference point system.



Provincial Supply Chain Management

PREFERENCE POINTS CLAIM FORM

Page: 5 of 7

Table 1: Specific goals for the tender and points claimed are indicated per the table below.

(Note to organs of state: Where either the 90/10 or 80/20 preference point system is applicable, corresponding points must also be indicated as such.

Note to tenderers: The tenderer must indicate how they claim points for each preference point system.)

| The specific goals allocated points in terms of this tender | Number of points allocated (90/10 system) (To be completed by the organ of state) | Number of points allocated (80/20 system) (To be completed by the organ of state) | Number of points claimed (90/10 system) (To be completed by the tenderer) | Number of points claimed (80/20 system) (To be completed by the tenderer) |
|---|--|--|--|--|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |



Provincial Supply Chain Management

PREFERENCE POINTS CLAIM FORM

Page: 6 of 7

DECLARATION WITH REGARD TO COMPANY/FIRM

4.3. Name of company/firm

4.4. Company registration number

4.5. TYPE OF COMPANY/ FIRM

[TICK APPLICABLE BOX]

| | |
|--|--|
| | Partnership/Joint Venture / Consortium |
| | One-person business/sole proprietor |
| | Close corporation |
| | Public Company |
| | Personal Liability Company |
| | (Pty) Limited |
| | Non-Profit Company |
| | State Owned Company |

4.6. I, the undersigned, who is duly authorised to do so on behalf of the company/firm, certify that the points claimed, based on the specific goals as advised in the tender, qualifies the company/ firm for the preference(s) shown and I acknowledge that:

- i) The information furnished is true and correct;
- ii) The preference points claimed are in accordance with the General Conditions as indicated in paragraph 1 of this form;
- iii) In the event of a contract being awarded as a result of points claimed as shown in paragraphs 1.4 and 4.2, the contractor may be required to furnish documentary proof to the satisfaction of the organ of state that the claims are correct;
- iv) If the specific goals have been claimed or obtained on a fraudulent basis or any of the conditions of contract have not been fulfilled, the organ of state may, in addition to any other remedy it may have –



Provincial Supply Chain Management

PREFERENCE POINTS CLAIM FORM

Page: 7 of 7

- (a) disqualify the person from the tendering process;
- (b) recover costs, losses or damages it has incurred or suffered as a result of that person's conduct;
- (c) cancel the contract and claim any damages which it has suffered as a result of having to make less favourable arrangements due to such cancellation;
- (d) recommend that the tenderer or contractor, its shareholders and directors, or only the shareholders and directors who acted on a fraudulent basis, be restricted from obtaining business from any organ of state for a period not exceeding 10 years, after the *audi alteram partem* (hear the other side) rule has been applied; and
- (e) forward the matter for criminal prosecution, if deemed necessary.

| | |
|------------------------------------|--|
| SIGNATURE(S) OF TENDERER(S) | |
| SURNAME AND NAME | |
| DATE | |
| ADDRESS | |