

REQUEST FOR QUOTATION FORM FOR GOODS AND SERVICES

AGRICULTURAL RESEARCH COUNCIL(ARC)

1. You are kindly requested to submit a written quotation to Agricultural Research Council as per below or attached specification.

Description	Quantity	UOM
- LB, Quick-Preparation Agar Media (with Amp) Pre-mixed, pre-sterilized E. coli growth agar medium that can be prepared in just five minutes without autoclaving, containing everything you need to prepare low-salt LB agar plates, with ampicillin - LB, Quick-Preparation Liquid Media (with Amp) Pre-mixed, pre-sterilized E. coli growth liquid medium that can be prepared in just five minutes without autoclaving, containing everything you need to prepare low-salt LB liquid medium with ampicillin - LB, Quick-Preparation Agar Blue screening Media (with Amp) Pre-mixed, pre-sterilized E. coli growth agar medium with X-gal/IPTG for blue/white colony screening that can be prepared in just five minutes without autoclaving, containing everything you need to prepare blue/white screening agar plates with ampicillin	05 EACH	REQ-054170
10 ul Filter tips X 30 10 x 96/PK DNAase/RNAase free Filtered tips that are able to fit on Finn pipettes and gilson pipettes without any leakage 200ul Filter tips X05 10 x 96/PK DNAase/RNAase free Filtered tips that are able to fit on Finn pipettes and gilson pipettes without any leakage 1000ul Filter tips X10 10 x 96/PK DNAase/RNAase free Filtered tips that are able to fit on Finn pipettes and gilson pipettes without any leakage		REQ-059629

10 ul Filter tips X30 10 x 96/PK DNAase/RNAase free Filtered tips that are able to fit on Finn pipettes and gilson pipettes without any leakage 200ul Filter tips X15 10 x 96/PK DNAase/RNAase free Filtered tips that are able to fit on Finn pipettes and gilson pipettes without any leakage 1000ul Filter tips X20 10 x 96/PK DNAase/RNAase free Filtered tips that are able to fit on Finn pipettes and gilson pipettes without any leakage			REQ-059632
FenWal single blood pack X15 With citrate phosphate glucose adenine anticoagulant 500 ml, 15 units AFR0193/pack		15 EACH	REQ-054137
- K2 EDTA plastic blood collection tubes 100/pkt X05 10ml 16 x 100mm; 100 tubes, purple top - Microscope slides 100/pkt X04 Frosted, 25x75x1mm - Serum-gel plastic collection tubes 100/pkt X05 8 ml 16 x 100ml, 100/pk, yellow top			REQ-054138
75cm Cell culture flasks, 100/box X03 Tissue Culture flasks, rectangular, canted neck, vent cap, 3 x 100/box 175cm Cell culture flasks, 50/box X03 Flat flask,, angled neck, plug seal cap, sterile			REQ-055269
Type-II HRM PCR Kit 400 reactions)			
Features	Specifications	01 EACH	REQ-054034
Applications	Detection of SNPs, mutations and mutation sc		
Product use	Functionally validated and developed for reliat detection of genetic differences using HRM		
Reaction type	PCR amplification		
Enzyme activity	5'-> 3' Exonuclease activity		
With or without ROX	Without ROX		
Sequence specific Probe	Not necessary, EvaGreen dye for detection ind in the Mastermix		
Sample/target type	Genomic DNA		
With/without hotstart	With		
Mastermix	Yes		
Real-time or endpoint	Both		

RPMI1640 RPMI 1640 Medium is unique from other media because it contains the reducing agent glutathione and high concentrations of vitamins. RPMI 1640 Medium contains biotin, vitamin B12, and PABA, which are not found in Eagle's Minimal Essential Medium or Dulbecco's Modified Eagle Medium.	07 EACH	REQ-05403 2
DMEM/F-12 with GlutaMAX™ supplement DMEM/F-12 (Dulbecco's Modified Eagle Medium/Nutrient Mixture F-12). Concentration: 1X. Sterility: Sterile-filtered. Serum Supplementation: Standard Serum supplementation. Shipping Condition: Room Temperature. Volume: 500 ml. Form: Liquid. High Glucose, Sodium Pyruvate. With: GlutaMAX-I and Phenol Red. Without: HEPEs. DMEM/F-12 with GlutaMAX™ supplement minimizes toxic ammonia build-up and improves cell viability and growth. DMEM/F-12 contains no proteins, lipids, or growth factors, uses sodium bicarbonate buffer system and therefore requires a 5–10% CO2 environment to maintain physiological pH.	40 EACH	REQ-05407 4
DNeasy Blood & Tissue kits Applications PCR, real-time PCR, genotyping Elution volume 100–200 µl Time per run or per prep 20 minutes – 1 hour Main sample type Blood, tissue Format 96-well plate, spin column Sample amount 100 µl/25 mg Processing Manual Yield 6 µg/30 µg Technology Silica technology Purification of total RNA, miRNA, poly A+ mRNA, DNA or protein DNA	01 EACH	REQ-05402 6
- Universal Insect Replacement fly glue boards X10 - Replacement fly glue board X15 - Glue Board X01 - Brush On sticky Trap coating, 1L X 02		REQ-05382 4
Antibiotics Penicillin-Streptomycin-Amphotericin B Mixture 100ml	06 EACH	REQ-05402 8
- Tough-tags Microcentrifuge tube labels Autoclavable/freezable, Chemically resistant. Size: 0.94 in. x 0.5 in. wide, Fits 0.5 - 0.65 ml tubes - Tough-tags Microcentrifuge tube labels Autoclavable/freezable, Chemically resistant. Size: 1.28 in. x 0.5 in wide	01 EACH	REQ-05404 5

Fits 1.5 - 2 ml tubes		
175cm2 tissue culture treated flasks - vent cap X 05 175cm2 cell culture flask,vented cap,sterile.treated 75cm2 tissue culture treated flasks - vent cap X02 75cm2 tissue culture treated flasks – vented cap		REQ-054072
CLOSING DATE 19 DECEMBER 2021 @13:00		
KINDLY ENSURE THAT YOUR COMPANY IS REGISTERED ON CENTRALISED SUPPLIER DATABASE IF NOT ,KINDLY REGISTER ON THIS WEBSITE www.csd.gov.za		

- The above specified goods/services should be delivered/rendered to:

Name of Institute : ARC-OVI

Address : 100 OLD SOUTPAN ROAD(M35)

: ONDERSTEPOORDT

: PRIVATE BAG X05

:0110
- The particulars of the guarantee that will apply to the goods quoted for, with particular regard to the period and extent of the warranty must be clearly stated. Where services are required, service providers must submit documentation pertaining the relevant experience.
- Your written quotation should be sent to:

Mosidi Ramosunya

Tel :012 529 9315

Cell:

Fax:

E Mail: ramosunyam@arc.agric.za
- All price quotations that have a rand value of R30,000-00 to R50 000 000.00 including VAT, will be evaluated by applying the 80/20 principle as prescribed by the Preferential Procurement Policy Framework Act 5 of 2000 and its Regulations. The lowest acceptable price will score 80 points, the 20 BBBEE points will be allocated as follows:

B-BBEE Status Level of Contributor	80/20
1	20
2	18

3	14
4	12
5	8
6	6
7	4
8	2
Non-compliant Contributor	0

6. Standard conditions:

- 6.1 The validity of the quotations must be indicated.
- 6.2 Prices quoted should be in South African Rand and inclusive of VAT costs such as delivery, insurance, taxes, etc.
- 6.3 No price adjustments or amendment of the delivery particulars contained in paragraph 2 will be considered by the ARC.
- 6.4 The supplier accepts full responsibility for the proper execution and fulfilment of the goods/services quoted for.
- 6.5 ARC reserves the right to accept or reject any special terms and conditions that may qualify the goods/services to be provided.
- 6.6 Quotes should be submitted on an official letterhead and duly signed
- 6.7 Goods and services should be supplied/rendered upon receipt of a purchase order from the ARC
- 6.8 The General Conditions of Contract issued by National Treasury are applicable.
- 6.9 The ARC supply chain management code of conduct is applicable.
- 6.10 SBD Forms must be signed and returned together with the quotation if your price is above R10 000.00, failure to comply will result to disqualification of your quotation.
- 6.11 Only quotation from suppliers who are requested to quote shall be evaluated and considered.
- 6.12 Your quotation must indicate the delivery date.
- 6.13 The ARC reserve the right to do due diligence on the quotations.
- 6.14 The ARC reserve the right to benchmark prices quoted.

Thank you in anticipation

Supply Chain Management: ARC

Date: 09/10/2020